




# Smart fluid handling to take you forward, faster • • • •







#### WHERE TO BUY

Find your local distributor by visiting cpcworldwide.com or call CPC's Customer Service at 1-800-444-2474 or 651-645-0091. You can also send an email to info@cpcworldwide.com.

developing innovative high quality products that make fluid handling safe and easy.

#### QUALITY

CPC meets or exceeds our customers' expectations. Everyone is involved, from our suppliers to our distribution network, and most importantly, our employees. CPC measures and continually improves our standards of product quality, support services and overall customer and employee satisfaction. CPC's Quality System conforms to ISO 9001:2008 and ISO 13485 standards. Our cleanroom is ISO Class 7 certified.


#### **NEW PRODUCTS**

Exciting additions to CPC's product offering include:

NS1 Series Connectors	Page 24
• NS4 Series Connectors – Elbow	Page 60
• LQ4 Series Connectors	Page 62
• LQ6 Series Connectors	Page 66
BQ45GL Series Connectors	Page 86
Hybrid Connectors	Page 88
MPC/MPX Back-to-Back Insert Adapters	Page 109
• Genderless AseptiQuik® G Connectors	Page 122
AseptiQuik® G STC Connectors	Page 132

# **OUR CAPABILITIES**





CPC's engineered, custom solutions improve the overall functionality and design of the equipment and processes in which they are used. Special product features include precise hose barbs for superior grip, built-in shutoff valves to prevent product spills and an easy-to-use, push-button thumb latch for quick connecting and disconnecting.

CPC's patented Hybrid Connector technology platform eliminates misconnections by integrating electrical, fluid and gas lines into a single connection point. RFID-enabled IdentiQuik® couplings redefine what's possible for controlling, protecting and streamlining fluid handling processes. IdentiQuik technology allows the robust transfer of fluid and information to identify device type or connected media, verify line connection accuracy or capture fluid and product data. Specific applications for CPC connection technologies include:

#### • LIFE SCIENCES

CPC's products are used across a broad range of medical devices and equipment, including surgical, dialysis, blood pressure monitoring, patient therapy devices, and in vitro diagnostics. A de facto standard in several applications, CPC connectors improve patient safety, prevent misconnections and increase the overall usability and effectiveness of medical devices. Medical device OEMs around the world choose CPC connectors because of our proven performance in the operating room, the life science lab and the home healthcare setting.

CPC is the leader in single-use connection technology for bioprocessing applications, offering a wide variety of solutions including sterile connect and sterile disconnect. Our innovative designs offer the flexibility to easily combine multiple components and systems including process containers, tubing manifolds, transfer lines, bioreactors and other bioprocess equipment. CPC's robust connectors maintain media sterility and integrity while improving production yields, reducing costs and decreasing time to market for biopharmaceutical manufacturers. Genderless connector technology simplifies the design process, reduces supply chain complexity and eases manufacturing integration for biopharmaceutical manufacturers.

#### • SPECIALTY INDUSTRIAL

CPC's couplings and fittings are installed in thousands of product applications. Our broad offering appeals to product manufacturers in a wide range of industrial and consumer-based markets, including analytical instrumentation, printing and ink management, engine systems, biohazard detection equipment, electronic cooling and portable hydration systems.

#### • CHEMICAL HANDLING

CPC's extensive options for chemically resistant quick disconnect fitments, closures, couplings, and chemical dispensing systems provide non-spill operation and improve the safety of chemical management environments. Applications include semiconductor, LCD, solar, pharmaceutical manufacturing and chemical management, pump, filter and de-ionized water connections, cleaning, soap and detergent dispensing, bulk ink management and consumer-based applications.

Closed-loop dispensing systems enable fast, easy and economic chemical extraction from drums, IBCs and jerry cans. These systems increase workplace safety and protect workers and the environment from hazardous leaks and spills.





# ENGINEERED SOLUTIONS WORLDWIDE

CPC innovates on the leading edge of what is possible and collaborates to develop the best solution for each customer. Many of CPC's standard couplings and fittings began as custom engineered solutions that found broader use in applicable markets. CPC delivers unmatched product design and quality with local engineering support in the United States, Europe and Asia. To date, CPC has produced over 10,500 custom and standard products to meet the fluid handling challenges of various industries worldwide.



# INTRODUCTION



# **Global Distribution**

Find your local distributor by visiting cpcworldwide.com/ContactUs or call CPC's Customer Service at 1-800-444-2474 or 651-645-0091.

# Quality

Quality is a key component of CPC's success. We are ISO 9001:2008 and ISO 13485:2003 certified and our cleanroom is ISO Class 7 certified.

# How to Choose the Right Connector

#### INTRODUCTION

With so many connection options, it can often be an overwhelming task to decide which connector is best suited for an application. By understanding your application requirements and selecting the correct connection type, you will have better performance and sealing results.

#### **ASSESSING THE APPLICATION**

Understanding your application is the key to selecting the proper connection. Use the following guide to simplify your selection.

# **FLOW**

TUBING MEDIA

**PRESSURE** 

TUBING

SHUTOFF

OPTIONS

MOUNTING

OPTIONS

SPECIAL REQUIREMENTS

CONNECTIONS

for the effect of shutoff valves and tubing connections on your calculations.

What size tubing, both inner and outer diameter, are you using?

What is your required flow and pressure drop? Be sure to allow

The viscosity and corrosiveness of the fluid going through the connection need to be considered. Make sure the media is chemically compatible with ALL coupling materials - including the seals, o-rings and springs. Double-check media compatibility.

**TEMPERATURE**Know your minimum and maximum temperature range.
Standard temperature ranges are from -40° to 200° F
(-40° to 93° C) depending on connection material.

fittings and tube ID and OD for compression.

What is the maximum pressure your connection will need to withstand during operation? Quick disconnects rated to 250 psi (17 bar) will handle most low pressure applications.

Type: Hose barb, compression fittings and push-to-connect are the most common termination styles.

Size: You will need to know ID for hose barbs, OD for push-in-

Do you need automatic or integral shutoff valves? Shutoff options are single, double and non-spill.

How is the connection going to be configured into your application? Common mounting options include pipe thread, panel mount, in-line or elbow.

Sterilization, NSF listed, USP Class VI approved materials, special packaging, color coding, keying, lot traceability, etc.

# Visit us online for interactive web tools











# **ALSO ON THE WEB**

- Part Specifier
- Chemical Compatibility
- Ask Our Engineers
- CAD Models
- Watch Product Videos



#### cpcworldwide.com

Quick Disconnect Couplings and Connectors.







# Material Options

The type of media flowing through a connection can affect the strength, surface appearance, color and performance of the connection. Here are some guidelines for the different types of material.

#### **THERMOPLASTICS**

ABS

Economical, medical-grade thermoplastic that withstands gamma and e-beam sterilization. It is commonly used in medical devices. ABS is an amorphous material with good physical properties and high resistance to chemical attack.

ACETAL

Strong, lightweight and economical material used for a wide variety of chemical and mechanical components. Acetal offers high strength and rigidity over a broad temperature range, low wear, toughness and resistance to repeated use.

**POLYAMIDE** (NYLON)

Very resistant to wear and abrasion, good mechanical properties even at elevated temperatures, low permeability to gases and good chemical resistance.

PEEK (POLYETHERETHERKETONE)

Highly temperature resistant, engineered thermoplastic with excellent chemical and fatigue resistance. It exhibits superior mechanical and electrical properties.

**POLYCARBONATE** 

Resistant to some chemicals, withstands sterilization and is typically transparent. It is commonly used in medical devices and offers impact resistance, outstanding dimensional stability and good optical properties.

**POLYETHYLENE** 

Low-cost, chemically resistant thermoplastic. It is opaque and can withstand reasonably high temperature.

**POLYPROPYLENE** 

Excellent general purpose resin that is highly resistant to chemical attack from solvents and chemicals in harsh environments. In general, polypropylene is resistant to environmental stress cracking and it can be exposed to challenging environments.

**POLYSULFONE** 

Rigid material with excellent strength, good chemical resistance; withstands repeated sterilization and higher temperatures than other thermoplastics. Its high hydrolytic stability allows its use in medical applications requiring repeated high temperature autoclave and steam sterilization.

PPS

(POLYPHENYLENE SULFIDE)

Broadest resistance to chemicals as an advanced engineering plastic.

#### **FLUOROPOLYMERS**

PTFE

Chemically resistant to all chemicals and solvents with the exception of some molten metals, molten sodium hydroxide, elemental fluorine and certain fluorinating agents. PTFE offers chemical resistance and stability at high temperature.

**PVDF** 

Tough engineering thermoplastic with a balance of physical and chemical properties that qualifies it for high performance in a wide range of applications. It is mechanically strong and tough, has good ductility and has a broad, useful temperature range.

#### **ALLOYS**

**ALUMINUM** 

Lightweight metal with an available hard anodized finish for durability. Aluminum is non-toxic, non-magnetic, non-sparking and is known for its high strength to weight ratio.

CHROME-PLATED BRASS

Rugged metallic material with an attractive appearance, chrome-plated brass is excellent for higher pressure and temperature.

**DIE-CAST ZINC** 

Durable and lightweight (about 20% less than comparable brass) material that withstands high pressure and temperature.

**STAINLESS** STEEL

Offers excellent rust resistance and is often used for connector components such as valve springs.

O-RING SELECTION

Selecting the correct o-ring is directly related to your application. One o-ring can offer you better chemical resistance, another can offer better heat resistance, a third can offer cold flexibility. Here are some guidelines for the different types of o-rings.

BUNA-N

The most common o-ring material is Buna-N due to its solvent, oil and water resistance. It has a temperature range of -40° to 250° F (-40° to 121° C).

**EPDM** 

Ethylene-propylene-diene rubber (EPDM, also sometimes referred to as EPR) is a chemically-resistant family of compounds. CPC uses high quality peroxide-cured EPDMs that provide exceptional resistance to temperatures when using a wide range of chemicals.

FFKM (CHEMRAZ®, SIMRIZ®, KALREZ®)

Broadest range of chemical resistance of any elastomeric material, combining the resilience and sealing force of an elastomer with chemical resistance approaching that of PTFE.

FKM

Well-known for its outstanding resistance to heat, oxidation, weathering and ozone. The temperature range is -15° to 400° F (-26° to 204° C).

PFA & FEP **ENCAPSULATED** SEALS SILICONE

Encapsulated seals combine the resiliency of the elastomer with the superior chemical resistance of the fluoropolymer to achieve a seal that is lower cost than a pure fluoroelastomer FFKM seal.

Seals have good temperature resistance. The temperature range is -70°F to 400°F (-57°C to 204°C) with special compounds that can reach up to 450°F (232°C). Silicone can also be supplied with Class VI requirements for life sciences applications.

# THE RIGHT NUMBER OF BARBS

Many things work together





to determine the quality of the connection between a hose barb and the tubing it connects. The sharpness of the barb, surface finish and the barb angle all contribute to the overall quality of the connection. Failing to optimize these technical aspects will result in a poor connection, regardless of the number of hose barbs

CPC offers a variety of hose barb styles and connection options for tubing ranging from 1/16" to 1" ID.

# **ELASTOMERS**

TPE

Thermoplastic Elastomer (TPE) is a blend of additives and copolymers in a special formulation that forms extremely durable bonds to a substrate, while offering the traditional properties of soft-touch overmold.

TPV

Thermoplastic Vulcanizate (TPV) is an alloy of polypropylene thermoplastic and fully vulcanized EPDM rubber. TPV is typically resistant to water, acids and bases.

6 cpcworldwide.com ● 800-444-2474

# **CUSTOM CAPABILITIES**

# **Connector solutions**

from CPC add value to your product by making your fluid and air connections easy to use and more reliable, increasing your product's modularity and serviceability and providing an overall cleaner, faster, safer and smarter way to make a connection. If you can't find the perfect solution from one of our thousands of standard products or you'd like to simply off-load the connector portion of your project to allow you to focus on your core technologies, CPC Engineered Solutions can help.

#### HYBRID COUPLINGS

Connect both fluid and electronics – at the same time in a single connection.



#### RFID-ENABLED INTELLIGENT FLUID CONNECTIONS

Redefine what's possible for controlling, protecting and streamlining fluid handling processes. Identify misconnections, control usage and protect your brand.



#### METAL NON-SPILL COUPLINGS



#### **CUSTOM MATERIALS**

A wide range of materials are available to meet application needs.



#### **MULTI-TUBE COUPLINGS**

Connect multiple lines with one easy connection.



#### **CUSTOM TUBE CONNECTIONS**

Connect to a wide variety of types and sizes of tubing.



are exclusively produced These proprietary products Contact your local distributor for more

#### **MOUNTING OPTIONS**

Application-specific connections including threads, caps, snap-fits and more.



# Consider a Custom-Designed Connector When:

- >> A quick disconnect will add value to your product, make it easier to use and more reliable
- >> Your specifications cannot be met by an existing standard CPC product
- >> Unique requirements, budgets or timing warrant your designer's collaboration with CPC's Engineered Solutions Team

Manifold







Polycarbonate



Metal non-snill



Cap with vent



Industrial



2-way valve



Dual-line connector





Air/fluid/electronic

tube set





RFID panel In-line filter



8 cpcworldwide.com • 800-444-2474

#### **LEGEND**




Shut-Off

Note: Product images may not be to scale

Shut-Off

# General Purpose Couplings






alternative to luer-type connections. Material: Acetal, polypropylene, ABS, chrome-plated brass **Tubing ID Sizes:** 1/16" to 1/8" (1.6mm to 3.2mm)

Through






**Tubing ID Sizes:** 3/32" to 1/8" (2.4mm to 3.2mm)



PAGE | DESCRIPTION



Material: Glass-filled polypropylene, PEEK®, polypropylene
Tubing: Microbore to 1/8" (Microbore to 3.2mm)





**Tubing ID Sizes:** 1/16" to 1/4" (1.6mm to 6.4mm)



**30** PMC12: Offered with a variety of configurations and chemical resistance for demanding applications; gamma sterilizable. **Material:** Polypropylene

**Tubing ID Sizes:** Microbore to 1/4" (Microbore to 6.4mm)



**34** MC: Durable and able to withstand higher pressure and temperature; easy one-hand connection and disconnection.

**Material:** Chrome-plated brass

**Tubing ID Sizes:** 1/8" to 1/4" (3.2mm to 6.4mm)




**38** NS212: Twist-to-connect design features non-spill valves designed to provide fast, safe and virtually leak-free fluid line connections. Material: Glass-filled polypropylene
Tubing ID Sizes: 1/8" to 1/4" (3.2mm to 6.4mm)



**40 PLC:** Widest selection of sizes and configurations offered; resistant to most mild chemical solutions. **Material:** Acetal

Tubing ID Sizes: 1/4" to 3/8" (6.4mm to 9.5mm)





44 PLC12: Materials of construction offer broad chemical resistance for demanding applications; gamma sterilizable. Material: Polypropylene

Tubing ID Sizes: 1/4" to 3/8" (6.4mm to 9.5mm)



# General Purpose Couplings

#### PAGE | DESCRIPTION



**48** LC: Durable and able to withstand higher pressure and temperature; easy one-hand connection and disconnection.


Material: Chrome-plated brass **Tubing ID Sizes:** 1/4" to 3/8" (6.4mm to 9.5mm)

 $\rightarrow + \leftarrow \rightarrow + \leftarrow \rightarrow + \leftarrow$ 



**52 APC:** Features plastic thumb latch, fewer moving parts and a smooth contoured design to deliver ease-of-use and excellent flow.

Tubing ID Sizes: 1/4" to 3/8" (6.4mm to 9.5mm)





**56** BREAKAWAY®: Provides safe and easy fluid transfer with protection from costly product loss and equipment damage.

**Material**: Acetal

Tubing ID Sizes: 1/4" and 3/8" (6.4mm and 9.5mm)





**58 EFC12:** High-efficiency valve design provides high flow capability; bulkhead panel mount option facilitates tight seals against tank walls

**Material:** Polypropylene

**Tubing ID Sizes:** 1/4" and 3/8" (6.4mm and 9.5mm)





**60** NS4: Non-spill coupling that virtually eliminates spills, minimizes downtime and enhances operator safety.

Material: Glass-filled polypropylene, ABS **Tubing ID Sizes:** 1/8" to 3/8" (3.2mm to 9.5mm)

**≫**-«⊙



**62** LQ4: Liquid cooling non-spill coupling guarantees a secure, reliable connection and dripless disconnect.

**Material**: Chrome-plated brass **Tubing ID Sizes:** 1/4" to 3/8" (6.4mm to 9.5mm)




**64** NS6: Durable, yet lightweight construction that features non-spill valves and is compatible with many chemicals. Material: Glass-filled polypropylene
Tubing ID Sizes: 3/8" and 1/2" (9.5mm and 12.7mm)



# **LEGEND**

Shut-Off



Note: Product images may not be to scale

Shut-Off

# General Purpose Couplings



PAGE | DESCRIPTION

**66** LQ6: Liquid cooling non-spill coupling guarantees a secure, reliable connection and dripless disconnect. *Material*: Chrome-plated brass **Tubing ID Sizes:** 3/8" to 1/2" ID (9.5mm to 12.7mm ID)

Through





**68** NSH: Pressure-balanced non-spill design with 100% metal-free, spring-free flow path. Material: Glass-filled polypropylene
Tubing ID Sizes: 3/8" to 3/4" (9.5mm to 19.0mm)





**70 HFC12:** Efficient valve design leads to high flow and exceptionally low spillage; shrouded thumb latch is easy to grip and simple to operate. **Material:** Polypropylene

**Tubing ID Sizes:** 3/8" to 3/4" (9.5mm to 19.0mm)





74 HFC35 & 57: Withstand harsh environments and offered with or without UV-resistant materials to withstand harmful rays without affecting

*Material:* Polysulfone (white), UV-resistant polysulfone (black) Tubing ID Sizes: 3/8" to 3/4" (9.5mm to 19.0mm)











**78 FFC35:** Features non-valved and unobstructed flow path to increase flow and minimize turbulence.

**Material:** Polysulfone Tubing ID Sizes: 3/4" (19.0mm)

# Specialty Products



PAGE | DESCRIPTION

**80 TENTUBETM:** Allows connection and disconnection of up to ten lines with one slide latch; tubing orientation ensured by physical keying. Material: Nylon, acetal, polypropylene Tubing ID Sizes: 1/16" to 1/8" (1.6mm to 3.2mm)

 $\rightarrow \leftarrow \rightarrow \leftarrow$ 

# Specialty Products (cont.)



**82 SIXTUBETM:** Snap-in panel mount design and flexibility to connect and disconnect six separate lines with or without valves. Material: Acetal, polypropylene

**Tubing ID Sizes:** 1/16" to 1/8" (1.6mm to 3.2mm)

 $\rightarrow \leftarrow \rightarrow \leftarrow$ 

TWIN TUBETM: One easy-to-use quick disconnect for two separate

Material: Acetal, ABS

**Tubing ID Sizes:** 1/16" to 1/8" (1.6mm to 3.2mm)

 $\rightarrow \leftarrow$ 



**84** MULTI-MOUNT: Connects three to five lines at once; keyed to prevent

**Material:** Acetal, chrome-plated brass

**Tubing ID Sizes:** 1/8" to 3/8" (3.2mm to 9.5mm)

 $\rightarrow \leftarrow \rightarrow \leftarrow \rightarrow \leftarrow$ 



**86 BQ45GL:** Quick disconnect coupling or fitting integrated into a cap designed for 45GL bottles. *Material:* Polyethylene

 $\rightarrow \leftarrow \rightarrow \leftarrow$ 







**HYBRID CONNECTOR:** One integrated connector combining 2 fluid lines with non-spill couplings and 4 electrical lines to simplify the user

Material: Glass-filled polypropylene



**90** MINI HYBRID CONNECTOR: A compact, single connection point combining 1 fluid line with 2 electrical lines with a breakaway mechanism for ease of use.

Material: Glass-filled polypropylene

 $\rightarrow + \leftarrow$ 






PAGE | DESCRIPTION

**92** INTELLIGENT FLUID CONNECTIONS: Redefine what's possible for controlling, protecting and streamlining fluid handling with RFID-enabled couplings.

Material: Acetal, polypropylene, glass-filled polypropylene












# **LEGEND**



Through Shut-Off Shut-Off Note: Product images may not be to scale

# Life Sciences - Medical & Bioprocessing



**102** SMC: Twist-to-connect design manufactured and packaged in a cleanroom. *Material:* Medical-grade polycarbonate **Tubing ID Sizes:** 1/16" to 1/8" (1.6mm to 3.2mm)

 $\rightarrow \leftarrow$ 

**104** SRC: A unique small bore connector that eliminates the potential for misconnections with luer fittings. Material: Medical-grade polypropylene **Tubing ID Sizes:** 1/8" to 3/16" (3.2mm to 4.8mm)

 $\rightarrow + \leftarrow$ 



**106** MPC: Easy-to-use and secure connection for critical fluid applications; includes pressure sealing caps and plugs with optional locking sleeves.

\*\*Material:\*\* Medical-grade ABS, medical-grade polycarbonate and polysulfone\*\*

**Tubing ID Sizes:** 1/8" to 3/8" (3.2mm to 9.5mm)

 $\rightarrow +\leftarrow$ 



**109** BACK-TO-BACK ADAPTERS: Allow end users to connect single-use systems that may feature identical connections at the end of their tubing. Material: Medical-grade polycarbonate and polysulfone, USP Class VI, ADCF



110 MPX: Larger flow, easy-to-use and secure connection for critical fluid applications, includes pressure sealing caps and plugs with optional locking

Material: Medical-grade polycarbonate and polysulfone, USP Class VI, ADCF **Tubing ID Sizes:** 3/8" to 1/2" (9.5mm to 12.7mm)



112 SANITARY SERIES: Attaches directly to 3/4", 1" and 1-1/2" sanitary terminations to provide greater flexibility for integrating components into single-use or hybrid process systems.

Material: Medical-grade polysulfone, USP Class VI, ADCF Termination Sizes: 3/4", 1" and 1-1/2" sanitary

 $\rightarrow \leftarrow$ 



**114 SANIQUIKTM:** Integral sanitary termination attaches to hard-plumbed systems with tri-clover clamps; permits quick and easy connection to single-use bag systems, manifolds or tube sets. Material: 316L stainless steel

Termination Sizes: 3/4" and 1-1/2" sanitary



116 MPU: Larger flow twist-to-connect design features easy-to-use locking mechanism that guards against accidental disconnects. Material: Medical-grade polysulfone, USP Class VI, ADCF **Tubing ID Sizes: 3/4" (19.0mm)** 

# Life Sciences - Medical & Bioprocessing





118 HFC39: Features sterile disconnect functionality with automatic shutoff valve, preventing external organisms from entering the media flow path upon

Material: Medical-grade polysulfone, USP Class VI, ADCF Termination Sizes: 1/4", 3/8" and 1/2" HB (6.4mm, 9.5mm, 12.7mm)

 $\rightarrow \leftarrow \rightarrow \leftarrow$ 



**120** ASEPTIQUIK® S: Genderless design provides quick and easy sterile connections for small-flow applications.

Material: Medical-grade polycarbonate, USP Class VI, ADCF

Termination Sizes: 1/8", 1/4" and 3/8" HB (3.2mm, 6.4mm and 9.5mm), 3/4" sanitary and MPC insert

 $\rightarrow +\leftarrow$ 



**122** ASEPTIQUIK® G: Robust genderless connectors provide quick and easy sterile connections and simplify systems integration.

Material: Medical-grade polycarbonate, USP Class VI, ADCF **Termination Sizes**: 1/4", 3/8", 1/2", 3/4" ID hose barb (6.4mm, 9.5mm, 12.7mm, 19.0mm) and 3/4" sanitary

 $\rightarrow \leftarrow$ 



124 ASEPTIQUIK® C: Provides a quick and easy sterile connection, even in

Material: Medical-grade polycarbonate, USP Class VI, ADCF Termination Sizes: 3/8" and 1/2" HB (9.5mm and 12.7mm) and 3/4" sanitary



**126** ASEPTIQUIK® DC: All-in-one, single-use connection technology offering both a sterile connect and a sterile disconnect.

Material: Medical-grade polycarbonate, USP Class VI, ADCF Termination Sizes: 1/4". 3/8" and 1/2" HB (6.4mm, 9.5mm, 12.7mm)

 $\rightarrow +\leftarrow$ 



**128** ASEPTIQUIK® X: Robust large format design enables quick sterile transfer for high flow single-use applications.

Material: Medical-grade polycarbonate, USP Class VI, ADCF Termination Sizes: 3/4" and 1" ID hose barb (19.0mm and 25.4mm) and 1-1/2" sanitary

 $\rightarrow +\leftarrow$ 



130 STEAM-THRU® CONNECTIONS: Allow quick and easy sterile connection via SIP between biopharmaceutical processing equipment and disposable bag and tube assemblies.

Material: Medical-grade polysulfone, USP Class VI, ADCF
Termination Sizes: 3/8" and 1/2" HB (9.5mm and 12.7mm) and 3/4" sanitary



132 ASEPTIQUIK® STC: Allows a gendered or genderless AseptiQuik sterile connection to be steamed on to stainless equipment via SIP. Material: Medical-grade polycarbonate and polysulfone, USP Class VI, ADCF

 $\rightarrow \leftarrow$ 

**14** cpcworldwide.com • 800-444-2474

#### **LEGEND**



Note: Product images may not be to scale

# High Purity Couplings

PAGE | DESCRIPTION

#### 138 CHEMQUIK® CQH06/CQV06: All plastic design for ultra-pure media and high chemical resistance.

Through

Material: Natural, virgin polypropylene and PVDF Tubing Sizes: 1/4" to 1/2" flare; 3/8" to 3/4" NPT



# **140 CHEMQUIK CQGO6:** Spring-free and metal-free flow path for high flow capacity and non-spill, high purity connections.

Material: Natural, virgin polypropylene Tubing Sizes: 3/8" to 3/4" flare; 3/8" to 3/4" NPT





# CHEMQUIK DUAL CONTAINMENT SYSTEM: Flare nuts

and panel mount fittings that facilitate double containment of critical chemical

Material: Natural, virgin polypropylene

# **Dispensing Connectors**

PAGE | DESCRIPTION

Thread Sizes:



# 142 DRUMQUIK® PRO & DRUMQUIK PUR: Closed chemical dispensing system with a reusable coupler and a recyclable dip-tube assembly

for chemical extraction from drums, jerry cans and IBCs. *Material:* Food grade, virgin polypropylene and polyethylene

Drum Inserts: 2" American buttress, BCS 56x4 and 2" NPS (G8), 70mm for BCS 70x5 and BCS70x6

Couplers: 1/2" NPT, 3/4" hose barb and 3/4" BSPP

# 145 DRUMQUIK ADAPTORS AND ACCESSORIES: CPC

offers many accessories intended for use with DrumQuik® PRO and DrumQuik PUR systems. These include vent check valves, foot valves, fittings and more.

# Dispensing Connectors (cont.)

PAGE | DESCRIPTION

(CON-18)

# **148** DRUMQUIK® 3-PORT UDA: Used with 3/4" male NPT coupling terminations to adapt dip-tubes to standard drum closures for closed system

Material: Food grade, virgin polypropylene
Thread Sizes: UDA 3-Port: 3/4" male NPT; UDA: none

#### **BOTTLE ADAPTOR KITS:** Intended for use on reagent

bottles with SP400-38mm threads. Material: Natural, virgin polypropylene Thread Sizes: SP400-38mm

# DRUMQUIK ASIAN DRUM ADAPTORS, PLUGS AND CAPS

**150** UDC: Provides closed connection to bag-in-box packaging with 38mm threaded and snap-in necks.

**Material:** Polypropylene and acetal

**Tubing ID Sizes:** 3/8" to 3/4" (9.5mm to 19.0mm)





# **152** PUNCTURE SEAL: Provides closed connection to Hedwin Cubitainer® bag-in-box container systems and flexible bags or bottles with SP400 38mm

Material: Polypropylene (coupling) and polyethylene (cap) Tubing ID Sizes: 1/4" and 3/8" (6.4mm and 9.5mm)



# Fittings, Luers, and Blood Pressure Connectors

PAGE | DESCRIPTION



156 FITQUIK® CONNECTORS: High-quality fittings for leak-free tubing connections. These precision-molded fittings are designed to eliminate tubing leak points in applications such as medical devices, analytical instrumentation or air-driven equipment.

# Additional Products and Information

PAGE | DESCRIPTION



**162** ACCESSORIES: Offering includes a variety of ancillary components from fittings and luers to PTF nuts and dust caps.


163 STERILIZATION AND DISINFECTANT METHODS


164 REGULATORY AND COMPLIANCE



The SMC and SMF1 are among CPC's smallest couplings. These twistto-connect couplings provide a reliable and more secure alternative to luer-type connections. They also allow for the tubing to rotate freely when connected. This important feature prevents both kinked tubing and accidental disconnection during use. The SMC Series is also available with optional RFID (Radio Frequency Identification) capability (page 98) and in polycarbonate material (page 102).

FEATURES	BENEFITS
Twist to connect	Prevents accidental disconnects
Free coupling rotation	Eliminates kinked tubing
Quarter turn latch clicks when connected	Prevents damage due to over torquing





# Specifications • • •

#### PRESSURE:

Vacuum to 100 psi, 6.9 bar

#### TEMPERATURE:

Acetal, Chrome-Plated Brass: -40°F to 180°F (-40°C to 82°C)

#### Polypropylene:

32°F to 180°F (0°C to 82°C)

-40°F to 160°F (-40°C to 71°C)

#### MATERIALS:

Main Components: ABS, acetal, polypropylene, chrome-plated brass, stainless steel

Locking sleeves: Acetal

Valves: Acetal

Valve spring: 316 stainless steel

**0-rings:** Buna-N with acetal or chrome-plated

brass, EPDM with polypropylene

#### COLOR:

Main components: Natural white (acetal and ABS), almond (polypropylene), chrome and black

# **TUBING SIZES:**

1/16" to 1/8" ID, 1.6mm to 3.2mm and 5.0mm

WARNING: Pressure, temperature, chemicals, and operating environment can affect the performance of couplings. It is the customer's responsibility to test the suitability of CPC's products in their own application conditions.



These graphs are intended to give you a general idea of the performance capabilities of each product line. The shaded area of each graph represents the operating range of the product family, i.e., upper and lower values are shown. Therefore, depending on the exact coupling configurations selected, you can reasonably expect values to fall within the

**Liquid Flow Rate Information for Couplings**The chart below shows the flow rate for CPC couplings. Each coupling was tested with water at 70°F (21°C). To determine flow rates for specific coupling configurations use the formula below.

# C<sub>v</sub> VALUES FOR SUBMINIATURE COUPLINGS


BODIES INSERTS	SMM01	SMM02
SMF01	.03	.03
SMFD01	.03	.03
SMF02	.03	.19
SMFD02	.03	.08
SMPT02	.03	.19
SMFD02	.03	.08


Q=C <sub>v</sub>	ΔP
v-V	S

- Q = Flow rate in gallons per
- $\mathbf{C}_{\mathbf{v}} = \text{Average coefficient}$ across various flow rates (see chart)
- $\Delta P$ = Pressure drop across coupling (psi)
- S = Specific gravity of liquid

VOTES:			
	4041414141		
		MAKAKA	
		44444	

# **SMC DIMENSIONS**





# **Panel Dimensions**

	PANEL OPENING	PANEL THICKNESS MAX. – MIN.	PANEL NUT HEX	PANEL NUT THREAD
COUPLING BODIES	7/16	.21 – .03	1/2	7/16-24UNF
COUPLING INSERTS	7/16	.21 – .03	1/2	7/16-24UNF

# **Coupling Bodies** • ACETAL



PIPE	THREAD
IN-LII HOSE	NE BARB

**TERMINATION** 

IN-LINE





IEL M	•••••	

<i>Tubing/Thread Size</i> 1/8" NPT 1/8" NPT	METRIC EQ.	STRAIGHT THRU SMPT02 SMPT02BLK	SHUTOFF SMPTD02 SMPTD02BLK	.51 .51	.86 .86	7/16 7/16
1/16" ID 1/16" ID 1/8" ID 1/8" ID	1.6mm ID 1.6mm ID 3.0mm ID 3.0mm ID 3.2mm ID 3.2mm ID	SMF01 SMF01MBLK SMFM3 SMFM3MBLK SMF02 SMF02MBLK	SMFD01 SMFD01MBLK SMFDM3 SMFDM3MBLK SMFD02 SMFD02MBLK	.48 .48 .48 .48 .48	.75/.90 .75/.90 1.00 1.00 .90	
1/16" ID 1/16" ID 1/8" ID 1/8" ID	1.6mm ID 1.6mm ID 3.0mm ID 3.0mm ID 3.2mm ID 3.2mm ID	SMFPM01 SMFPM01BLK SMFPMM3 SMFPMM3BLK SMFPM02 SMFPM02BLK	SMFPMD01 SMFPMD01BLK SMFPMDM3 SMFPMDM3BLK SMFPMD02 SMFPMD02BLK	.58 .58 .58 .58 .58	.75/.90 .75/.90 1.00 1.00 .90	1/2 1/2 1/2 1/2 1/2 1/2 1/2
1/16" ID 1/8" ID	1.6mm ID 3.2mm ID	iSMFT0103 iSMFT0203	iSMFDT0103 iSMFDT0203	.76 .76	.76/.91 .91/.91	

#### POLYPROPYLENE

HOSE BARB



TERMINATION IN-LINE PIPE THREAD HOSE BARB (Acetal Locking Sleeve)





TUBING/THREAD SIZE 1/8" NPT	METRIC EQ.	STRAIGHT THRU SMPT0212	.51	.86	<b>C</b> 7/16
1/16" ID	1.6mm ID	SMF0112	.48	.75/.90	······································
1/8" ID	3.2mm ID	SMF0212	.48	.90	
1/16" ID	1.6mm ID	SMFPM0112	.58	.75/.90	1/2
1/8" ID	3.2mm ID	SMFPM0212	.58	.90	1/2
1/16" ID	1.6mm ID	iSMFT0100	.76	.76	•••••••••••
1/8" ID	3.2mm ID	iSMFT0200	.76	.91	

# SMF1 Disposable Bodies • ABS



TERMINATION IN-LINE HOSE BARB

TUBING/THREAD SIZE	METRIC FO	STRAIGHT THRU	A	
		***************************************		
1/8" ID	3.2mm ID	SMF10297	.38	.8
1/8" ID	3.2mm ID	SMF10297MBLK	.42	.8
	5mm ID	SMF1M597	.42	.6
	5mm ID	SMF1M597MBLK	.42	.6

CPC has added a proprietary lubrication step to our SMF1 manufacturing process. This revolutionary technique allows the lube to be placed only where you need it, ensuring consistent connection forces throughout the lifetime of the connectors even when aggressive cleaning techniques are used.

# **Coupling Inserts**


# ACETAL



#### • CHROME-PLATED BRASS



#### • POLYPROPYLENE



# SMC Complete In-Line Couplings

All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters. MBLK = molded black material. BLK = dyed black material.

#### ACETAL



TUBING/THREAD SIZE	METRIC EQ.	STRAIGHT THRU	SHUTOFF	B
1/16" ID	1.6mm ID	SMC01	SMCD01	1.32/1.
1/16" ID	1.6mm ID	SMC01MBLK	SMCD01MBLK	1.32/1.4
	3.0mm ID	SMCM3	SMCDM3	1.62
	3.0mm ID	SMCM3MBLK	SMCDM3MBLK	1.62
1/8" ID	3.2mm ID	SMC02	SMCD02	1.61
1/8" ID	3.2mm ID	SMC02MBLK	SMCD02MBLK	1.61

# POLYPROPYLENE



TUBING/ THREAD SIZE	METRIC EQ.
./16" ID	1.6mm ID
./8" ID	3.2mm ID

STRAIGHT THRU
SMC0112
SMC0212

# SNAPQUIK SERIES CONNECTOR

# Specifications • • •

#### PRESSURE:

Vacuum to 15 psi, 1.0 bar

#### TEMPERATURE:

32°F to 122°F (0°C to 50°C)

#### MATERIALS:

Main components: Acetal

O-ring: Buna-N

# COLOR:

Natural White

#### **TUBING SIZES:**


3/32" and 1/8" ID


WARNING: Pressure, temperature, chemicals, and operating environment can affect the performance of couplings. It is the customer's responsibility to test the suitability of CPC's products in their own application conditions.

# SnapQuik® Series Connectors are ideal for small bore applications. The compact design delivers high flow fluid transfer in a small profile. Featuring a unique internal latching mechanism that creates a smooth, streamlined exterior that is easy to clean and use. The user-friendly break away design of SnapQuik prevents misconnections with luers, protects equipment and ensures patient safety.

#### **FEATURES BENEFITS** Break away latch Intuitive, rapid disconnection Internal latch design Streamlined exterior is easy to clean Audible click at connection Assures reliable connection is made Does not mate to luer fittings Improved patient safety


he lifetime of the connectors, cleaning techniques are





These graphs are intended to give you a general idea of the performance capabilities of each product line. The shaded area of the graph represents the operating range of the product family, i.e. upper and lower values are shown. Therefore, depending on the exact coupling configurations selected, you can reasonably expect values to fall within the shaded area.

# SNAPQUIK DIMENSIONS





# **Coupling Bodies • ACETAL**



**TUBING SIZE** 3/32" ID 1/8" ID

METRIC EQ. 2.4mm 3.2mm

BSQ17015 BSQ1702

# **Coupling Inserts • ACETAL**



**TUBING SIZE** 3/32" ID 1/8" ID

METRIC EQ. 2.4mm 3.2mm

STRAIGHT THRU BSQ22015 BSQ2202

A 0.42 0.42

0.36 0.36

B 1.14 1.14

B 1.09

1.09

All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters.

# DID YOU KNOW?

#### AN ALTERNATIVE TO LUERS

An alternative to luers, SnapQuik® will not mate with luers designed to the ISO 594 standards



# NS1 SERIES CONNECTOR

# Specifications • • •

#### PRESSURE:

Vacuum to 15 psi, 1.0 bar

#### TEMPERATURE:

32°F to 120°F (0°C to 49°C)

#### MATERIALS:

#### Main Components:

Glass-filled polypropylene, PEEK®, polypropylene

Thumb latch: PEEK®
Springs: 316 stainless steel

COLOR: Gray

O-rings: EPDM

#### **TUBING SIZES:**

Microbore to 1/8" ID (Microbore to 3.2mm)

#### SPILLAGE:

0.02 cc per disconnect rated @ 15 psi


**INCLUSION:** 0.01 cc per connect

**WARNING**: Pressure, temperature, chemicals, and operating environment can affect the performance of couplings. It is the customer's responsibility to test the suitability of CPC's products in their own application conditions.


**NS1 Series couplings** compact design is packed full of features including non-spill shutoff valves that ensure drip-free disconnects, and eliminate the need to purge air from the system. The smallest plastic non-spill coupling on the market, NS1 is perfect for applications with tubing from 1/8-inch and smaller. Its rotating design allows tubing to move freely without kinks. Ideal for applications that need to be connected and disconnected often, the ergonomic, user-friendly interface enables a secure, reliable connect and disconnect.

FEATURES	BENEFITS
Non-spill design	Closed system for drip-free disconnect
Ultra low air inclusion	Minimizes contaminants, eliminates the need to purge air bubbles from the system
Intuitive interface	User friendly ergonomic design enables a secure, reliable connection
Rotating design	Allows tubing to move freely or relax without kinking or leaking
Chemically resistant materials	Reliable strength and performance in a range of applications

➤ This graphs is intended to give you a general idea of the performance capabilities of each product line. The shaded area of each graph represents the operating range of the product family, i.e., upper and lower values are shown. Therefore, depending on the exact coupling configurations selected, you can reasonably expect values to fall within the shaded area.



# **NS1 SERIES DIMENSIONS**




# Panel Dimensions

	PANEL	PANEL THICKNESS	PANEL	PANEL NUT
	OPENING	MAX. – MIN.	NUT HEX	THREAD
COUPLING INSERTS	see drawing	.2905	5/8	1/2-24UNS



# **Coupling Bodies • POLYPROPYLENE**



# **Coupling Inserts • POLYPROPYLENE**

	DESCRIPTION IN-LINE HOSE BARB	<b>TUBING SIZE</b> 1/8" ID	METRIC EQ. 3.2 MM ID	SHUTOFF NS1D220212	A 0.44	B 1.04
O THE OWNER OF THE OWNER OWNER OF THE OWNER	IN-LINE WITH A 1/4-28 FLAT BOTTOM PORT	1/8" OD OR SMALLER*	3.2 MM OD OR SMALLER*	NS1D26042812	0.44	0.92
E III	IN-LINE WITH A 1/4-28 MALE THREAD			NS1D24042812	0.49	0.99
	PANEL MOUNT HOSE BARB	1/8" ID	3.2 MM ID	NS1D420212	0.70	1.56
	PANEL MOUNT WITH A 1/4-28 FLAT BOTTOM PORT	1/8" OD OR SMALLER*	3.2 MM OD OR SMALLER*	NS1D48042812	0.70	1.45

\*Tubing size dependent on nut and ferrule used with 1/4-28 port.

# Accessories



MATERIAL EPDM

PART NO. 731303 FOR 1/4-28 NUTS AND FERRULES, SEE PAGE 33 OR VISIT CPCWORLDWIDE.COM


*Ferrules* 




The 1/8" flow PMC coupling covers a wide variety of general-purpose applications. Featuring the CPC thumb latch, the PMC is easier to use than ball-andsleeve designs. One-hand connection/disconnection and integral terminations make the PMC the choice for ease of use and manufacture.

This coupling is NuSeal® compatible (see page 27 for details).

FEATURES	BENEFITS	
Hose barb shroud	Extra protection for the 1/16" hose barb	
Integral terminations	Fewer leak points, shorter assembli faster installations	
Compatible	Mates with PMC12 and MC Series couplings (see pages 30 and 34)	
Clicks when connected	Assurance of a reliable connection	





# Specifications • • •

# PRESSURE:

Vacuum to 120 psi, 8.3 bar

#### TEMPERATURE:

-40°F to 180°F (-40°C to 82°C)

# MATERIALS:

Main components and valves: Acetal Thumb latch: Stainless steel

Valve spring: 316 stainless steel

External springs and pin: Stainless steel O-rings: Buna-N

# COLOR:

Natural white, others available


#### **TUBING SIZES:**

1/16" to 1/4" ID, 1.6mm to 6.4mm ID

WARNING: Pressure, temperature, chemicals, and operating environment can affect the performance of couplings. It is the customer's responsibility to test the suitability of CPC products in their own application conditions.



These graphs are intended to give you a general idea of the performance capabilities of each product line. The shaded area of each graph represents the operating range of the product family, i.e., upper and lower values are shown. Therefore, depending on the exact coupling configurations selected, you can reasonably expect values to fall within the shaded area.



# NU-SEAL® QUICK DISCONNECT COUPLING

# Nu-Seal® Specifications

#### PRESSURE:

Vacuum to 20 psi, 1.4 bar

#### **MATERIALS:**

Main components:

Low density polyethylene

#### COLOR:

Translucent

#### **TUBING SIZE:**

1/8" ID, 3.2mm ID

# Nu-Seal® quick disconnect coupling inserts feature a

patented integral seal and are USP Class VI compliant. These economic inserts are ideal for applications that require limited connect/disconnect cycles.

FEATURES	BENEFITS
Patented integral seal	Economical, allows for disposability and provides leak-free connections
Low density polyethylene	Sterilizable by EtO, E-beam & gamma
Chemically resistant	Works in many critical USP Class VI material applications
Compatible	Mates with MC, PMC, PMC12 coupling



**TUBING SIZE** 1/8" ID

PART NO. DSM2202 HEIGHT

LENGTH

**TERMINATION** IN-LINE HOSE BARB

METRIC EQ. 3.2 mm ID

.98

# Liquid Flow Rate Information for Couplings


The chart below shows the flow rate for CPC couplings. Each coupling was tested with water at  $70^{\circ}\text{F}$ (21°C). To determine flow rates for specific coupling configurations use the formula at the right.



- Q = Flow rate in gallons per minute Average coefficient across various
- flow rates (see chart)
- S = Specific gravity of liquid
- C<sub>v</sub> VALUES FOR 1/8" FLOW PMC COUPLINGS

INSERTS BODIES	PMC 2004	<i>PMCD</i> 2004	PMC 2006	PMCD 2006	PMC 2202	PMCD 2202	<i>PMC</i> 2204	PMCD 2204	PMC 2402	<i>PMCD</i> 2402	PMC 2404	PMCD 2304	PMC 2602	<i>PMCD</i> 2304	PMC 2104	PMCD 2304	PMC 2203	<i>PMCD</i> 2203	PMC 2201	PMCD 2201
PMC1002	.40	.18	.50	.19	.25	.16	.50	.19	.50	.20	.51	.19	.50	.50	.38	.24	.30	.17	.03	.03
PMCD1002	.27	.18	.31	.18	.24	.16	.28	.20	.26	.20	.29	.18	.26	.26	.27	.24	.25	.17	.03	.03
PMC1004	.40	.21	.50	.24	.26	.18	.50	.24	.50	.20	.51	.24	.50	.50	.38	.26	.30	.19	.03	.03
PMCD1004	.29	.19	.32	.23	.25	.17	.30	.23	.27	.21	.28	.23	.27	.28	.29	.24	.25	.18	.03	.03
PMC1204	.40	.18	.50	.18	.25	.16	.40	.18	.40	.16	.36	.18	.40	.40	.38	.21	.30	.17	.03	.03
PMCD1204	.21	.17	.22	.17	.20	.16	.22	.17	.21	.17	.20	.17	.21	.22	.21	.18	.21	.16	.03	.03
PMC1602	.23	.15	.28	.18	.19	.14	.27	.15	.27	.15	.28	.18	.27	.27	.23	.16	.20	.14	.03	.03
PMCD1602	.19	.15	.19	.15	.17	.14	.19	.15	.18	.15	.18	.15	.18	.19	.19	.15	.18	.14	.03	.03
PMC1604	.33	.23	.44	.24	.24	.18	.44	.23	.44	.20	.38	.24	.38	.44	.33	.26	.26	.19	.03	.03
PMCD1604	.23	.17	.26	.21	.22	.16	.26	.21	.26	.19	.25	.21	.21	.26	.23	.24	.22	.16	.03	.03
PMC1703	.25	.20	.30	.20	.20	.17	.30	.20	.30	.19	.28	.20	.28	.30	.25	.18	.21	.17	.03	.03
PMCD1703	.20	.17	.20	.17	.19	.15	.21	.17	.19	.17	.20	.17	.19	.20	.20	.16	.19	.16	.03	.03
PMC1701	.03	.03	.03	.03	.03	.03	.03	.03	.03	.03	.03	.03	.03	.03	.03	.03	.03	.03	.03	.02
PMCD1701	.03	.03	.03	.03	.03	.03	.03	.03	.03	.03	.03	.03	.03	.03	.03	.03	.03	.03	.02	.02

# PMC DIMENSIONS



# Panel Dimensions

	PANEL OPENING	PANEL THICKNESS MAX. – MIN.	PANEL NUT HEX	PANEL NUT THREAD
COUPLING BODIES	see drawing	.50 — .05	5/8	1/2-24UNS
COUPLING INSERTS	1/2	.30 – .05	5/8	1/2-24UNS




# **Coupling Bodies • ACETAL**

60	TERMINATION IN-LINE PIPE THREAD	<b>TUBING/THREAD SIZE</b> 1/8" NPT 1/8" BSPT 1/4" NPT 1/4" BSPT	METRIC EQ.	STRAIGHT THRU PMC1002 PMC1002BSPT PMC1004 PMC1004BSPT	SHUTOFF PMCD1002 PMCD1002BSPT PMCD1004 PMCD1004BSPT	.79 .79 .79 .79	B 1.00 1.00 1.10 1.10
OF	PANEL MOUNT FERRULELESS POLYTUBE FITTING, PTF†	5/32" OD, .10" ID 1/4" OD, .17" ID	4.0mm OD, 2.0mm ID 4.0mm OD, 2.5mm ID 6.4mm OD, 4.3mm ID	PMC12M42 PMC12025 PMC1204	PMCD12M42 PMCD12025 PMCD1204	.77 .77 .77	1.65 1.70 1.72
	PANEL MOUNT Hose Barb	1/16" ID 1/8" ID 3/16" ID 1/4" ID	1.6mm ID 3.2mm ID 4.8mm ID 6.4mm ID	PMC1601 PMC1602 PMC1603 PMC1604	PMCD1601 PMCD1602 PMCD1603 PMCD1604	.77 .77 .77 .77	1.40 1.65 1.85 1.85
	PANEL MOUNT FEMALE THREAD	10-32 UNF		PMC181032	PMCD181032	.77	1.25
POP .	IN-LINE FERRULELESS POLYTUBE FITTING, PTF†	5/32" OD, .10" ID	4.0mm OD, 2.0mm ID 4.0mm OD, 2.5mm ID 5.0mm OD, 3.0mm ID 6.4mm OD, 4.3mm ID	PMC13M42 PMC13025 PMC13M5 PMC1304	PMCD13M42 PMCD13025 PMCD13M5 PMCD1304	.81 .81 .81 .81	1.65 1.65 1.72 1.75
	IN-LINE HOSE BARB	1/16" ID 1/8" ID 3/16" ID 1/4" ID	1.6mm ID 3.2mm ID 4.8mm ID 6.4mm ID	PMC1701 PMC1702 PMC1703 PMC1704	PMCD1701 PMCD1702 PMCD1703 PMCD1704	.81 .81 .81 .81	1.41 1.45/1.65 1.65/1.85 1.65/1.85
80	IN-LINE FEMALE THREAD	10-32 UNF		PMC191032	PMCD191032	.81	1.25

All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters. †NOTE: CPC's Ferruleless PTF (polytube fitting) terminations do not require ferrules to achieve a secure connection and are therefore easier to use and reuse. PTF fittings are designed for semi-rigid tubing, i.e., polyethylene, nylon, polyurethane, etc.

# **Coupling Inserts • ACETAL**




All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters. Couplings are pictured with valves unless otherwise noted. †NOTE: CPC's Ferruleless PTF (polytube fitting) terminations do not require ferrules to achieve a secure connection and are therefore easier to use and reuse. PTF fittings are designed for semi-rigid tubing, i.e., polyethylene, nylon, polyurethane, etc.


# PMC12 SERIES CONNECTOR



The 1/8" flow polypropylene PMC12 offers many of the same configuration options as the PMC. The polypropylene material adds greater chemical resistance and is gamma sterilizable. The PMC12 also mates to small diameter rigid tubing. Available with a 1/4-28 flat bottom port and 1/4-28 UNF threads, these couplings eliminate the need to thread and re-thread common compression nuts each time a tubing connection is made.

FEATURES	BENEFITS
Polypropylene material	Chemically resistant and gamma- sterilizable
EPDM o-ring	Greater chemical resistance
CPC thumb latch	One-hand connection and disconnection
Integral terminations	Fewer leak points, shorter assemblies, faster installations





# Specifications • • •

#### PRESSURE:

Vacuum to 120 psi, 8.3 bar

#### TEMPERATURE:

32°F to 160°F (0°C to 71°C)

# **MATERIALS:**

Main components and valves: Polypropylene

Thumb latch: Stainless steel Valve spring: 316 stainless steel External spring and pin: Stainless steel

O-rings: EPDM

#### STERILIZATION:


**Gamma:** Up to 50 kGy irradiation


COLOR: Almond


### **TUBING SIZES:**

Microbore to 1/4" ID, Microbore to 6.4mm ID

WARNING: Pressure, temperature, chemicals, and operating environment can affect the performance of couplings. It is the customer's responsibility to test the suitability of CPC products in their own application conditions. Use the graph to the right as a guide.







# Liquid Flow Rate Information for Couplings

C<sub>v</sub> VALUES FOR 1/8" FLOW PMC12 COUPLINGS

The chart below shows the flow rate for CPC couplings. Each coupling was tested with water at  $70^{\circ}\text{F}$ (21°C). To determine flow rates for specific coupling configurations use the formula at the right.




- Q = Flow rate in gallons per minute
- Average coefficient across various flow rates (see chart)
- S = Specific gravity of liquid

INSERTS		PMCD12			I
BODIES	2004	2004	2006	2006	2

INSERTS BODIES	PMC12 2004	PMCD12 <b>2004</b>	PMC12 <b>2006</b>	PMCD12 <b>2006</b>	PMC12 <b>2202</b>	PMCD12 <b>2202</b>	PMC12 <b>2204</b>	PMCD12 <b>2204</b>	PMC12 <b>2402</b>	PMCD12 <b>2402</b>	PMC12 <b>2404</b>	PMCD12 <b>2304</b>	PMC12 <b>2602</b>	PMCD12 <b>2304</b>	PMC12 <b>2104</b>	PMCD12 <b>2304</b>	PMC12 <b>2203</b>	PMCD12 <b>2203</b>	PMC12 <b>2201</b>	PMCD12 <b>2201</b>
PMC100212	.40	.18	.50	.19	.25	.16	.50	.19	.50	.20	.51	.19	.50	.50	.38	.24	.30	.17	.03	.03
PMCD100212	.27	.18	.31	.18	.24	.16	.28	.20	.26	.20	.29	.18	.26	.26	.27	.24	.25	.17	.03	.03
PMC100412	.40	.21	.50	.24	.26	.18	.50	.24	.50	.20	.51	.24	.50	.50	.38	.26	.30	.19	.03	.03
PMCD100412	.29	.19	.32	.23	.25	.17	.30	.23	.27	.21	.28	.23	.27	.28	.29	.24	.25	.18	.03	.03
PMC120412	.40	.18	.50	.18	.25	.16	.40	.18	.40	.16	.36	.18	.40	.40	.38	.21	.30	.17	.03	.03
PMCD120412	.21	.17	.22	.17	.20	.16	.22	.17	.21	.17	.20	.17	.21	.22	.21	.18	.21	.16	.03	.03
PMC160212	.23	.15	.28	.18	.19	.14	.27	.15	.27	.15	.28	.18	.27	.27	.23	.16	.20	.14	.03	.03
PMCD160212	.19	.15	.19	.15	.17	.14	.19	.15	.18	.15	.18	.15	.18	.19	.19	.15	.18	.14	.03	.03
PMC160412	.33	.23	.44	.24	.24	.18	.44	.23	.44	.20	.38	.24	.38	.44	.33	.26	.26	.19	.03	.03
PMCD160412	.23	.17	.26	.21	.22	.16	.26	.21	.26	.19	.25	.21	.21	.26	.23	.24	.22	.16	.03	.03
PMC170312	.25	.20	.30	.20	.20	.17	.30	.20	.30	.19	.28	.20	.28	.30	.25	.18	.21	.17	.03	.03
PMCD170312	.20	.17	.20	.17	.19	.15	.21	.17	.19	.17	.20	.17	.19	.20	.20	.16	.19	.16	.03	.03
PMC170112	.03	.03	.03	.03	.03	.03	.03	.03	.03	.03	.03	.03	.03	.03	.03	.03	.03	.03	.03	.02
PMCD170112	.03	.03	.03	.03	.03	.03	.03	.03	.03	.03	.03	.03	.03	.03	.03	.03	.03	.03	.02	.02

# PMC12 DIMENSIONS




# Panel Dimensions

	PANEL OPENING	PANEL THICKNESS MAX. – MIN.	PANEL NUT HEX	PANEL NUT THREAD
COUPLING BODIES	see drawing	.50 — .05	5/8	1/2-24UNS
COUPLING INSERTS	see drawing	.30 — .06	5/8	1/2-24UNS



# **Coupling Bodies • POLYPROPYLENE**




# PMC12 1/4-28 Coupling Bodies • POLYPROPYLENE



All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters. †NOTE: CPC's Ferruleless PTF (polytube fitting) terminations do not require ferrules to achieve a secure connection and are therefore easier to use and reuse. PTF fittings are designed for semi-rigid tubing, i.e., polyethylene, nylon, polyurethane, etc.

# **Coupling Inserts • POLYPROPYLENE**



# PMC12 1/4-28 Coupling Inserts • POLYPROPYLENE



# Nuts

	1/16" and 1.8mm 1/8" and 3mm 1/16" and 1.8mm 1/16" and 3mm	DESCRIPTION F 1/4-28 Polypropylene Nut (natural) 1/4-28 Polypropylene Nut (natural) 1/4-28 Acetal Nut (black) 1/4-28 Acetal Nut (black)	
--	---	---	--

# **Ferrules**



All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters. Couplings are pictured with valves unless otherwise noted.


†NOTE: CPC's Ferruleless PTF (polytube fitting) terminations do not require ferrules to achieve a secure connection and are therefore easier to use and reuse. PTF fittings are designed for semi-rigid tubing, i.e., polyethylene, nylon, polyurethane, etc.




# CPC's MC Series chrome-plated brass couplings are

built tough and made to last in the most demanding applications. Ideal for use with higher temperatures or pressures, MC Series couplings feature one-hand operation for swift and easy connects and disconnects. These couplings offer the flexibility of multiple configurations and terminations and mate with both PMC acetal and PMC12 polypropylene couplings.

FEATURES	BENEFITS
Brass material	Durable and able to withstand higher pressure and temperature
High temperature capability	Versions rated to 400°F (204°C)
CPC thumb latch	One-hand connection and disconnection
Compatible	MC mates with PMC and PMC12 Series couplings





# Specifications • • •

#### PRESSURE:

Vacuum to 250 psi, 17.3 bar

#### TEMPERATURE:

-40°F to 180°F (-40°C to 82°C) (High temperature versions available with ratings to 400°F)

#### **MATERIALS:**

Main components: Chrome-plated brass

Thumb latch: Stainless steel

Valves: Acetal

Valve spring: 316 stainless steel

External springs and pin: Stainless steel

**O-rings:** Buna-N

FINISH: Chrome

#### **TUBING SIZES:**

1/8" to 1/4" ID, 3.2mm to 6.4mm ID

WARNING: Pressure, temperature, chemicals, and operating environment can affect the performance of couplings. It is the customer's responsibility to test the suitability of CPC products in their own application conditions. Use the graph to the right as a guide.

Note:
High temperature versions available with ratings to 400 F. Call customer service for more information.

Also available in NSF listed versions, please visit our website for part number information.

These graphs are intended to give you a general idea of the performance capabilities of each product line. The shaded area of each graph represents the operating range of the product family, i.e., upper and lower values are shown. Therefore, depending on the exact coupling configurations selected, you can reasonably expect values to fall within the shaded area.



# Liquid Flow Rate Information for Couplings

C<sub>v</sub> VALUES FOR 1/8" FLOW MC COUPLINGS


The chart below shows the flow rate for CPC couplings. Each coupling was tested with water at 70°F (21°C). To determine flow rates for specific coupling configurations use the formula at the right.



- Q = Flow rate in gallons per minute
- C<sub>v</sub> = Average coefficient across various flow rates (see chart)
- $\Delta extbf{P}= ext{ Pressure drop across coupling (psi)}$
- S = Specific gravity of liquid

BODIES INSERTS	MC12 2004	MCD12 2004	MC12 2006	MCD12 2006	MC12 2202	MCD12 2202	MC12 2204	MCD12 2204	MC12 2402	MCD12 2402	MC12 2404	MCD12 2304	MC12 2602	MCD12 2304	MC12 2104	MCD12 2304	MC12 2203	MCD12 2203
MC100212	.40	.18	.50	.19	.25	.16	.50	.19	.50	.20	.51	.19	.50	.50	.38	.24	.30	.17
MCD100212	.27	.18	.31	.18	.24	.16	.28	.20	.26	.20	.29	.18	.26	.26	.27	.24	.25	.17
MC100412	.40	.21	.50	.24	.26	.18	.50	.24	.50	.20	.51	.24	.50	.50	.38	.26	.30	.19
MCD100412	.29	.19	.32	.23	.25	.17	.30	.23	.27	.21	.28	.23	.27	.28	.29	.24	.25	.183
MC120412	.40	.18	.50	.18	.25	.16	.40	.18	.40	.16	.36	.18	.40	.40	.38	.21	.30	.173
MCD120412	.21	.17	.22	.17	.20	.16	.22	.17	.21	.17	.20	.17	.21	.22	.21	.18	.21	.16
MC160212	.23	.15	.28	.18	.19	.14	.27	.15	.27	.15	.28	.18	.27	.27	.23	.16	.20	.14
MCD160212	.19	.15	.19	.15	.17	.14	.19	.15	.18	.15	.18	.15	.18	.19	.19	.15	.18	.14
MC160412	.33	.23	.44	.24	.24	.18	.44	.23	.44	.20	.38	.24	.38	.44	.33	.26	.26	.19
MCD160412	.23	.17	.26	.21	.22	.16	.26	.21	.26	.19	.25	.21	.21	.26	.23	.24	.22	.16
MC170312	.25	.20	.30	.20	.20	.17	.30	.20	.30	.19	.28	.20	.28	.30	.25	.18	.21	.17
MCD170312	.20	.17	.20	.17	.19	.15	.21	.17	.19	.17	.20	.17	.19	.20	.20	.16	.19	.16

# MC DIMENSIONS



# **Panel Dimensions**

	PANEL OPENING	PANEL THICKNESS MAXMIN.	PANEL NUT HEX	PANEL NUT THREAD
COUPLING BODIES (except male thread)	see drawing	.50 — .05	5/8	1/2-24UNS
MALE THREAD	see drawing	.21 – .06	5/8	1/2-24UNS
COUPLING INSERTS	see drawing	.30 – .06	5/8	1/2-24UNS




# **Coupling Bodies • CHROME-PLATED BRASS**

1	TERMINATION IN-LINE PIPE THREAD	TUBING/THREAD SIZE 1/8" NPT 1/8" BSPT 1/4" NPT 1/4" BSPT	E METRIC EQ.	STRAIGHT THRU MC1002 MC1002BSPT MC1004 MC1004BSPT	SHUTOFF MCD1002 MCD1002BSPT MCD1004 MCD1004BSPT	.76 .76 .76 .76	1.00 1.00 1.10 1.10
OT THE	PANEL MOUNT FERRULELESS POLYTUBE FITTING, PTF†	5/32" OD, .10" ID 1/4" OD, .17" ID	4.0mm OD, 2.5mm ID 6.4mm OD, 4.3mm ID	MC12025 MC1204	MCD12025 MCD1204	.76 .76	1.67 1.77
OF THE	PANEL MOUNT Hose Barb	1/8" ID 3/16" ID 1/4" ID	3.2mm ID 4.8mm ID 6.4mm ID	MC1602 MC1603 MC1604	MCD1602 MCD1603 MCD1604	.76 .76 .76	1.67 1.89 1.89
OTTO	PANEL MOUNT Female Thread	10-32 UNF		MC181032	MCD181032	.76	1.22
	PANEL MOUNT Male Thread	1/8" NPT		MC1502	MCD1502	.76	1.20
	IN-LINE FERRULELESS POLYTUBE FITTING, PTF†		6.4mm OD, 4.3mm ID	MC1304	MCD1304	.76	1.77
The same	IN-LINE HOSE BARB	1/8" ID 3/16" ID 1/4" ID	3.2mm ID 4.8mm ID 6.4mm ID	MC1702 MC1703 MC1704	MCD1702 MCD1703 MCD1704	.76 .76 .76	1.67 1.89 1.89

All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters. †NOTE: CPC's Ferruleless PTF (polytube fitting) terminations do not require ferrules to achieve a secure connection and are therefore easier to use and reuse. PTF fittings are designed for semi-rigid tubing, i.e., polyethylene, nylon, polyurethane, etc. NOTE: Elbow configurations are available. Contact CPC for more information.

# **Coupling Inserts • CHROME-PLATED BRASS**



All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters. Couplings are pictured with valves unless otherwise noted.

tNote: CPC's Ferruleless PTF (polytube fitting) terminations do not require ferrules to achieve a secure connection and are therefore easier to use and reuse. PTF fittings are designed for semi-rigid tubing, i.e., polyethylene, nylon, polyurethane, etc. NOTE: Elbow configurations are available. Contact CPC for more information.



# Specifications • • •

#### PRESSURE:

Vacuum to 45 psi, 3.1 bar

#### TEMPERATURE:

32°F to 120°F (0°C to 49°C)

#### MATERIALS:

Main components: Glass-filled polypropylene

Valve spring: 316 stainless steel

Seal material:

EPDM (FKM/Simriz® FFKM optional)

**COLOR:** Gray with black accent standard

#### **TUBING SIZES:**


1/8" and 1/4" ID, 3.2mm and 6.4mm ID

**SPILLAGE:** <0.025 cc/disconnect


INCLUSION: 0.009 cc/connect


WARNING: Pressure, temperature, chemicals, and operating environment can affect the performance of couplings. It is the customer's responsibility to test the suitability of CPC's products in their own application conditions.

**NS212 Series couplings** were built on the company's proven platform of non-spill valves designed to provide fast, safe and leak-free fluid line connections. The NS212 is an easy, twist-to-connect coupling with an integrated locking mechanism and double-sided non-spill shutoff valves. It is also chemically resistant and ideal for a diverse range of applications including printing and ink management, analytical instrumentation, electronic cooling and chemical handling. NS212 couplings provide high-flow fluid transfer in a small footprint.




▼ These graphs are intended to give you a general idea of the performance capabilities of each product line. The shaded area of each graph represents the operating range of the product family, i.e., upper and lower values are shown. Therefore, depending on the exact coupling configurations selected, you can reasonably expect values to fall within the shaded area.





# NS212 DIMENSIONS



# **Panel Dimensions**

	PANEL OPENING	PANEL THICKNESS MAXMIN.	+
COUPLING BODIES	see drawing	.3005	
	1		Ø .688" +.010 000



# **Coupling Bodies** • POLYPROPYLENE



# **Coupling Inserts • POLYPROPYLENE**



TERMINATION
IN-LINE
PIPE THREAD

IN-LINE
FERRULESS
POLYTUBE FITTING,
PTF†

IN-LINE
HOSE BARB

TUBING SIZE 1/8" NPT	METRIC EQ.	SHUTOFF NS2D240212	.64	B 1.27
1/4" OD, .17 ID	6.4mm OD, 4.3mm ID	NS2D200412	.64	1.39
1/8" ID	3.2mm ID	NS2D220212		1.32
1/4" ID	6.4mm ID	NS2D220412	.64	1.52

The NS212 does not mate with NS2. Contact CPC for additional compatibility information.

# Accessories

**DESCRIPTION PANEL MOUNT GASKET REPLACEMENT:** For sealing panel mount bodies listed above

MATERIAL PART NO.
Buna-N 1830300


All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters. PTF fittings are designed for semi-rigid tubing, i.e., polyethylene, nylon, polyurethane, etc. †NOTE: CPC's Ferruleless PTF (polytube fitting) terminations do not require ferrules to achieve a secure connection and are therefore easier to use and reuse.


# PLC SERIES CONNECTOR



The 1/4" flow PLC Series is proven worldwide in thousands of applications and offers the widest selection of sizes and configurations. PLC couplings are injection molded from acetal thermoplastic and are resistant to most mild chemical solutions. One-hand connection/disconnection, plus integral terminations make the PLC Series the choice for ease of use and manufacture.

FEATURES	BENEFITS
CPC thumb latch	One-hand connection and disconnection
Integral terminations	Fewer leak points, shorter assemblies, faster installations
Clicks when connected	Assurance of a reliable connection
Compatible	Mates with LC and PLC12 Series couplings





# Specifications • • •

#### PRESSURE:

Vacuum to 120 psi, 8.3 bar

#### TEMPERATURE:

-40°F to 180° F (-40°C to 82°C)

#### MATERIALS:

**O-rings:** Buna-N

Main components and valves: Acetal Thumb latch: Stainless steel Valve spring: 316 stainless steel External springs and pin: Stainless steel

**COLOR:** Natural white, others available

#### **TUBING SIZES:**

1/4" to 3/8" ID. 6.4mm to 9.5mm ID

#### JG® TUBING SPECIFICATIONS:

Tube tolerances: 1/4" OD, +0.001/-0.004 | 3/8" OD,


Tube Types: Plastic tube: Polyethylene, nylon, polyurethane. For soft or thin wall tubing with JG® terminations, tube supports are recommended.

Metal tube: Brass, copper and mild steel

Also available in NSF listed versions, please visit our website for part number information.



> These graphs are intended to give you a general idea of the performance capabilities of each product line. The shaded area of each graph represents the operating range of the product family, i.e., upper and lower values are shown. Therefore, depending on the exact coupling configurations selected, you can reasonably expect values to fall within the shaded area.



# **Decoding Coupling Options**

Often times, it can be an overwhelming task to decipher which coupling best fits your application needs. The chart below helps to decode popular CPC couplings.

1/4" FLOW OPTIONS	PRIMARY MATERIAL	O-RING	APPLICATION BENEFITS		
PLC	Acetal	Buna-N	General purpose plastic connector with metal latch		
PLC12	Polypropylene	EPDM	Expanded chemical compatibility		
LC	Chrome-Plated Brass	Buna-N	Higher operating temperature and pressure		
APC	Acetal	Buna-N	Plastic thumb latch		
BREAKAWAY®	Acetal	Buna-N	Breakaway latch		
NS4	Glass-filled polypropylene	EPDM	Drybreak, non-spill connect and disconnect, even under pressure		

For further information please contact CPC Customer Service at 1-800-444-2474 or 651-645-0091.

# Liquid Flow Rate Information for Couplings

C<sub>v</sub> VALUES FOR 1/4" FLOW PLC COUPLINGS

The chart below shows the flow rate for CPC couplings. Each coupling was tested with water at 70°F (21°C). To determine flow rates for specific coupling configurations use the formula at the right.



- Q = Flow rate in gallons per minute
- Cv = Average coefficient across various flow rates (see chart)
- $\Delta P$  = Pressure drop across coupling (psi)
- S = Specific gravity of liquid


BODIES INSERTS	PLC 20004	PLCD 20004	PLC 20006	PLCD 20006	PLC 22004	PLCD 22004	PLC 22006	PLCD 22006	PLC 24004	PLCD 24004	PLC 24006	PLCD 24006	PLC 26004
PLC10004	.40	.36	1.05	.58	.83	.56	1.40	.82	1.40	.75	1.40	.77	.83
PLCD10004	.36	.31	.73	.48	.66	.41	.82	.50	.80	.45	.77	.45	.81
PLC10006	.40	.36	1.05	.60	.83	.56	1.40	.81	1.40	.76	1.40	.76	.83
PLCD10006	.37	.31	.81	.47	.70	.43	1.02	.51	.98	.46	.99	.48	.98
PLC12006	.38	.36	.84	.63	.74	.56	1.14	.75	1.14	.70	1.14	.72	.74
PLCD12006	.38	.33	.78	.49	.68	.44	.84	.49	.81	.43	.82	.44	.81
PLC16004	.38	.37	.87	.54	.95	.51	1.00	.70	.95	.64	1.00	.66	.95
PLCD16004	.37	.31	.61	.44	.57	.41	.78	.44	.78	.43	.75	.46	.78
PLC16006	.38	.37	1.00	.57	.95	.53	1.40	.80	1.40	.71	1.40	.73	1.40
PLCD16006	.38	.32	.71	.49	.63	.42	.89	.51	.96	.45	.92	.49	.97

# PLC DIMENSIONS

# A = Height/Diameter B = Total Length

# **Panel Dimensions**

	PANEL OPENING	PANEL THICKNESS MAXMIN.	PANEL NUT HEX	PANEL NUT THREAD
COUPLING BODIES	see drawing	.50 — .05	13/16	11/16-24UNEF
COUPLING INSERTS	see drawing	.30 — .05	13/16	11/16-24UNEF



# **Coupling Bodies • ACETAL**

D = Elbow Radial Length

6	TERMINATION IN-LINE PIPE THREAD	TUBING/THREAD SIZE 1/4" NPT 1/4" BSPT 3/8" NPT 3/8" BSPT	METRIC EQ.	STRAIGHT THRU PLC10004 PLC10004BSPT PLC10006 PLC10006BSPT	SHUTOFF PLCD10004 PLCD10004BSPT PLCD10006 PLCD10006BSPT	.90 .90 .90	B 1.17 1.15 1.15 1.15
(T)	PANEL MOUNT FERRULELESS POLYTUBE FITTING, PTF†	1/4" OD, .17" ID 3/8" OD, .25" ID	6.4mm OD, 4.3mm ID 8.0mm OD, 6.0mm ID 9.5mm OD, 6.4mm ID 10.0mm OD, 8.0mm ID	PLC12004 PLC120M8 PLC12006 PLC120M10	PLCD12004 PLCD120M8 PLCD12006 PLCD120M10	.88 .88 .88	1.87 1.95 1.95 1.95
	PANEL MOUNT Hose Barb	1/4" ID 5/16" ID 3/8" ID	6.4mm ID 7.9mm ID 9.5mm ID	PLC16004 PLC16005 PLC16006	PLCD16004 PLCD16005 PLCD16006	.88 .88 .88	1.95 1.95 1.95
(2)	PANEL MOUNT JG® PUSH-TO- CONNECT	1/4" OD	6.4mm OD	PLC11004	PLCD11004	.88	1.95
	IN-LINE FERRULELESS POLYTUBE FITTING, PTF†	1/4" OD, .17" ID 3/8" OD, .25" ID	6.4mm OD, 4.3mm ID 8.0mm OD, 6.0mm ID 9.5mm OD, 6.4mm ID 10.0mm OD, 8.0mm ID	PLC13004 PLC130M8 PLC13006 PLC130M10	PLCD13004 PLCD130M8 PLCD13006 PLCD130M10	.94 .94 .94 .94	1.82 1.95 1.95 1.95
60	IN-LINE HOSE BARB	1/4" ID 5/16" ID 3/8" ID	6.4mm ID 7.9mm ID 9.5mm ID	PLC17004 PLC17005 PLC17006	PLCD17004 PLCD17005 PLCD17006	.94 .94 .94	1.80 /1.95 1.95 1.80 /1.95
	IN-LINE JG® PUSH-TO- CONNECT	1/4" OD 3/8" OD (See page 162 for JG ac	6.4mm OD 9.5mm OD ccessories)	PLC14004 PLC14006	PLCD14004 PLCD14006	.94 .94	1.95 2.09

All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters. †NOTE: CPC's Ferruleless Polytube Fitting terminations do not require ferrules to achieve a secure connection, which makes them easier to use and reuse. PTF fittings are designed for semi-rigid tubing, i.e., polyethylene, nylon, polyurethane, etc. NOTE: JG is a registered trademark of John Guest USA, Inc.

# **Coupling Inserts • ACETAL**



All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters. Couplings are pictured with valves unless otherwise noted.


†NOTE: CPC's Ferruleless Polytube Fitting terminations do not require ferrules to achieve a secure connection, which makes them easier to use and reuse. PTF fittings are designed for semi-rigid tubing, i.e., polyethylene, nylon, polyurethane, etc. NOTE: JG is a registered trademark of John Guest USA, Inc.


# PLC12 SERIES CONNECTOR



**The 1/4" flow polypropylene PLC12** offers many of the same configuration options as the PLC. The polypropylene material adds greater chemical resistance for more demanding applications. PLC12 couplings are also gamma sterilizable. PLC12 coulings are also gamma sterilizable and available with optional RFID (Radio Frequency Identification) capability (see page 96).

FEATURES	BENEFITS					
Polypropylene material	Chemically resistant and gamma sterilizable					
EPDM o-rings	Greater chemical resistance					
CPC thumb latch	One-hand connection and disconnection					
Integral terminations	Fewer leak points, shorter assemblies, faster installations					





# Specifications • • •

#### PRESSURE:

Vacuum to 120 psi, 8.3 bar

#### TEMPERATURE:

32°F to 160°F (0°C to 71°C)

#### MATERIALS:

Main components and valves:

Polypropylene

Thumb latch: Stainless steel

Valve spring: 316 stainless steel

External springs and pin: Stainless steel

**0-rings:** EPDM

#### **STERILIZATION:**

Gamma: Up to 50 kGy irradiation

COLOR: Almond


#### **TUBING SIZES:**

1/4" to 3/8" ID, 6.4mm to 9.5mm ID

WARNING: Pressure, temperature, chemicals, and operating environment can affect the performance of couplings. It is the customer's responsibility to test the suitability of CPC's products in their own application conditions. Use the graph at the right as a guide.



These graphs are intended to give you a general idea of the performance capabilities of each product line. The shaded area of each graph represents the operating range of the product family, i.e., upper and lower values are shown. Therefore, depending on the exact coupling configurations selected, you can reasonably expect values to fall within the shaded area.



# Liquid Flow Rate Information for Couplings

The chart below shows the flow rate for CPC couplings. Each coupling was tested with water at 70°F (21°C). To determine flow rates for specific coupling configurations use the formula at the right.

# $Q=C_{V}\sqrt{\frac{\Delta P}{S}}$

Q = Flow rate in gallons per minute

C<sub>v</sub> = Average coefficient across various flow rates (see chart)


 $\Delta P$  = Pressure drop across coupling (psi) S = Specific gravity of liquid

C<sub>v</sub> Values for 1/4" flow Plc12 couplings

BODIES INSERTS	PLC12 20004	PLCD12 20004	PLC12 20006	PLCD12 20006	PLC12 22004	PLCD12 22004	PLC12 22006	PLCD12 22006	PLC12 24004	PLCD12 24004	PLC12 24006	PLCD12 24006	PLC12 26004
PLC1000412	.40	.36	1.05	.58	.83	.56	1.40	.82	1.40	.75	1.40	.77	.83
PLCD1000412	.36	.31	.73	.48	.66	.41	.82	.50	.80	.45	.77	.45	.81
PLC1000612	.40	.36	1.05	.60	.83	.56	1.40	.81	1.40	.76	1.40	.76	.83
PLCD1000612	.37	.31	.81	.47	.70	.43	1.02	.51	.98	.46	.99	.48	.98
PLC1200612	.38	.36	.84	.63	.74	.56	1.14	.75	1.14	.70	1.14	.72	.74
PLCD1200612	.38	.33	.78	.49	.68	.44	.84	.49	.81	.43	.82	.44	.81
PLC1600412	.38	.37	.87	.54	.95	.51	1.00	.70	.95	.64	1.00	.66	.95
PLCD1600412	.37	.31	.61	.44	.57	.41	.78	.44	.78	.43	.75	.46	.78
PLC1600612	.38	.37	1.00	.57	.95	.53	1.40	.80	1.40	.71	1.40	.73	1.40
PLCD1600612	.38	.32	.71	.49	.63	.42	.89	.51	.96	.45	.92	.49	.97


											Т																								$\neg$						
M	N7	ES	c.			_					+						Н									+				Н									+		
/40	וע	L	,								+					+	Н													Н											
-	+		_			-	+	+	_	-	+	-	+	-	_	+			+			-	_		-	+			+		-	-	+		-	+					
+	+	+	+			+	+	+		+	+	+		-	-	+	Н		+			+	+		+	+		+		Н		+	+		+	+			+	+	
-	+	+	+			-	+	+	-	-	+	-	-	-	-	+	Н		+			+	+		+	+	-	+	+	Н	-	-	+		+	+		-	+	+	
	-	-						-			+		-	-		+	$\vdash$								-	-	-			Н			+		-	-	-			-	
-	-	-	+			-	+	+	-	-	+	-	-	-	-	+	Н		+			-	+		-	+	-	-	+	Н		-	+		+	+		-	-	-	
-	_	-				-	-	-		_	+	_	-	-		+	Н		+			-			-	_	-	-	-	Н		_	+		-	-		-	_	_	_
_	_	_				_	_	_		_	4			_		-			_						_	_				Ш			-		4	_			_	_	
	_										4						Ш		_							_				Ш					_	_					
											Т			T			П																								
	$\top$													T			П													П	$\top$								$\top$		
	$\top$	$\top$					$\top$	$\top$		$\exists$	$\top$		$\Box$				П					$\top$			$\top$					П	$\top$				$\top$	$\top$		$\exists$	$\top$	7	
$\dashv$	$\top$	$\top$				$\dashv$	$\top$				$\top$						П					$\top$			$\top$					П	$\top$				$\top$	$\top$		$\exists$	$\top$		
	_																													Н											
_	+					_					+																														
	+							+			+		+			+			+							+				Н			+						+	+	
	-	+					+	+		_	+			-		+			+						+	+		+		Н			+		+	+			_		_
-	+		-			-		+	-	-	+	-	+	-		+	Н		+				-		+	+	-	-		Н			+		+	+			-	+	
-	+	+	+			-	+	+	-	-	+	-	-	-	-	+	Н		+			+	+		+	+	-	+	+	Н	-	-	+		+	+		-	+	+	
-	+	+	+			-	+	+	-	-	+	-	-	-	-	+	Н		+			+	+		+	+	-	+	+	Н	-	-	+		+	+		-	+	+	
-	-	-	-			-	-	+		-	+	-	-	-		+	Н		+			-	-		-	+	-	-	+	Н		-	+		+	+	-		-	-	
-	-	-	+			-	+	+	-	-	+	-	-	-	-	+	Н		+			-	+		-	+	-	-	+	Н		-	+		+	+		-	-	-	
-	-	-	+			-	+	+	-	-	+	-	-	-	-	+	Н		+			-	+		-	+	-	-	+	Н		-	+		+	+		-	-	-	
-	-	-	_			-	-	-		-	+	_	-	-		+	Н		-			-	_		-	-	-	-	-	Н		_	+		-	-	-		_	-	
-	_	-				-	-	+		_	+	_	-	-		+	Н		+			-			-	_	-	-	-	Н		_	+		-	-		-	_	_	_
_	_	_				_	_	_	-	_	4		-	4	_	-	Ш		-						-	_	$\perp$						+		4	_		_	_	4	
4	4	_			Ш	4	4	$\perp$		_	4			4	_	-	Ш		-			_				-	4	4		$\square$		$\perp$	1		4	4		_	_	4	
4						4					_			4			Ш					4				1				Ш			1		4	_					
						_					_			4												X									4	_					
											1																														
											1																			H											
																														Z											
																		7		V		1													7						
	đ	*				J	X			4		7		1	W	٨.	M		V	M	74	XI.	ΠV	$\square$	V	X		7	X	И	7				7						
	×					ZÍ.									All	Ħ	V	X	A			V	A		XI.	11	11										H				
X	7	1				1	Ť	X		X		X			HÌ	1	1X	1	HÌ	1/1	M	A	17		¥	A		X	1			1						7	1	7	
	4		4	×		H	×			đ	×			1	P	1	H	X 17	₩	X	///	///	#	X	Alv	4	1/1	X	1	Н	X	1			7				7		
ä	*					1			K		7		H	7		#	H	$\forall$	+	H	H	¥	7	H	7	X	44	X		H	4	4	+	d	ď				7	1	
			+			=						+	1/1	X	X	4	H	$\sqrt{}$	/ /	$\forall$	X	Ak	+	Н	$\forall$	$\times$									+				-		/

# PLC12 DIMENSIONS




# Panel Dimensions

	PANEL OPENING	PANEL THICKNESS MAXMIN.	PANEL NUT HEX	PANEL NUT THREAD
COUPLING BODIES	see drawing	.30 – .02	13/16	11/16-24UNEF
COUPLING INSERTS	see drawing	.30 — .02	13/16	11/16-24UNEF



# **Coupling Bodies • POLYPROPYLENE**




# **PLC12 Coupling Readers • POLYPROPYLENE**



All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters. †NOTE: CPC's Ferruleless Polytube Fitting terminations do not require ferrules to achieve a secure connection, which makes them easier to use and reuse. PTF fittings are designed for semi-rigid tubing, i.e., polyethylene, nylon, polyurethane, etc. NOTE: JG is a registered trademark of John Guest USA, Inc.

# **Coupling Inserts • POLYPROPYLENE**



All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters. Couplings are pictured with valves unless otherwise noted.


†NOTE: CPC's Ferruleless Polytube Fitting terminations do not require ferrules to achieve a secure connection, which makes them easier to use and reuse. PTF fittings are designed for semi-rigid tubing, i.e., polyethylene, nylon, polyurethane, etc. NOTE: JG is a registered trademark of John Guest USA, Inc.


# LC SERIES CONNECTOR



**CPC's LC Series chrome-plated brass** couplings are built tough and made to last in the most demanding applications. Ideal for use with higher temperature or pressure, the LC Series features a one-hand operation for swift and easy connects and disconnects. These couplings offer the flexibility of multiple configurations and terminations to mate with both the PLC acetal and PLC12 polypropylene couplings.

FEATURES	BENEFITS					
Brass material	Durable construction withstands higher pressure and temperature					
Chrome plating	Attractive appearance					
High temperature capability	Versions rated to 400°F (204°C)					
CPC thumb latch One-hand connection and disco						
Compatible	LC mates with PLC Series couplings					





# Specifications • • •

#### PRESSURE:

Vacuum to 250 psi, 17.3 bar

#### TEMPERATURE:

-40°F to 180°F (-40°C to 82°C) (High temperature versions available with ratings to 400°F)

### **MATERIALS:**

Main components: Chrome-plated brass

Thumb latch: Stainless steel

Valves: Acetal

Valve springs: 316 stainless steel

External springs and pin: Stainless steel

**0-rings**: Buna-N

FINISH: Chrome

#### **TUBING SIZES:**

1/4" to 3/8" ID. 6.4mm to 9.5mm ID

WARNING: Pressure, temperature, chemicals, and operating environment can affect the performance of couplings. It is the customer's responsibility to test the suitability of CPC's products in their own application conditions.

# Note:

High temperature versions available with ratings to 400 F. Call customer service for more information.

Also available in NSF listed versions for foodbased applications, please visit our website for part number information.



These graphs are intended to give you a general idea of the performance capabilities of each product line. The shaded area of each graph represents the operating range of the product family, i.e., upper and lower values are shown. Therefore, depending on the exact coupling configurations selected, you can reasonably expect values to fall within the shaded area.

# DID YOU KNOW?

**O-RING SELECTION IS A KEY DECISION** in determining which connector will perform best in your specific application. Understanding the material characteristics and how they can be affected by both the media being transferred and the environment in which the connector is being used is important.

**BUNA-N** is the most common o-ring seal due to its solvent, oil and water resistance. Temperature range is -40°F to 250°F.

**FKM** is best suited for applications where chemical and strong acid resistance and/or high temperatures are a requirement. The temperature range is -15°F to 400°F.

**EPDM** has excellent resistance to polar solvents. This o-ring is not compatible with petroleum. The temperature range is -70°F to 300°F.

**SILICONE** seals have good temperature resistance. The temperature range is -70°F to 400°F and with special compounds can reach 175°F to 450°F. Silicone can also be supplied with Class VI requirements for medical/bioprocessing applications.

For further information about which o-ring is best suited for your application, please contact CPC Customer Service at 1-800-444-2474 or 651-645-0091.



# Liquid Flow Rate Information for Couplings

C<sub>v</sub> VALUES FOR 1/4" FLOW LC COUPLINGS


The chart below shows the flow rate for CPC couplings. Each coupling was tested with water at 70°F (21°C). To determine flow rates for specific coupling configurations use the formula at the right.

 $Q=C_{V}\sqrt{\frac{\Delta P}{S}}$ 

- Q = Flow rate in gallons per minute
- C<sub>v</sub> = Average coefficient across various flow rates (see chart)
- $\Delta P$  = Pressure drop across coupling (psi)
- S = Specific gravity of liquid


BODIES INSERTS	LC 20004	LCD 20004	LC 20006	LCD 20006	LC 22004	LCD 22004	LC 22006	LCD 22006	LC 24004	LCD 24004	LC 24006	LCD 24006	LC 26004
LC10004	.40	.36	1.05	.58	.83	.56	1.40	.82	1.40	.75	1.40	.77	.83
LCD10004	.36	.31	.73	.48	.66	.41	.82	.50	.80	.45	.77	.45	.81
LC10006	.40	.36	1.05	.60	.83	.56	1.40	.81	1.40	.76	1.40	.76	.83
LCD10006	.37	.31	.81	.47	.70	.43	1.02	.51	.98	.46	.99	.48	.98
LC12006	.38	.36	.84	.63	.74	.56	1.14	.75	1.14	.70	1.14	.72	.74
LCD12006	.38	.33	.78	.49	.68	.44	.84	.49	.81	.43	.82	.44	.81
LC16004	.38	.37	.87	.54	.95	.51	1.00	.70	.95	.64	1.00	.66	.95
LCD16004	.37	.31	.61	.44	.57	.41	.78	.44	.78	.43	.75	.46	.78
LC16006	.38	.37	1.00	.57	.95	.53	1.40	.80	1.40	.71	1.40	.73	1.40
LCD16006	.38	.32	.71	.49	.63	.42	.89	.51	.96	.45	.92	.49	.97

# LC DIMENSIONS



# **Product Dimensions**

	PANEL OPENING	PANEL THICKNESS MAXMIN.	PANEL NUT HEX	PANEL NUT THREAD
<b>COUPLING BODIES</b>	see drawing	.50 — .05	13/16	11/16-24UNEF
COUPLING INSERTS	see drawing	.090– .300	13/16	11/16-24UNEF




# **Coupling Bodies • CHROME-PLATED BRASS**

0-1	TERMINATION IN-LINE PIPE THREAD	TUBING/THREAD SIZE 1/4" NPT 1/4" BSPT 3/8" NPT 3/8" BSPT	METRIC EQ.	STRAIGHT THRU LC10004 LC10004BSPT LC10006 LC10006BSPT	SHUTOFF LCD10004 LCD10004BSPT LCD10006 LCD10006BSPT	.88 .88 .88	1.15 1.15 1.15 1.15 1.15
D-BR.	PANEL MOUNT FERRULELESS POLYTUBE FITTING, PTF†	1/4" OD, .17" ID 3/8" OD, .25" ID	6.4mm OD, 4.3mm ID 9.5mm OD, 6.4mm ID	LC12004 LC12006	LCD12004 LCD12006	.88 .88	1.87 2.00
S. Barre	PANEL MOUNT HOSE BARB	1/4" ID 5/16" ID 3/8" ID	6.4mm ID 7.9mm ID 9.5mm ID	LC16004 LC16005 LC16006	LCD16004 LCD16005 LCD16006	.88 .88 .88	2.00 2.00 2.00
	PANEL MOUNT Male Thread	1/4" NPT		LC15004	LCD15004	.88	1.50
	IN-LINE FERRULELESS POLYTUBE FITTING, PTF†	1/4" OD, .17" ID 3/8" OD, .25" ID	6.4mm OD, 4.3mm ID 9.5mm OD, 6.4mm ID	LC13004 LC13006	LCD13004 LCD13006	.88 .88	1.87 2.00
S. A. S.	IN-LINE HOSE BARB	1/4" ID 5/16" ID 3/8" ID	6.4mm ID 7.9mm ID 9.5mm ID	LC17004 LC17005 LC17006	LCD17004 LCD17005 LCD17006	.88 .88 .88	2.00 2.00 2.00

All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters. †NOTE: CPC's Ferruleless Polytube Fitting terminations do not require ferrules to achieve a secure connection, which makes them easier to use and reuse. PTF fittings are designed for semi-rigid tubing, i.e., polyethylene, nylon, polyurethane, etc.

# **Coupling Inserts • CHROME-PLATED BRASS**




All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters. Couplings are pictured with valves unless otherwise noted.


†NOTE: CPC's Ferruleless Polytube Fitting terminations do not require ferrules to achieve a secure connection, which makes them easier to use and reuse. PTF fittings are designed for semi-rigid tubing, i.e., polyethylene, nylon, polyurethane, etc.

# APC SERIES CONNECTOR

The 1/4" flow APC Series couplings feature a plastic thumb latch, fewer moving parts and a smooth contoured design to deliver ease-of-use and excellent flow in a compact size. The APC Series includes couplings with valves and thumb latch shrouds. A more rugged design on some models adds durability to increase the APC's application range.

FEATURES	BENEFITS
Plastic thumb latch	No exposed metal parts
Shrouded thumb latch	Protects against accidental disconnects
Integral terminations	Fewer leak points, shorter assemblies, faster installations
Compatibility	Mates with most EFC couplings





These graphs are intended to give you a general idea of the performance capabilities of each product line. The shaded area of each graph represents the operating range of the product family, i.e., upper and lower values are shown. Therefore, depending on the exact coupling configurations selected, you can reasonably expect values to fall within the shaded area.

# Specifications • • •

#### PRESSURE:

Vacuum to 120 psi, 8.3 bar

#### TEMPERATURE:

-40°F to 180°F (-40°C to 82°C)

#### MATERIALS:

Main components and valves: Acetal

Thumb latch: Acetal

Valve springs: 316 stainless steel

External springs and pin: Stainless steel

**0-rings**: Buna-N

**COLOR:** Natural white with light gray latch

#### **TUBING SIZES:**

1/4" to 3/8" ID. 6.4mm to 9.5mm ID

**WARNING:** Pressure, temperature, chemicals, and operating environment can affect the performance of couplings. It is the customer's responsibility to test the suitability of CPC's products in their own application conditions.



# DID YOU KNOW?

# THE SHELF LIFE OF A CONNECTOR IS DEPENDENT UPON THE CONNECTOR

MATERIALS AND THE STORAGE CONDITIONS. CPC manufactures quick disconnect couplings and connections made from a variety of metal, plastic and elastomeric raw materials. In general, it is best to store any CPC product in an environment shielded from light at a temperature between 50°F and 90°F.

The shelf life of CPC couplings and connections is estimated to be 3 years minimum if the storage conditions comply with those shown above.

**Note:** For product containing lubricated o-rings, the lubricant will dry out over time. Prior to use, it may be advisable to relubricate external o-rings. For product containing silicone o-rings, if lube is allowable, use of a synthetic lube is recommended.

# Liquid Flow Rate Information for Couplings

The chart below shows the flow rate for CPC couplings. Each coupling was tested with water at 70°F (21°C). To determine flow rates for specific coupling configurations use the formula at the right.




- Q = Flow rate in gallons per minute
- C<sub>v</sub> = Average coefficient across various flow rates (see chart)
- ΔP= Pressure drop across coupling (psi)
- S = Specific gravity of liquid

<b>VALUES F</b>	OR APC CO	JUPLINGS
-----------------	-----------	----------

BODIES INSERTS	APC 20004	APCD 20004	APC 20006	APCD 20006	APC 22004	APCD 22004	APC 22006	APCD 22006	APC 24004	APCD 24004	APC 24006	APCD 24006
APC10004	0.4	0.4	1.2	0.9	1.0	0.9	2.0	0.9	1.9	0.8	1.7	0.8
APCD10004	0.36	0.4	1.0	0.7	1.0	0.7	0.8	0.7	1.2	0.6	0.8	0.6
APC10006	0.4	0.4	1.2	0.9	1.2	0.7	1.8	0.9	2.0	0.9	1.7	0.8
APCD10006	0.37	0.4	1.3	0.8	0.6	0.5	1.1	0.7	1.3	0.7	1.0	0.7
APC13004SH	0.38	0.5	0.5	0.4	0.4	0.4	0.5	0.4	0.5	0.4	0.4	0.4
APCD13004SH	0.38	0.5	0.4	0.4	0.4	0.4	0.4	0.4	0.4	0.4	0.4	0.4
APC13006SH	0.38	0.4	1.0	0.7	0.9	0.6	1.1	0.7	1.1	0.7	1.0	0.7
APCD13006SH	0.37	0.4	0.8	0.6	0.7	0.5	0.8	0.6	0.9	0.6	0.8	0.6
APC17004SH	0.38	0.4	0.8	0.6	1.0	0.7	0.9	0.7	0.8	0.6	0.7	0.6
APCD17004SH	0.38	0.4	0.7	0.5	0.8	0.6	0.7	0.6	0.7	0.5	0.6	0.5
APC17006SH	0.4	0.4	1.2	0.9	1.4	0.9	1.4	0.9	2.0	0.9	1.9	0.8
APCD17006SH	0.4	0.4	1.2	0.8	1.3	0.8	1.2	0.7	1.2	0.7	1.1	0.7




# **APC DIMENSIONS**




# Panel Dimensions

	PANEL	PANEL THICKNESS	PANEL	PANEL NUT
	OPENING	MAXMIN.	NUT HEX	THREAD
COUPLING BODIES	see drawing	.05– .02	13/16	11/16-24UNF



# **Coupling Bodies • ACETAL**



All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters. †NOTE: CPC's Ferruleless Polytube Fitting terminations do not require ferrules to achieve a secure connection, which makes them easier to use and reuse. PTF fittings are designed for semi-rigid tubing, i.e., polyethylene, nylon, polyurethane, etc.

# Accessories

DESCRIPTION
PANEL MOUNT GASKET REPLACEMENT: For sealing
panel mount bodies listed above

MATERIAL<br/>Buna-NPART NO.1830300

# **Coupling Inserts • ACETAL**

TERMINATION IN-LINE

PIPE THREAD

IN-LINE FERRULELESS POLYTUBE FITTING,

IN-LINE HOSE BARB

ELBOW FERRULELESS POLYTUBE FITTING,

PTF†

**ELBOW** 

HOSE BARB











TUBING/THREAD SIZE 1/4" NPT 1/4" BSPT 3/8" NPT 3/8" BSPT	METRIC EQ.	STRAIGHT THRU APC24004 APC24004BSPT APC24006 APC24006BSPT	SHUTOFF APCD24004 APCD24004BSPT APCD24006 APCD24006BSPT	.72 .72 .72 .72	B 1.43/1.77 1.63/1.77 1.50/1.64 1.51/1.65	O
,	. ,	APC20004 APC200M8 APC20006 APC200M10	APCD20004 APCD200M8 APCD20006 APCD200M10	.72 .72 .72 .72	1.86/1.99 1.79/1.93 1.78/1.91 1.79/1.93	•••••
	6.4mm ID 9.5mm ID	APC22004 APC22006	APCD22004 APCD22006	.60/.72 .63/.72	1.40/2.10 1.40/1.92	•••••
,	6.4mm OD, 4.3mm ID 9.5mm OD, 6.4mm ID	APC21004 APC21006	APCD21004 APCD21006	1.22 1.27	1.32/1.45 1.33/1.46	.91 .96
27 1 10	6.4mm ID 9.5mm ID	APC23004 APC23006	APCD23004 APCD23006	.63 .63	1.20/1.45 1.20/1.50	.96 .96

All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters Couplings are pictured with valves unless otherwise noted.

tNOTE: CPC's Ferruleless Polytube Fitting terminations do not require ferrules to achieve a secure connection, which makes them easier to use and reuse. PTF fittings are designed for semi-rigid tubing, i.e., polyethylene, nylon, polyurethane, etc.


54 cpcworldwide.com • 800-444-2474


# **BREAKAWAY SERIES CONNECTOR**



CPC's BreakAway® Series couplings provide safe and easy fluid transfer with protection from costly product loss and equipment damage. These high quality, leak-free connectors are designed to disconnect when pull force is applied without depressing a latching mechanism, allowing users to make a clean, fast and easy break at critical moments.

FEATURES	BENEFITS
Reliable break-away design	Rapid, automatic disconnection
Automatic dual shutoff valves	Safe and clean fluid handling
Streamlined, flush finish	Easy-to-use for maximum productivity





# Specifications • • •

# PRESSURE:

0 to 20 psi, 1.4 bar

#### TEMPERATURE:

0°F to 150°F (-18°C to 66°C)

#### MATERIALS:

Main components and valves: Acetal

Valve springs: 316 stainless steel

O-rings: Buna-N **Retention ring: PPS** 

COLOR: White or black


#### **TUBING SIZES:**


1/4" and 3/8" ID, 6.4mm and 9.5mm ID

WARNING: Pressure, temperature, chemicals, and operating environment can affect the performance of couplings. It is the customer's responsibility to test the suitability of CPC's products in their own application conditions.


These graphs are intended to give you a general idea of the performance capabilities of each product line. The shaded area of the graphs represent the operating range of the product family, i.e., upper and lower values are shown. Therefore, depending on the exact coupling configurations selected, you can reasonably expect values to fall within the shaded area.

Force to disconnect is approximately 11 pounds pull force at 0 psi line pressure and 8 pounds pull force at 20 psi. In some applications and with certain types of tubing, a hose clamp may be required.





# **BREAKAWAY DIMENSIONS**



# Panel Dimensions

	PANEL OPENING	PANEL THICKNESS MAXMIN.	
IPLING	see drawing	.5002	
IIES			Ø .6



# **Coupling Bodies**

#### • WHITE ACETAL


	TERMINATION IN-LINE HOSE BARB	<b>TUBING SIZE</b> 1/4" ID 3/8" ID	<b>METRIC EQ.</b> 6.4mm ID 9.5mm ID	STRAIGHT THRU	SHUTOFF BACD17004 BACD17006	0.92 0.92	B 2.29 2.29	
0	PANEL MOUNT HOSE BARB	1/4" ID 3/8" ID	6.4mm ID 9.5mm ID		BACD16004 BACD16006	0.92 0.92	2.29 2.29	

#### • BLACK ACETAL

	TERMINATION IN-LINE HOSE BARB	<b>TUBING SIZE</b> 1/4" ID 3/8" ID	METRIC EQ. 6.4mm ID 9.5mm ID	STRAIGHT THRU	SHUTOFF BACD17004MBLK BACD17006MBLK	A 0.92 0.92	B 2.29 2.29	
Real Property of the Property	PANEL MOUNT Hose Barb	1/4" ID 3/8" ID	6.4mm ID 9.5mm ID		BACD16004MBLK BACD16006MBLK	0.92 0.92	2.29 2.29	

# **Coupling Inserts**

# WHITE ACETAL



#### • BLACK ACETAL




All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters. MBI K = molded black material.

# **EFC12 SERIES CONNECTOR**




The 9/32" flow EFC12 Series couplings feature a high efficiency valve design that provides a greater flow capability than any other coupling its size. Chemically resistant polypropylene material makes it ideal for harsh environments. The EFC12 Series adds a bulkhead panel mount option for tight seals against tank walls and drums.

# FEATURES BENEFITS High efficiency valve More flow than PLC Series in a compact size Plastic thumb latch Fewer moving parts Chemically resistant and gamma sterilizable



Compatible



Mates with most APC couplings

# Specifications • • •

#### PRESSURE:

Vacuum to 105 psi, 7.2 bar

#### TEMPERATURE:

32°F to 160°F (0°C to 71°C)

#### **MATERIALS:**

Main components and valves: Polypropylene

Thumb latch: Polypropylene
Valve spring: 316 stainless steel
Panel mount gasket: EPDM
External springs: 302 stainless steel

**0-rings:** EPDM

#### COLOR:


Gray with dark gray latch

#### **TUBING SIZES:**

1/4" and 3/8" ID, 6.4mm and 9.5mm ID

**WARNING:** Pressure, temperature, chemicals, and operating environment can affect the performance of couplings. It is the customer's responsibility to test the suitability of CPC's products in their own application conditions.a

These graphs are intended to give you a general idea of the performance capabilities of each product line. The shaded area of each graph represents the operating range of the product family, i.e., upper and lower values are shown. Therefore, depending on the exact coupling configurations selected, you can reasonably expect values to fall within the shaded area.



# Liquid Flow Rate Information for Couplings

The chart below shows the flow rate for CPC couplings. Each coupling was tested with water at 70°F (21°C). To determine flow rates for specific coupling configurations use the formula at the right.

# $Q=C_v\sqrt{\frac{\Delta P}{S}}$

Q = Flow rate in gallons per minute


C<sub>V</sub> = Average coefficient across various flow rates (see chart) ΔP = Pressure drop across coupling (psi)

S = Specific gravity of liquid

<b>C.</b>	VAII	IFS F	:NR	EFC12	COUP	LINGS
UΝ	INL	LUI	UIV	LIUIZ	0001	LINUS

BODIES INSERTS	EFC12 20004	EFCD12 20004	EFC12 20006	EFCD12 20006	EFC12 22004	EFCD12 22004	EFC12 22006	EFCD12 22006	EFC12 24004	EFCD12 24004
EFCD10412	0.51	0.51	0.51	0.51	0.50	0.45	0.50	0.50	0.51	0.51
EFCD10612	0.61	0.51	1.13	0.72	0.50	0.45	0.81	0.69	0.51	0.72
EFCD16412	0.51	0.51	0.51	0.51	0.50	0.45	0.50	0.50	0.51	0.51
EFCD16612	0.61	0.51	1.13	0.72	0.50	0.45	0.81	0.69	0.51	0.72
EFCD17412	0.51	0.51	0.51	0.51	0.50	0.45	0.50	0.50	0.51	0.51
EFCD17612	0.61	0.51	1.13	0.72	0.50	0.45	0.81	0.69	0.51	0.72

# **EFC12 DIMENSIONS**



# Panel Dimensions

	PANEL OPENING	PANEL THICKNESS MAXMIN.	PANEL NUT HEX	PANEL NUT THREAD	
COUPLING BODIES	see drawing	.25 — .03	13/16	11/16-24UN	
Ø .72"	.78"	EFC Series Gasket Thickness: . Mounting Hole: .720 Coupling Spacing: I Greenlee® 720 Keyv 730BB-3/4 Hole Pu	)" diameter 1.25" min. vay Punch a	nd	

# **Coupling Bodies • POLYPROPYLENE**

September 1	TERMINATION IN-LINE PIPE THREAD	TUBING/THREAD SIZE 1/4" NPT 3/8" NPT	METRIC EQ.	STRAIGHT THRU	SHUTOFF EFCD10412 EFCD10612	.93 .93	B 2.29 2.29	0	
S) III	BULKHEAD PANEL MOUNT HOSE BARB	1/4" ID 3/8" ID	6.4mm ID 9.5mm ID		EFCD16412 EFCD16612	.93 .93	2.23 2.23		
0/100	IN-LINE Hose Barb	1/4" ID 3/8" ID	6.4mm ID 9.5mm ID		EFCD17412 EFCD17612	.93 .93	2.23 2.23		

# **Coupling Inserts • POLYPROPYLENE**

All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters. MBLK = molded black material.



# Accessories

DESCRIPTION PANEL MOUNT GASKET REPLACEMENT	For sealing panel mount bodies listed above	<i>MATERIAL</i> Buna-N	<b>PART NO.</b> 1830300

# **NS4 SERIES CONNECTOR**



NS4 Series couplings feature non-spill valves in a compact size, at a great price. Use the NS4 when even a few drops pose problems regarding safety, media cost or environmental regulations. These innovative couplings are lightweight, chemically resistant and easy to use. The non-spill design effectively eliminates spills, minimizes downtime and enhances operator safety. NS4 coulings are also available with optional RFID (Radio Frequency Identification) capability (see page 97).

FEATURES	BENEFITS		
Non-spill design	Disconnect under pressure with no spills		
Color coding	Instant visual differentiation of media line		
Glass-filled polypropylene	Durable and compatible with many chemicals		
Medical-grade ABS	Gamma sterilizable		

These graphs are intended to give you a general idea of the performance capabilities of each product line. The shaded area of each graph represents the operating range of the product family, i.e., upper and lower values are shown. Therefore, depending on the exact coupling configurations

selected, you can reasonably expect values to fall within the

shaded area

# **NS4 Pressure Range**

# Specifications • • •

#### PRESSURE:

Vacuum to 120 psi, 8.3 bar

#### TEMPERATURE:

32°F to 160°F (0°C to 71°C)

#### MATERIALS:

#### Main components and valves:

Glass-filled polypropylene with TPV\* overmold, ABS with TPE\* soft-touch overmold

**Thumb latch:** Glass-filled polypropylene, ABS

Valve spring: 316 stainless steel External spring: 316 stainless steel

O-rings: EPDM

#### COLOR:

**Polypropylene:** Gray with dark gray overmold standard; gray with red or blue overmold available†

ABS: White with teal overmold

#### **TUBING SIZES:**

1/8" to 3/8" ID, 3.2mm to 9.5mm ID

**LUBRICANTS:** Krytox® PFPE (inert)

# SPILLAGE:


< 0.10 cc per disconnect at all rated pressures

#### **INCLUSION:** 0.26 cc per connect


WARNING: Pressure, temperature, chemicals, and operating environment can affect the performance of couplings. It is the customer's responsibility to test the suitability of CPC's products in their own application conditions

\*The overmold material TPV (themoplastic vulcanizate) is used with the NS4 polypropylene couplings. TPV is an alloy of polypropylene thermoplastic and fully vulcanized EPDM rubber. TPV is typically resistant to water, acids and bases. The overmold material TPE (thermoplastic elastomer) is used with the NS4 ABS couplings. TPE is a blend of additives and copolymers in a special formulation that forms extremely durable bonds to the ABS substrate, while offering the traditional properties of soft-touch overmold.

tNOTE: Standard product is grav: color options require a minimum quantities. Please contact CPC for details.



# **NS4 DIMENSIONS**




# **Panel Dimensions**

PANEL PANEL THICKNESS OPENING MAX.-MIN. COUPLING .25 - .03see drawing INSERTS

EPDM

1879800

# **Coupling Bodies** • POLYPROPYLENE / ABS



# **Coupling Inserts • POLYPROPYLENE / ABS**

	TERMINATION IN-LINE PIPE THREAD	TUBING/THREAD SIZE 1/4" NPT 1/4" BSPT	METRIC EQ.	<b>ABS SHUTOFF</b> NS4D2400406 NS4D24004BSPT06	POLYPROPYLENE SHUTOFF NS4D24004 NS4D24004BSPT	.96 .96	B 1.95 1.95
	IN-LINE FERRULESS POLYTUBE FITTING, PTF†	1/4" OD, .17" ID 3/8" OD, 1/4" ID	6.4mm OD, 4.3mm ID 9.5mm OD, 6.4mm ID	NS4D2000606	NS4D20004 NS4D20006	.96 .96	1.95 2.09
and the	IN-LINE Hose Barb	1/8" ID* 1/4" ID 3/8" ID	3.2mm ID* 6.4mm ID 9.5mm ID	NS4D2200206 NS4D2200406 NS4D2200606	NS4D22002 NS4D22004 NS4D22006	.96 .96 .96	1.73 1.98 1.98
	PANEL MOUNT FERRULESS POLYTUBE FITTING, PTF†	1/4" OD, .17" ID 3/8" OD, 1/4" ID	6.4mm OD, 4.3mm ID 9.5mm OD, 6.4mm ID	NS4D4000606	NS4D40004 NS4D40006	1.17 1.17	2.30 2.43
The same	PANEL MOUNT Hose Barb	1/8" ID* 1/4" ID 3/8" ID	3.2mm ID* 6.4mm ID 9.5mm ID	NS4D4200206 NS4D4200406 NS4D4200606	NS4D42002 NS4D42004 NS4D42006	1.17 1.17 1.17	2.07 2.32 2.32
all y	PANEL MOUNT Hose Barb Elbow	1/4" ID 3/8" ID	6.4mm ID 9.5mm ID		NS4D43004 NS4D43006	1.19 1.19	2.39 2.39
And the same	ELBOW HOSE BARB	1/4" ID 3/8" ID	6.4mm ID 9.5mm ID		NS4D23004 NS4D23006	0.96 0.96	2.04 2.04
ÎD	IN-LINE HOSE BARB (RS-232 Interface)	1/8" ID 1/4" ID 3/8" ID	3.2mm 6.4mm 9.5mm		ins4dt2200200 ins4dt2200400 ins4dt2200600	For production dimension cpcworldv	ns visit
	Access	ories	DESCRIPTION		MATERIAL	PART N	0.

All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters. PTF fittings are designed for semi-rigid tubing, i.e., polyethylene, nylon, polyurethane, etc. †NOTE: CPC's Ferruleless PTF (polytube fitting) terminations do not require ferrules to achieve a secure connection and are therefore easier to use and reuse. \*NOTE: Patented 1/8" ID (3.2mm) hose barb

panel mount inserts listed above

PANEL MOUNT GASKET REPLACEMENT: For sealing


connection has a support shroud allowing a maximum tube OD of 1/4" (6.4mm). 60 cpcworldwide.com • 800-444-2474 Smart fluid handling to take you forward, faster. 61


# LQ4 SERIES CONNECTOR



**LQ4 Series quick disconnect couplings** provide ultra-reliable, dripless connections and disconnections that protect valuable electronics. Designed specifically for liquid cooling applications, LQ4 utilizes patent-pending valve technology that eliminates drips and is able to withstand long-term use.

FEATURES	BENEFITS
Non-spill design	Disconnect under pressure with no spills
Redundant sealing when connected	Extra protection from leaks
Color-coding	Instant visual difference in cooling lines
Friction-free valve	Valving system designed for extended periods in connected state
Ruggedness	Able to withstand long-term, repeated use
Rotation allows movement	Eliminates tube kinking





# Specifications • • •

#### PRESSURE:

Vacuum to 120 psi, 8.3 bar

#### TEMPERATURE:

0°F to 240°F (-17°C to 115°C)

#### MAXIMUM FLOW AT DISCONNECT:

3.0 gal/min, 11.3L/min

#### MATERIALS:

Main housing components:

Chrome-plated brass

Valves and thumb latch: Polysulfone

Valve spring (wetted): Stainless steel

External spring: Stainless steel

Seals: EPDM

#### COLOR:

Chrome with Cool Blue or Warm Red

#### **TUBING SIZES:**

1/4" to 3/8" ID, 6.4mm to 9.5mm ID


**LUBRICANTS:** Krytox® PFPE


# SPILLAGE:

0.025 cc per disconnect rated at 0 psi 0.055 cc per disconnect rated at 120 psi


INCLUSION: 0.025 cc per connect

**WARNING:** Pressure, temperature, chemicals, and operating environment can affect the performance of couplings. It is the customer's responsibility to test the suitability of CPC's products in their own application conditions.





# LQ4 DIMENSIONS



# **Coupling Bodies • CHROME-PLATED BRASS**

	TERMINATION IN-LINE HOSE BARB	TUBING/THREAD SIZE 1/4" ID 1/4" ID 3/8" ID 3/8" ID	METRIC EQ. 6.4mm ID 6.4mm ID 9.5mm ID 9.5mm ID	SHUTOFF LQ4D17004RED LQ4D17004BLU LD4D17006RED LQ4D17006BLU	0.98 0.98 0.98 0.98	2.49 2.49 2.49 2.49 2.49
a second	IN-LINE STRAIGHT THREAD PORT SAE	1/4 SAE-4 1/4 SAE-4 3/8 SAE-6 3/8 SAE-6		LQ4D30004RED LQ4D30004BLU LQ4D30006RED LQ4D30006BLU	0.98 0.98 0.98 0.98	2.29 2.29 2.29 2.29
	IN-LINE Pipe Thread	1/4" NPT 1/4" NPT 3/8" NPT 3/8" NPT		LQ4D10004RED LQ4D10004BLU LQ4D10006RED LQ4D10006BLU	0.98 0.98 0.98 0.98	2.31 2.31 2.31 2.31 2.31
A P	IN-LINE PTF	1/4" ID x 3/8" OD 1/4" ID x 3/8" OD		LQ4D13006RED LQ4D13006BLU	0.98 0.98	2.46 2.46

# **Coupling Inserts • CHROME-PLATED BRASS**

	TERMINATION IN-LINE HOSE BARB	<b>TUBING/THREAD SIZE</b> 1/4"  D 1/4"  D 3/8"  D 3/8"  D	METRIC EQ. 6.4mm ID 6.4mm ID 9.5mm ID 9.5mm ID	SHUTOFF LQ4D22004RED LQ4D22004BLU LQ4D22006BLU LQ4D22006RED	0.74 0.74 0.74 0.74	B 2.22 2.22 2.22 2.22 2.22
A.	IN-LINE STRAIGHT THREAD PORT SAE	1/4 SAE-4 1/4 SAE-4 3/8 SAE-6 3/8 SAE-6		LQ4D46004BLU LQ4D46004RED LQ4D46006BLU LQ4D46006RED	0.80 0.80 0.80 0.80	2.02 2.02 2.02 2.02
	IN-LINE Pipe Thread	1/4" NPT 1/4" NPT 3/8" NPT 3/8" NPT		LQ4D24004RED LQ4D24004BLU LQ4D24006RED LQ4D24006BLU	0.80 0.80 0.80 0.80	2.04 2.04 2.04 2.04
01 F	IN-LINE PTF	1/4" ID x 3/8" OD 1/4" ID x 3/8" OD		LQ4D20006RED LQ4D20006BLU	0.80 0.80	2.19 2.19

All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters.

# **NS6 SERIES CONNECTOR**

**NS6 Series couplings** couplings feature non-spill valves at a great price. Use the NS6 when even a few drops of spillage pose problems regarding safety, media cost or environmental regulations. These innovative couplings are lightweight, chemically resistant and easy to use. The non-spill design virtually eliminates spills, minimizes downtime and enhances operator safety. Soft touch overmold makes them comfortable in the hand and very attractive.

FEATURES	BENEFITS		
Non-spill design	Disconnect under pressure with no spills		
Color coding	Instant visual differentiation of media lines		
Glass-filled polypropylene	Durable and compatible with many chemicals		
CPC thumb latch	One-hand connection and disconnection		

shaded area.

# **NS6 Pressure Range** Temperature (°C) -17.8 -6.7 4.4 15.6 26.7 37.8 48.9 60.0 71.1 These graphs are intended to give you a general idea of the performance capabilities of each product line. The shaded area of each graph represents the operating range of the product family, i.e., upper and lower values are shown. Therefore, depending on the exact coupling configurations selected, you can reasonably expect values to fall within the

# Specifications • • •

#### PRESSURE:

Vacuum to 120 psi, 8.3 bar

#### TEMPERATURE:

32°F to 160°F (0°C to 71°C)

#### **MATERIALS:**

#### Main components and valves:

Glass-filled polypropylene with TPV\* soft touch

Thumb latch: Glass-filled polypropylene Valve spring (wetted): 316 stainless steel External spring: 316 stainless steel O-rings: EPDM

# COLOR:

Gray with dark gray overmold standard; gray with red or blue overmold available†

#### **TUBING SIZES:**

3/8" and 1/2" ID, 9.5mm and 12.7mm ID

**LUBRICANTS:** Krytox® PFPE (inert)


# SPILLAGE:

~0.03 cc per disconnect @ 0 psi, ~0.30 cc/ disconnect @ 120 psi

# INCLUSION: 0.42 cc per connect

\*The overmold material is a TPV (thermoplastic vulcanizate). This TPV is an alloy of polypropylene thermoplastic and fully vulcanized EPDM rubber. The material is typically resistant to water, acids and bases.

†NOTE: Standard product is gray; color options require a set-up charge and minimum quantities. Please contact CPC



# **NS6 DIMENSIONS**

# A = Height/Diameter **B** = Total Length

# **Panel Dimensions**

COUPLING

PANEL OPENING	PANEL THICKNESS MAXMIN.	
see drawing	.25 – .03	
	1.250"	

# **Coupling Bodies • POLYPROPYLENE**

920	TERMINATION IN-LINE PIPE THREAD	TUBING/THREAD SIZE 1/2" NPT 1/2" BSPT	METRIC EQ.	<b>SHUTOFF</b> NS6D10008 NS6D10008BSPT	1.31 1.31	B 3.01 3.01
O.O.	IN-LINE COMPRESSION	1/2" OD	12.7mm OD	NS6D13008	1.31	3.53
000	IN-LINE HOSE BARB	3/8" ID 1/2" ID	9.5mm ID 12.7mm ID	NS6D17006 NS6D17008	1.31 1.31	3.16 3.16

# **Coupling Inserts • POLYPROPYLENE**

	TERMINATION IN-LINE PIPE THREAD	<b>TUBING/THREAD SIZE</b> 1/2" NPT 1/2" BSPT	METRIC EQ.	SHUTOFF NS6D24008 NS6D24008BSPT	1.31 1.31	B 2.44 2.44
	IN-LINE COMPRESSION	1/2" OD	12.7mm OD	NS6D20008	1.31	3.01
· Die	IN-LINE HOSE BARB	3/8" ID 1/2" ID	9.5mm ID 12.7mm ID	NS6D22006 NS6D22008	1.31 1.31	2.59 2.59
No septiment of the sep	PANEL MOUNT COMPRESSION	1/2" OD	12.7mm 0D	NS6D40008	1.50	3.32
· De	PANEL MOUNT HOSE BARB	3/8" ID 1/2" ID	9.5mm ID 12.7mm ID	NS6D42006 NS6D42008	1.50 1.50	2.85 2.85

# Accessories

DESCRIPTION	MATERIAL	PART NO.
PANEL MOUNT GASKET REPLACEMENT: For sealing	EPDM	1884300
panel mount bodies listed above	FKM	1889600


All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters. Couplings are pictured with valves unless otherwise noted.


# LQ6 SERIES CONNECTOR



**LQ6 Series couplings** closed-system design provides ultra-reliable, dripless connections and disconnections. The first of their kind, the couplings utilize patent-pending valve technology specifically designed for liquid cooling applications.

FEATURES	BENEFITS		
Non-spill design	Disconnect under pressure with no spills		
Redundant sealing when connected	Extra protection from leaks		
Color-coding	Instant visual difference in cooling lines		
Friction-free valve	Valving system designed for extended periods in connected state		
Ruggedness	Able to withstand long-term, repeated use		
Rotation allows movement	Eliminates tube kinking		





# Specifications • • •

#### PRESSURE:

Vacuum to 120 psi, 8.3 bar

#### TEMPERATURE:

0°F to 240°F (-17°C to 115°C)

#### MAXIMUM FLOW AT DISCONNECT:

4.0 gal/min, 15.1 L/min

#### **MATERIALS:**

Main housing components:

Chrome-plated brass

Valves and thumb latch: Polysulfone

Valve spring (wetted): Stainless steel

External spring: Stainless steel

Seals: EPDM

#### COLOR:

Chrome with Cool Blue or Warm Red

#### **TUBING SIZES:**

3/8" to 1/2" ID, 9.5mm to 12.7mm ID


**LUBRICANTS:** Krytox® PFPE


#### SPILLAGE:

0.05 cc per disconnect rated at 0 psi 0.125 cc per disconnect rated at 120 psi


**INCLUSION:** 0.42 cc per connect

**WARNING:** Pressure, temperature, chemicals, and operating environment can affect the performance of couplings. It is the customer's responsibility to test the suitability of CPC's products in their own application conditions.





# LQ6 DIMENSIONS




A = Height/DiameterB = Total Length

# **Coupling Bodies • CHROME-PLATED BRASS**

	TERMINATION IN-LINE HOSE BARB	TUBING/THREAD SIZE 3/8" ID 3/8" ID 1/2" ID 1/2" ID	METRIC EQ. 9.5mm ID 9.5mm ID 12.7mm ID 12.7mm ID	SHUTOFF LQ6D17006RED LQ6D17006BLU LQ6D17008RED LQ6D17008BLU	1.185 1.185 1.185 1.185	2.73 2.73 2.73 2.73 2.73
OF THE	IN-LINE STRAIGHT THREAD PORT SAE	3/8 SAE-6 3/8 SAE-6 1/2 SAE-8 1/2 SAE-8		LQ6D30006RED LQ6D30006BLU LQ6D30008RED LQ6D30008BLU	1.185 1.185 1.185 1.185	2.48 2.48 2.50 2.50
	IN-LINE PIPE THREAD	3/8" NPT 3/8" NPT		LQ6D10006RED LQ6D10006BLU	1.185 1.185	2.56 2.56
	IN-LINE PTF	3/8" ID x 1/2" OD 3/8" ID x 1/2" OD		LQ6D13008RED LQ6D13008BLU	1.185 1.185	2.50 2.50

# **Coupling Inserts • CHROME-PLATED BRASS**



All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters.

# **NSH SERIES CONNECTOR**

NSH Series couplings couplings are general purpose versions of our pressure-balanced, non-spill, ChemQuik® CQG06 couplings. Molded polypropylene, EPDM seals and a 100% springless and metal-free flow path provide broad chemical resistance and exceptionally high flow capacity, allowing instant disconnects (and reconnects), even under pressure. Their non-spill design virtually eliminates spills, minimizes downtime and enhances operator safety in a very low cost package.

FEATURES	BENEFITS		
Non-spill design	Ultimate protection from chemicals and fumes		
Pressure-balanced design	Failsafe disconnect, even under pressure; easy to reconnect at high pressure		
Springless flow path design	Eliminates source of metallic contaminant		
Optional color coding	Helps prevent accidental misconnects		


These graphs are intended to give you a general idea of the performance capabilities of each product line. The shaded area of each graph represents the operating range of the

product family, i.e., upper and lower values are shown.

shaded area.

Therefore, depending on the exact coupling configurations

selected, you can reasonably expect values to fall within the



# Specifications • • •

#### PRESSURE:

Vacuum to 120 psi, 8.3 bar

#### TEMPERATURE:

32°F to 120°F (0°C to 49°C)

#### MATERIALS:

Main components: Glass-filled polypropylene Thumb latch: Glass-filled polypropylene

**0-rings:** EPDM

Springs (non-wetted): 316 stainless steel Compression nuts: White polypropylene

#### COLOR:

Gray with dark gray accent standard; red or blue thumb latch and insert tip available.

#### TUBING SIZES:

3/8" to 3/4" ID, 9.5mm to 19.0mm ID

LUBRICANTS: Krytox® PFPE (inert)

SPILLAGE: 0.03 cc/connect

INCLUSION: 0.42 cc per connect

#### PANEL MOUNT:


Optional adaptor kit (see page 69)

WARNING: Pressure, temperature, chemicals, and operating environment can affect the performance of couplings. It is the customer's responsibility to test the suitability of CPC products in their own application conditions. Use the graphs

Compatible with ChemQuik® Dual Containment System; DrumQuik® Universal Dispense System and Asian Drum

# **NSH Water Flow** 45.4 68.1 90.8 113.6 22.7

# **NSH DIMENSIONS**



# **Panel Dimensions**

COUPLING INSERTS


OPENING	MAX.—MIN.	
see drawing	.50 – .12	
	 _	
	-	1
	 <u></u>	
2.000"		
	_	

PANEL THICKNESS

# **Coupling Bodies • POLYPROPYLENE**

18 m	TERMINATION IN-LINE PIPE THREAD	TUBING/THREAD SIZE 3/8" FEMALE NPT 3/8" BSPT 1/2" NPT 3/4" NPT	METRIC EQ.	SHUTOFF NSHD19006 NSHD10006BSPT NSHD10008 NSHD10012	2.0 2.0 2.0 2.0 2.0	3.70 3.67 3.81 4.01
9	IN-LINE COMPRESSION	3/8" OD 1/2" OD	9.5mm OD 12.7mm OD	NSHD13006 NSHD13008	2.0 2.0	3.97 4.10
9	IN-LINE HOSE BARB	3/8" ID 1/2" ID 5/8" ID 3/4" ID	9.5mm ID 12.7mm ID 15.9mm ID 19.0mm ID	NSHD17006 NSHD17008 NSHD17010 NSHD17012	2.0 2.0 2.0 2.0	3.73 3.73 3.88 3.88

# **Coupling Inserts • POLYPROPYLENE**



All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters. \*Panel mount kit can be used with body for front mount or insert for front or rear mount.

# Accessories



PANEL MOUNT KIT: with panel seal Dust plug/cap for body or insert\*


PART NO. NSHPMKIT12 1800700 (not shown)


# HFC12 SERIES CONNECTOR



**HFC12 Series couplings** have flow comparable to many 1/2" flow couplings in a 3/8" body size. Compact and lightweight, HFC couplings replace bulky and heavy brass ball-and-sleeve couplings in a wide range of applications. An ergonomic design and a large, shrouded thumb latch pad produce a coupling that is easy to grip and simple to operate. An efficient valve design leads to high flow and low spillage.

FEATURES BENEFITS	
High efficiency valve More flow in a compact size	
Ergonomic design	Easy to grip, simple to operate
Polypropylene material	Chemically resistant and gamma sterilizable
Compatible	Mates with HFC35 and HFC57 couplings





# Specifications • • •

#### PRESSURE:

Vacuum to 60 psi, 4.2 bar

#### TEMPERATURE:

32°F to 160°F (0°C to 71°C)

#### **MATERIALS:**

Main components and valves: Polypropylene

Thumb latch: Polypropylene

Valve spring (wetted): 316 stainless steel External springs: 316 stainless steel

**0-rings:** EPDM

Panel mount gasket: EPDM Compression nut, gripper, ferrule:

Polypropylene

**COLOR:** Gray with dark gray latch

#### **TUBING SIZES:**


3/8" to 3/4" ID, 9.5mm to 19.0mm ID


**WARNING:** Pressure, temperature, chemicals, and operating environment can affect the performance of couplings. It is the customer's responsibility to test the suitability of CPC's products in their own application conditions.

For compression termination specifications: 3/8" OD or 1/2" OD, +0.010/-0.000



These graphs are intended to give you a general idea of the performance capabilities of each product line. The shaded area of each graph represents the operating range of the product family, i.e., upper and lower values are shown. Therefore, depending on the exact coupling configurations selected, you can reasonably expect values to fall within the shaded area.





# Liquid Flow Rate Information for Couplings

C<sub>v</sub> VALUES FOR HFC12 COUPLINGS


The chart below shows the flow rate for CPC couplings. Each coupling was tested with water at  $70^{\circ}$ F ( $21^{\circ}$ C). To determine flow rates for specific coupling configurations use the formula at the right.



- Q = Flow rate in gallons per minute
- C<sub>v</sub> = Average coefficient across various flow rates (see chart)
- $\Delta P$  = Pressure drop across coupling (psi)
- S = Specific gravity of liquid

BODIES INSERTS	HFC 22612	HFCD 22612	HFC 22812	HFCD 22812	HFC 23612	HFCD 23612	HFC 23812	HFCD 23812	HFC 24612	HFCD 24612	HFC 24812	HFCD 24812	HFCD 221212	HFC 221212
HFCD10612	1.27	1.27	1.62	1.51	1.14	1.14	1.46	1.36	1.80	1.58	1.70	1.65	-	-
HFCD10812	1.28	1.34	1.62	1.51	1.15	1.24	1.46	1.36	1.81	1.54	1.72	1.56	-	-
HFCD16612	1.07	1.00	1.17	1.14	0.96	0.90	1.05	1.03	1.33	1.26	1.30	1.24	-	-
HFCD16812	1.25	1.23	1.61	1.52	1.13	1.11	1.45	1.37	1.79	1.60	1.68	1.56	-	-
HFCD17612	1.07	1.00	1.17	1.14	0.96	0.90	1.05	1.03	1.33	1.26	1.30	1.24	-	-
HFCD17812	1.25	1.23	1.61	1.52	1.13	1.11	1.45	1.37	1.79	1.60	1.68	1.56	-	-
HFC171212	-	-	-	<u> </u>	<b>///</b>	\\\-\\\\-\\\\	VX-	(	-	-	-	-	3.94	-
HFCD171212	-	-	-		/ / /- T			// /-	X	-	-	-	-	2.04

# HFC12 DIMENSIONS



# **Coupling Bodies • POLYPROPYLENE**




All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters.

# Accessories



DESCRIPTION	MATERIAL	PART NO
PANEL MOUNT GASKET REPLACEMENT: For sealing	EPDM	621200
panel mount bodies listed above		
panel mount bodies listed above		

# **Coupling Inserts • POLYPROPYLENE**



All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters. Couplings are pictured with valves unless otherwise noted

# DID YOU KNOW?

WHEN SELECTING A VALVED COUPLING FOR YOUR APPLICATION, make sure to order a part number with a D just before the numeric portion of the part number. For example, HFC101212 does not have a shutoff valve; HFCD101212 would be the correct part number to order for a valved coupling.

To visually identify a part to determine whether or not it is valved, disconnect the coupling body and insert, and then look through the part. If you can see light all the way through, your part is non-valved. If there is anything obstructing the light, you probably have a valved part. Please contact Customer Service at 1-800-444-2474 or 651-645-0091 if you need further assistance.

# HFC35 & HFC57 SERIES CONNECTOR



**HFC35 Series polysulfone couplings** are able to withstand demanding applications and feature physical strength, chemical resistance and autoclavability.

**HFC57 Series polysulfone couplings** are UV resistant and able to withstand continued exposure to harmful UV rays without affecting performance.

Both can operate in harsh environments in temperatures from -40°F to 280°F (-40°C to 138°C) and pressures from vacuum to 125 psi (8.6 bar). An ergonomic design and a large, shrouded thumb latch pad produce a coupling that is easy to grip and simple to operate. An efficient valve design leads to high flow and exceptionally low spillage. Both series also add bulkhead panel mount, garden hose thread and Hastelloy® C spring options.

FEATURES	BENEFITS
Polysulfone (HFC35) Physically strong and autoclavable	
UV Polysulfone (HFC57) Physically strong and UV resistant	
High efficiency valve	More flow and less spillage in a compact size
Ergonomic design	Easy to grip, simple to operate
Compatible	Mates with HFC12 couplings

#### Specifications • • •

#### PRESSURE:

Vacuum to 125 psi, 8.6 bar

#### TEMPERATURE:

-40°F to 280°F (-40°C to 138°C)

#### MATERIALS:

#### Main components and valves:

HFC35 - Polysulfone

HFC57 - UV-resistant polysulfone

#### THUMB LATCH:

HFC35 - Polysulfone

HFC57 - UV-resistant polysulfone
Valve spring: 316 stainless steel
External springs: 316 stainless steel

O-rings: EPDM

Panel mount gasket: EPDM Female GHT gasket: FDA EPDM

#### COLORS:

HCF35 - Natural white, dark gray latch


HCF57 - Black, black latch


#### **TUBING SIZES:**


3/8" to 3/4" ID, 9.5mm to 19.0mm ID

**WARNING:** Pressure, temperature, chemicals, and operating environment can affect the performance of couplings. It is the customer's responsibility to test the suitability of CPC products in their own application conditions. Use the graphs on the following page as a guide.









These graphs are intended to give you a general idea of the performance capabilities of each product line. The shaded area of each graph represents the operating range of the product family, i.e., upper and lower values are shown. Therefore, depending on the exact coupling configurations selected, you can reasonably expect values to fall within the shaded area.

# DID YOU KNOW?

25 50 75 100 125 150 175 200 225 250 Flow (scfm)


#### THE HFC57 SERIES POLYSULFONE COUPLINGS ARE UV-RESISTANT AND ABLE TO WITHSTAND CONTINUED EXPOSURE TO HARMFUL RAYS WITHOUT AFFECTING PERFORMANCE?

These couplings offer all the same features as our HFC12 and HFC35 Series with the added advantage of UV-resistance to protect the strength and durability of your fluid connections. Use the HFC57 couplings in outdoor applications such as watercraft, lawn irrigation systems and solar heating on roofs and decks where a connector is subjected to the harsh rays of the sun.



# HFC35 DIMENSIONS

# **HFC35 Coupling Bodies** • POLYSULFONE



#### **HFC35 Coupling Inserts • POLYSULFONE**


	TERMINATION IN-LINE PIPE THREAD (non-valved shown)	TUBING/THREAD : 3/8" NPT 3/8" BSPT 1/2" NPT 1/2" BSPT 3/4" NPT 3/4" BSPT 3/4" GHT	SIZE METRIC EQ.	STRAIGHT THRU HFC24635 HFC24635BSPT HFC24835 HFC24835BSPT HFC241235 HFC241235BSPT HFC241235GHT	SHUTOFF HFCD24635 HFCD24635BSPT HFCD24835 HFCD24835BSPT HFCD241235 HFCD241235BSPT HFCD241235GHT	1.16 1.16 1.16 1.16 1.23 1.23 1.23	B 1.82/1.94 1.97/2.09 1.97/2.09 1.17/2.09 2.22/2.34 2.21/2.33 2.03/2.15	D
	IN-LINE PIPE THREAD (FEMALE) (non-valved shown)	3/4" GHT 3/4" BSPP		HFC261235GHT HFC261235BSPP	HFCD261235GHT HFCD261235BSPP	1.37 1.37	1.79/1.91 1.79/1.91	
5	IN-LINE HOSE BARB (non-valved shown)	3/8" ID 1/2" ID 5/8" ID 3/4" ID	9.5mm ID 12.7mm ID 15.9mm ID 19.0mm ID	HFC22635 HFC22835 HFC221035 HFC221235	HFCD22635 HFCD22835 HFCD221035 HFCD221235	1.00 1.00 1.00 1.00	1.91/2.03 1.91/2.03 2.03/2.15 2.03/2.15	
	ELBOW HOSE BARB (non-valved shown)	3/8" ID 1/2" ID	9.5mm ID 12.7mm ID	HFC23635 HFC23835	HFCD23635 HFCD23835	1.00 1.00	1.92/2.04 1.96/2.08	.93 .93

#### ♠ HFC35 & HFC57 Accessories ≫

<b>DESCRIPTION PANEL MOUNT GASKET REPLACEMENT:</b> For sealing panel mount bodies listed above	<i>MATERIAL</i> EPDM	<b>PART NO.</b> 621200
GHT GASKET REPLACEMENT: For sealing female threads	FDA EPDM	2339400


All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters.

# HFC57 DIMENSIONS



#### **Panel Dimensions**

	PANEL OPENING	PANEL THICKNESS MAXMIN.
COUPLING BODIES	see drawing	.2503



# **HFC57 Coupling Bodies** • UV POLYSULFONE

TERMINATION IN-LINE PIPE THREAD	TUBING/THREAD SIZE 3/8" BSPT 1/2" NPT 1/2" BSPT 3/4" BSPT 3/4" GHT	METRIC EQ.	STRAIGHT THRU HFC10657BSPT HFC10857 HFC10857BSPT HFC101257BSPT HFC101257GHT	SHUTOFF HFCD10657BSPT HFCD10857 HFCD10857BSPT HFCD101257BSPT HFCD101257GHT	1.44 1.44 1.44 1.44 1.44	2.73 2.88 2.88 3.12 2.94
IN-LINE Pipe Thread (Female)	3/4" GHT 3/4" BSPP		HFC191257GHT HFC191257BSPP	HFCD191257GHT HFCD191257BSPP	1.44 1.44	2.70 2.70
BULKHEAD Panel mount Hose Barb	3/8" ID 1/2" ID 5/8" ID 3/4" ID	9.5mm ID 12.7mm ID 15.9mm ID 19.0mm ID		HFCD16657 HFCD16857 HFCD161057 HFCD161257	1.82 1.82 1.82 1.82	2.82 2.82 2.94 2.94
IN-LINE HOSE BARB	3/8" ID 1/2" ID 5/8" ID 3/4" ID	9.5mm ID 12.7mm ID 15.9mm ID 19.0mm ID	HFC17657 HFC17857 HFC171057 HFC171257	HFCD17657 HFCD17857 HFCD171057 HFCD171257	1.44 1.44 1.44 1.44	2.82 2.82 2.94 2.94

# **HFC57 Coupling Inserts • UV POLYSULFONE**

TERMINATION IN-LINE PIPE THREAD (non-valved shown)	TUBING/THREAD SIZ 3/8" BSPT 1/2" BSPT 1/2" NPT 3/4" BSPT 3/4" GHT	-	STRAIGHT THRU HFC24657BSPT HFC24857BSPT HFC24857 HFC241257BSPT HFC241257GHT	SHUTOFF HFCD24657BSPT HFCD24857BSPT HFCD24857 HFCD241257BSPT HFCD241257GHT	1.16 1.16 1.16 1.23 1.16	1.84/1.96 1.99/2.07 1.95/2.11 2.19/2.31 2.15/2.27
IN-LINE PIPE THREAD (FEMALE) (non-valved shown)	3/4" GHT 3/4" BSPP		HFC261257GHT HFC261257BSPP	HFCD261257GHT HFCD261257BSPP	1.37 1.37	1.79/1.91 1.79/1.91
IN-LINE HOSE BARB (non-valved shown)	3/8" ID 1/2" ID 5/8" ID 3/4" ID	9.5mm ID 12.7mm ID 15.9mm ID 19.0mm ID	HFC22657 HFC22857 HFC221057 HFC221257	HFCD22657 HFCD22857 HFCD221057 HFCD221257	1.00 1.00 1.00 1.00	1.91/2.03 1.91/2.03 2.03/2.15 2.03/2.15

All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters. Couplings are pictured with valves unless otherwise noted.

# FFC35 SERIES CONNECTOR



# Specifications • • • PRESSURE:

Vacuum to 125 psi, 8.6 bar

#### TEMPERATURE:

-40°F to 280°F (-40°C to 138°C)

#### **MATERIALS:**

Main components: Polysulfone Thumb latch: Polysulfone

External spring: 316 stainless steel

O-rings: FDA Buna-N Panel mount gasket: EPDM Female GHT gasket: FDA EPDM

#### COLOR:

Natural white with dark gray latch

#### TUBING SIZES: 3/4" ID, 19.0mm ID

WARNING: Pressure, temperature, chemicals, and operating environment can affect the performance of couplings. It is the customer's responsibility to test the suitability of CPC products in their own application conditions. Use the graphs on the following page as a guide.


FFC35 Series 1/2" flow couplings feature double the flow of the HFC Series couplings. The FFC Series has a non-valved bore that increases flow and minimizes turbulence. The polysulfone material is able to withstand demanding applications and features physical strength, chemical resistance and autoclavability. An ergonomic design and a large, shrouded thumb latch pad produce a coupling that is easy to grip and simple to operate.

FEATURES BENEFITS		
o valves Exceptional flow		
Polysulfone material	Physically strong and autoclavable	
Ergonomic design	Easy to grip, simple to operate	

#### shaded area of each graph represents the operating range of the product family, i.e., upper and lower values are shown. Therefore, depending on the exact coupling configurations selected, you can reasonably expect values to fall within the shaded area. FFC Water Flow

These graphs are intended to give you a general idea of

the performance capabilities of each product line. The




#### NON-VALVED


CPC's FFC35 Series couplings feature a wide open, straight-thru flow path.







# FFC35 DIMENSIONS



#### **Panel Dimensions**


	PANEL OPENING	PANEL THICKNESS MAX.—MIN.
COUPLING BODIES	see drawing	.25– .03
Ø 1.21"	1.30"	Gasket Thickness: .06* Mounting Hole: 1.21* diameter Coupling Spacing: 1.92* minimum

#### **Coupling Bodies • POLYSULFONE**



# **Coupling Inserts • POLYSULFONE**

All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters. Couplings are pictured with valves unless otherwise noted.



TUBING/THREAD SIZE 1/2" NPT 1/2" BSPT 3/4" NPT 3/4" BSPT 3/4" GHT	METRIC EQ.	STRAIGHT THRU FFC24835 FFC24835BSPT FFC241235 FFC241235BSPT FFC241235GHT	1.16 1.16 1.23 1.23 1.16	1.95 1.95 2.14 2.14 2.01
3/4" GHT 3/4" BSPP		FFC261235GHT FFC261235BSPP	1.37 1.37	1.77 1.77
3/4" ID	19.0mm ID	FFC221235	1.00	2.02

# Accessories

<b>DESCRIPTION PANEL MOUNT GASKET REPLACEMENT:</b> For sealing panel mount bodies listed above	<i>Material</i> EPDM	<b>PART NO.</b> 621200	
GHT GASKET REPLACEMENT: For sealing female threads	FDA EPDM	2339400	



The Tentube<sup>TM</sup> coupling provides one easy-to-use quick disconnect for up to ten separate fluid lines. The Tentube is available with valves on the body side in panel mount or in-line configurations for maximum design flexibility. The Tentube can connect and disconnect up to ten lines with the use of a simple slide latch. Tubing orientation is ensured by keying the coupling body and insert. Different fluids and tubing sizes can be accommodated by using acetal or polypropylene inserts in three popular sizes.

FEATURES	BENEFITS	
Ten line connection	Fast, efficient operation	
Separate flow paths	Pressure and/or vacuum in one coupli	
Insert is keyed to body	Ensures correct line orientation	
Acetal or polypropylene inserts	Maximize design flexibility	

# Specifications • • •

#### PRESSURE:

Vacuum to 100 psi, 6.9 bar per line

#### TEMPERATURE:

**Acetal fitting inserts:** 

-40°F to 180°F (-40°C to 82°C)

**Polypropylene fitting inserts:** 32°F to 180°F (0°C to 82°C)

#### **MATERIALS:**

Main components: Nylon

Valves: Acetal

Valve springs: 316 stainless steel Fitting inserts: Acetal or polypropylene

O-rings: Buna-N with acetal or EPDM with polypropylene Panel mount adapter: Acetal Tube shroud: Acetal


#### COLOR:


Main components: Black with red locking latch Tube shroud: Black Panel mount adapter: Black

#### **TUBING SIZES:**


1/16" to 1/8" ID, 1.6mm to 3.2mm ID

These graphs are intended to give you a general idea of the performance capabilities of each product line. The shaded area of each graph represents the operating range of the product family, i.e., upper and lower values are shown. Therefore, depending on the exact coupling configurations selected, you can reasonably expect values to fall within the shaded area.





# TENTUBE DIMENSIONS



#### **Panel Dimensions**

	PANEL OPENING	PANEL THICKNESS MAXMIN.
COUPLING BODIES	see drawing	.4406

# Fitting Bodies • ACETAL/POLYPROPYLENE



TERMINATION COMPLETE COUPLING SET

COUPLING BODY WITH FEMALE FITTING BODIES

COUPLING INSERT WITH MALE FITTING INSERTS

<b>TUBING SIZE</b> 1/16" ID 1/8" ID	METRIC EQ. 1.6mm ID 3.0mm ID 3.2mm ID	ACETAL STRAIGHT THRU TT1001 TT10M3 TT1002	ACETAL SHUTOFF TTD1001 TTD10M3 TTD1002	POLYPROPYLENE STRAIGHT THRU TT100112
1/16" ID  1/8" ID Without fitting bodies	1.6mm ID 3.0mm ID 3.2mm ID	TFB1001 TFB10M3 TFB1002 TFB10	TFBD1001 TFBD10M3 TFBD1002	TFB100112 TFB100212
1/16" ID  1/8" ID Without fitting inserts	1.6mm ID 3.0mm ID 3.2mm ID	TMB1001 TMB10M3 TMB1002 TMB10		TMB100112 TMB100212

# Fitting Inserts • ACETAL/POLYPROPYLENE



TERMINATION IN-LINE HOSE BARB

IN-LINE HOSE BARB (FEMALE)

<b>TUBING SIZE</b> 1/16" ID 1/8" ID	<b>METRIC EQ.</b> 1.6mm ID 3.0mm ID 3.2mm ID	AGETAL STRAIGHT THRU TMF01 TMFM3 TMF02	ACETAL SHUTOFF	POLYPROPYLENE STRAIGHT THRU TMF0112 TMF0212
1/16" ID	1.6mm ID 3.0mm ID	TFF01 TFFM3	TFFD01 TFFDM3	TFF0112
1/8" ID	3.2mm ID	TFF02	TFFD02	TFF0212

Note: Acetal bodies and inserts are also available in white. To order, add MWHT to the end of the part number.

#### Accessories




DESCRIPTION
TUBE SHROUD: snaps into either half

PART NO.
TS10
.....TPM10

**PANEL MOUNT ADAPTER & NUT:** 7/16" max. panel thickness; requires 1.875" diameter hole

All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters, Couplings are pictured with valves unless otherwise noted.

# TWIN TUBE & SIXTUBE SERIES CONNECTOR



# Twin Tube Specifications

#### PRESSURE:

Vacuum to 120 psi, 8.3 bar per line

#### TEMPERATURE:

-40°F to 180°F (-40°C to 82°C)

#### MATERIALS:

Main components: Acetal, ABS

Valves: Acetal

Thumb latch: Stainless steel

External springs and pin: Stainless steel

O-rings: Buna-N

**COLOR:** Natural white with aqua latch


#### **TUBING SIZES:**


1/16" to 1/8" ID, 1.6mm to 3.2mm ID

**The Twin Tube<sup>TM</sup> coupling** provides one easy-to-use quick disconnect for two separate fluid lines. The non-valved Twin Tube maintains two individual flow paths in one coupling. Featuring the CPC thumb latch, the Twin Tube can disconnect two lines with the push of a button. The panel mount design keeps a low profile on the front of the equipment panel.

**The Sixtube<sup>TM</sup> coupling** provides one easy-to-use quick disconnect for up to six separate fluid lines. The Sixtube is available with valves on either side for maximum design flexibility. Featuring the CPC thumb latch, the Sixtube can disconnect up to six lines with the push of a button. The snap-in panel mount design keeps a low profile on the front of the equipment panel.

FEATURES	BENEFITS	
Two/six line connection	Fast, efficient operation	
Separate flow paths	Pressure and/or vacuum in one coupli	
Sixtube insert is keyed to body	Ensures correct line orientation	
Free coupling rotation (Twin Tube)	Eliminates kinked tubing	
CPC thumb latch	One-hand connection and disconnection	





#### Sixtube Specifications • • •

#### PRESSURE:

Vacuum to 100 psi, 6.9 bar per line

#### **TEMPERATURE:**

#### **Acetal fitting inserts:**

-40°F to 180°F (-40°C to 82°C)

Polypropylene fitting inserts:

32°F to 180°F (0°C to 82°C)

#### MATERIALS:

Main components: Acetal

Valves: Acetal

Fitting inserts: Acetal, polypropylene

Thumb latch: Acetal

Valve spring: 316 stainless steel

**0-rings:** Buna-N with acetal or EPDM with

polypropylene


COLOR: Natural white with aqua latch

#### **TUBING SIZES:**

1/16" to 1/8" ID, 1.6mm to 3.2mm ID

**NOTE:** Caution must be observed when ordering Sixtube assemblies. Body and insert halves which both contain female or male fittings will not couple. Contact CPC for questions and ordering assistance.

These graphs are intended to give you a general idea of the performance capabilities of each product line. The shaded area of each graph represents the operating range of the product family, i.e., upper and lower values are shown. Therefore, depending on the exact coupling configurations selected, you can reasonably expect values to fall within the shaded area.




# TWINTUBE & SIXTUBE DIMENSIONS

#### Twin Tube™ Body Panel Mount Dimensions



Minimum Panel Thickness: 0.02"

#### Sixtube™ Insert Panel Mount Dimensions



#### Twintube Coupling Bodies • ACETAL/ABS



TERMINATION PANEL MOUNT HOSE BARB **TUBING SIZE** 1/16" ID 1/8" ID

METRIC EQ. 1.6mm ID 3.2mm ID ACETAL STRAIGHT THRU PTC16010 PTC16020

STRAIGHT THRU

PTC1602096

# Twintube Coupling Inserts • ACETAL/ABS



TERMINATION IN-LINE HOSE BARB **TUBING SIZE** 1/16" ID 1/8" ID

METRIC EQ. 1.6mm ID 3.2mm ID ACETAL STRAIGHT THRU PTC22010

PTC22020

ABS STRAIGHT THRU

PTC2202096

SXFD1702

# Sixtube Coupling Bodies • ACETAL/POLYPROPYLENE



TERMINATION
WITH MALE FITTING
INSERTS

WITH FEMALE FITTING BODIES

WITHOUT FITTINGS

٧G		
• • • •		

 TUBING SIZE
 METRIC EQ.

 1/16" ID
 1.6mm ID

 3.0mm ID
 3.2mm ID

 1/16" ID
 1.6mm ID

 3.0mm ID
 3.0mm ID

 3.2mm ID
 3.2mm ID

ACETAL ABS STRAIGHT THRU SHUTOFF SXM1701 SXM17M3 SXM1702 SXF1701 SXFD1701 SXF17M3 SXFD17M3

SX17

SXF1702

**ACETAL** 

SXM4201

SXM42M3

STRAIGHT THRII

# Sixtube Coupling Inserts • ACETAL/POLYPROPYLENE



TERMINATION WITH MALE FITTING INSERTS

WITH FEMALE FITTING BODIES

WITHOUT FITTINGS

 TUBING SIZE
 METRIC EQ.

 1/16" ID
 1.6mm ID

 3.0mm ID
 3.2mm ID

 1/16" ID
 1.6mm ID

1/8" ID

3.0mm ID 3.2mm ID SXM4202 SXF4201 SXF42M3 SXF4202

SXFD4201 SXFD42M3 SXFD4202

ABS

SHUTOFF

SXF420212

Accessories

DESCRIPTION
SIXTUBE PRESSURE PLUGS

**PART NO.** TMPMWHT

**POLYPROPYLENE** 

STRAIGHT THRU

SXM170112

SXM170212

SXF170112

SXF170212

**POLYPROPYLENE** 

STRAIGHT THRI

SXM420112

SXM420212

SXF420112

# **MULTI-MOUNT SERIES CONNECTOR**



Multi-Mount couplings provide one easy-to-use coupling for connecting from three to five lines at once. Multi-mount couplings are available in a wide variety of materials including acetal and chrome-plated brass. Multi-mounts are keyed to prevent mismatched connections and can be configured with or without valves for maximum design flexibility.

FEATURES	BENEFITS
Three to five line connection	Fast, efficient operation
Separate flow paths	Pressure and/or vacuum in one coupling
Insert is keyed to body	Ensures correct line orientation
Acetal or chrome-plated brass	Maximize design flexibility

# Specifications • • •

#### PRESSURE:

Minimum: Vacuum

Maximum: Brass = 500 psi / no. of stations or Acetal = 240 psi / no. of stations

#### TEMPERATURE:

-40°F to 180°F (-40°C to 82°C)

#### MATERIALS:

Main components: Acetal or chrome-plated brass

Thumb latch: Stainless steel

Mounting Plate: Black anodized aluminum

Valves: Acetal

Valve spring: 316 stainless steel External springs and pin: Stainless steel

**0-ring:** Buna-N

#### **TUBING SIZES:**

1/8" to 3/8" ID. 3.2mm to 9.5mm ID

Multiple mount inserts are designed for mounting in 1/4"

Note: To order in chrome-plated brass, eliminate the "P" at the beginning of the part number.

# Multiple mount plates Step 4 Multiple mount coupling inserts Step 2 Standard panel mount coupling bodies with thumb latch for exterior stations Step 3 Multiple mount coupling bodies for interior stations

# Four-Step Multi-Mount Ordering Process

**STEP 1:** Determine flow size and termination type.


**STEP 2:** Order two "multiple mount plates" (one for the coupling bodies and one for the coupling inserts) with the appropriate number of connection stations for your specific application (available in 3, 4, or 5 connections, or construct your own). Refer to next page.

**STEP 3:** Order two "standard panel mount coupling bodies with thumb latches" for the exterior stations.

**STEP 4:** Order an appropriate number of "multiple mount coupling bodies without thumb latches" to accommodate the remaining interior stations.

**STEP 5:** Order an appropriate number of "multiple mount inserts" to connect with standard and multiple mount bodies.

# **MULTI-MOUNT SERIES DIMENSIONS**




# LM SERIES 1/4" Flow

MM SERIES 1/8" Flow

Multiple mount plates are keyed to prevent misconnections. >>

Multiple mount minimum coupling spacing: CPC offsets the last hole on the mounting plates so they cannot be miscon-



Multiple Mount Plates Information in RED indicates 1/8\* flow couplings. These are not compatible with 1/4\* flow couplings in non-shaded sections.

# Retaining Rings

NUMBER OF STATIONS	1/8" FLOW SIZE	LENGTH	1/4" FLOW SIZE	LENGTH
3	CP103	2.60	LCP103	3.00
4	CP104	3.35	LCP104	3.35
5	CP105	4.10	LCP105	4.90

COUPLING	PART NO.
SERIES	(REPLACEMENTS)
MC/PMC	101700
LC/PLC	278100

# **Coupling Bodies • ACETAL**



# **Coupling Inserts • ACETAL**

POLYTUBE FITTING,



TUBING/THREAD SIZE	METRIC EQ.	STRAIGHT THRU	SHUTOFF	A	В	D
5/32" OD, .10" ID	4.0mm OD, 2.5mm ID	PMM20025	PMMD20025	.58	1.86/1.96	
1/4" OD, .17" ID	6.4mm OD, 4.3mm ID	PMM2004	PMMD2004	.58	1.77/1.87	
1/4" OD, .17" ID	6.4mm OD, 4.3mm ID	PLM20004	PLMD20004	.72	2.10/2.23	
3/8" OD, .25" ID	9.5mm OD, 6.4mm ID	PLM20006	PLMD20006	.72	2.02/2.15	
1/8" ID	3.2mm ID	PMM2202	PMMD2202	.56	1.87/1.97	
3/16" ID	4.8mm ID	PMM2203	PMMD2203	.58	2.08/2.18	
1/4" ID	6.4mm ID	PMM2204	PMMD2204	.58	1.83/1.93	
1/4" ID	6.4mm ID	PLM22004	PLMD22004	.72	2.20/2.33	
5/16" ID	7.9mm ID	PLM22005	PLMD22005	.72	2.20/2.33	
3/8" ID	9.5mm ID	PLM22006	PLMD22006	.72	2.03/2.16	
1/4" OD17" ID	6.4mm OD, 4.3mm ID	PMM2104	PMMD2104	.53	1.41/1.51	.77

All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters. Couplings are pictured with valves unless otherwise noted. †NOTE: CPC's Ferruleless Polytube Fitting terminations do not require ferrules to achieve a secure connection, which makes them easier to use and reuse. PTF fittings are designed for semi-rigid tubing, i.e., polyethylene, nylon, polyurethane, etc.

# **BQ45GL SERIES CONNECTOR**



BQ45GL Series Bottle Caps include a quick disconnect coupling or FitQuik® connector integrated into a cap designed for 45GL bottles. The BQ45GL series feature a hydrophobic packaging vent, polyethylene gasket and a hose barb for mounting a dip tube. Adding quick disconnects to media bottles makes changing laboratory consumable products safer and easier.

FEATURES	BENEFITS
Factory installed quick disconnects	Ensures cap properly seals to bottle (even without sealing rings)
Hydrophobic packaging vent	Encourages proper fluid handling protocol allowing easy disconnections to replace or refill consumables in the correct location
Polyethylene foam gasket	Allows dispensing from container while preventing accidental spillage
Integral hose barb	Allows for mounting of 1/4" ID DIP tube

#### Specifications • • •



#### TEMPERATURE:

32°F to 120°F (0°C to 49°C)

#### **MATERIALS:**

Base Cap: High density polyethylene

Gasket: Low density polyethylene

Packaging Vent: Hydrophobic PTFE

**Quick Disconnects:** Polypropylene, with 316 stainless steel springs and EPDM o-rings

FitQuik Connector: Polypropylene


COLOR: Blue

#### **TUBING SIZES:**

FitQuik connector: 1/4" OD DIP tube interface: 1/4" ID

WARNING: Pressure, temperature, chemicals, and operating environment can affect the performance of couplings. It is the customer's responsibility to test the suitability of CPC products in their own application conditions.

# **BQ45GL DIMENSIONS**





# **Bottle Caps** • POLYETHYLENE




#### GL45 HDPE CAP WITH POLYETHYLENE FOAM SEALING GASKET, HYDROPHOBIC VENT AND INTEGRAL VALVED PMC12 BODY QUICK DISCONNECT INTERFACE



GL45 HDPE CAP WITH POLYETHYLENE FOAM SEALING GASKET, HYDROPHOBIC VENT AND INTEGRAL VALVED PMC12 INSERT QUICK DISCONNECT INTERFACE

G	L45 HDPE CAP WITH POLYETHYLENE FOAM
S	EALING GASKET, HYDROPHOBIC VENT AND
11	NTEGRAL POLYPROPYLENE ELBOW FITTING
(1	WITH 1/4" HOSE BARB)

	PART NO. BQ45GL1001	DIP TUBE SIZE 1/4" ID	CONNECTOR SIZE 1/4" OD	0.68"	B 0.83"	2.04"
	BQ45GL2001	1/4" ID	1/4" OD	1.08"	0.83"	2.04"
1	BQ45GL5001	1/4" ID	1/4" OD	0.63"	0.83"	2.04"



# HYBRID SERIES CONNECTOR

**Hybrid Connectors** integrate air lines, electrical connections and fluid lines into a single, compact connection point for remote tools on medical devices and equipment. Eliminating the potential for misconnecting tubing and wires, CPC's Hybrid Connector platform options reduce multiple connections and allow clinicians or patients to quickly and easily change or replace modular tools, umbilicals or hand pieces.

**Large Format Hybrid Connector** — Now available with four times the flow capacity, this connector features multiple dry-break couplings and easy-to-use twist latch for secure, repeatable connection of multiple lines.

#### **FEATURES BENEFITS** Single connection point for pneumatic, fluid, Increases safety and makes systems data and power connections between a device and more user friendly its tool/hand piece Non-spill shutoff valves with dry-break Eliminate fluid drips upon disconnect technology and minimize air inclusion Reduces multiple components Turnkey solution Minimizes engineering design and manufacturing effort Off-the-shelf assembly Custom options available Flexibility to add additional fluid connector or electrical contacts

#### Specifications • • •

#### PRESSURE:

Vacuum to 45 psi, 3.1 bar

#### TEMPERATURE:

32°F to 120°F (0°C to 49°C)

#### MATERIALS:

**Non-wetted coupler components:** Glass-filled polypropylene

**Wetted coupler components:** Glass-filled polypropylene

Springs: Stainless steel

**Electrical contacts:** Gold-plated copper alloy

**0-rings:** EPDM

Fluid tubing: Polyurethane (PU)
Umbilical tubing: PVC

COLOR: Gray

#### **TERMINATIONS:**

**Panel:** Two 1/8" (3.2mm) hose barbs on panel and four 20 ga wires with 12" (304.3mm) pigtail


**Umbillical Hose:** 8' (2.44M) umbilical with two 1/8" (3.2mm) ID PU tubing and four 20ga wires with 3" (76.2mm) pigtail

#### NUMBER OF FLUID LINES: 2


#### **ELECTRICAL SPECIFICATIONS**

Voltage: Up to 24VDC Current: Up to 3 Amps Wire gauge: 20AWG Number of contacts: 4

# HYBRID SERIES DIMENSIONS



#### Panel Dimensions




#### **Coupling Hose** • POLYPROPYLENE



#### **Coupling Panel** • POLYPROPYLENE





# MINI HYBRID SERIES CONNECTOR



**Mini Hybrid Connector** – A compact, single connection point with a smaller form factor that is ideal for use in low-pressure applications where ease of use and ergonomics are important. This connector makes it easy to connect power, signal and gases in one simple motion.

FEATURES	BENEFITS	
Easy-to-use breakaway design	Quick, automatic disconnect	
Small package design	Minimize device size	
Streamlined exterior	Clean, attractive appearance	
Proven design	Accelerate product development	

# Specifications • • •

#### PRESSURE:

Vacuum to 25 psi, 1.7 bar

#### TEMPERATURE:

40°F to 80°F (4°C to 27°C)

#### **MATERIALS:**

External coupler components: PET Air coupler components: Acetal Springs: Stainless steel

Panel Mount Nut: Zinc-plated steel

Electrical contacts: Gold-plated copper alloy

O-rings: Buna-N
Air line: PVC

COLOR: Gray

#### **EXTERNAL UMBILICAL HOSE:** 4 feet

(121 cm) PVC

WIRE: 22AWG PVC jacket

#### **TERMINATIONS:**

Panel: 1/16" hose barb and two 22AWG wires

with 11" pigtail


**Umbilical:** 4' (121 cm) umbilical with one 1/16" ID PVC tube and two 22AWG wires with 6" pigtail

Number of air lines: 1

#### **ELECTRICAL SPECIFICATIONS**

Voltage: Up to 24VDC Current: Up to 3 Amps Wire: 22AWG with PVC jacket Number of contacts: 2

# MINI HYBRID SERIES DIMENSIONS



#### Panel Dimensions




# **Coupling Hose** • PET AND ACETAL



#### Coupling Panel • PET AND ACETAL







# INTELLIGENT FLUID CONNECTIONS



# INTELLIGENT FLUID CONNECTIONS



RFID-based IdentiQuik® quick disconnect couplings redefine what's possible for controlling, protecting and streamlining fluid handling processes. A patented design embeds RFID-based electronic intelligence into a variety of our most popular and industry-leading coupling configurations, enabling the exchange of rewritable data once the coupling halves come within a few centimeters of one another. This intelligent technology allows you to protect your brand, minimize misconnections, prolong equipment life and increase efficiency.

IdentiQuik couplings are ideal for applications with a product involving: onboard electronics and parts or media that are consumable or expensive; the need to control the process to ensure the use of the right tool at the right time.

IdentiQuik couplings allow the robust transfer of fluid and information in a limitless number of applications across multiple industries. Examples of potential applications include:

- >> Connection validation
- >> Media validation
- >> Operational control settings
- >> Track and trace
- >> Usage control

# THOUSANDS OF CONFIGURATIONS TO MEET YOUR NEEDS

IdentiQuik couplings are available in a vast array of sizes and configurations that includes all our most popular product lines, opening endless new possibilities for your fluid handling processes.

Use this table to quickly find the right IdentiQuik coupling for your application. Can't find the ideal solution from our standard offering? Let our engineered solutions team design one for you.

# Flow Capacity



# Spillage Upon Disconnect



# OUR HIGHEST-FLOW-CAPACITY COUPLING

IdentiQuik® Universal Dispensing Coupler (UDC)



#### IdentiQuik Universal Dispensing Couplers (UDC) provide

a universal connection to a 38mm fitment neck, allowing you to make instant connections to bulk packaging systems. The automatic flush face valves minimize costly or dangerous spillage, and an ergonomic design and large, shrouded thumb latch pad make the coupling easy to grip and simple to operate.

#### APPLICATIONS FOR UDC COUPLINGS INCLUDE:

**UDC Coupling Readers** • POLYPROPYLENE

- Bag-in-box (BIB) packaging systems
- Rigid bulk packaging systems
- Flexible bulk packaging systems
- Many more

#### Specifications • • •



#### PRESSURE:

20" Hg Vacuum to 15 psi, 1 bar (body only)

TEMPERATURE: 32°F to 150°F (0°C to 65°C)

#### MATERIALS:

#### Main body components and valves:

Polypropylene and Acetal Thumb latch: Polypropylene **Springs:** 316 stainless steel

O-rings (food grade): Silicone, Buna-N and EPDM

**INTERFACE:** Multi-drop 2-ft. (60 cm) cable with RS-232 or USB connector

# OPERATING VOLTAGE:

RS-232: 8V − 24V • USB: 5V


**POWER CONSUMPTION:** 1.2W maximum

#### RF COMMUNICATION RANGE:

Approximately 1"

#### RFID TAG:

ISO 15693, 13.56MHz, 112 bytes programmable



# **UDC Caps** • ACETAL



#### POLYPROPYLENE

All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters Couplings are pictured with valves unless otherwise noted



ATION I 38MM	
-ON ALVED	

DESCRIPTION	METRIC EQ.	SHUTOFF	A	B
	38mm	iUDCDTSN0000	1.62	1.76
	Note: Snap-in closure	e designed to fit Scholle 900 bag neck	S.	
Without thread seal	38mm	iUDCDTCN3800	1.70	1.48
With thread seal	38mm	iUDCDTCS3800	1.70	1.48

# CHEMICAL-RESISTANT AND HIGHLY VERSATILE

IdentiQuik® PLC12 Coupling



IdentiQuik PLC12 Couplings are injection-molded from polypropylene thermoplastic and are resistant to many chemical solutions.

This popular coupling features one-hand connection and disconnection, and offers the flexibility of multiple configurations and terminations. Its polypropylene provides superior chemical resistance, and the EPDM o-ring adds expanded chemical compatibility for more demanding applications.

PLC12 couplings can be ordered in a vast range of sizes and configurations to meet the needs of your application.

#### **APPLICATIONS FOR PLC12 COUPLINGS INCLUDE:**

- Flavor concentrates
- Laser cooling
- Cosmetic ingredients
- Many more

#### Specifications • • •

#### PRESSURE:

Vacuum to 120 psi, 8.3 bar

#### TEMPERATURE:

32°F to 150°F (0°C to 65°C)

#### MATERIALS:

Main components and valves: Polypropylene

Thumb latch: Stainless steel Valve spring: 316 stainless steel External springs and pin: Stainless steel

O-rings: EPDM

COLOR: Almond

#### TUBING SIZE:

1/4" to 3/8", 6.4mm to 9.5mm

#### INTERFACE:

Multi-drop 2-ft. (60 cm) cable with RS-232 or USB connector

#### **OPERATING VOLTAGE:**

RS-232: 8V - 24V USB: 5V

#### POWER CONSUMPTION:

1.2W maximum

#### RF COMMUNICATION RANGE:

Approximately 1"

#### RFID TAG:

ISO 15693, 13.56MHz, 112 bytes programmable

# **Coupling Readers • POLYPROPYLENE**



# **Coupling Inserts • POLYPROPYLENE**



All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters. Couplings are pictured with valves

# FAST, SAFE AND LEAK-FREE FLUID CONNECTIONS

IdentiQuik® NS4 Non-Spill (NS) Coupling



IdentiQuik NS4 Non-Spill (NS) Couplings are designed with double-sided non-spill shutoff valves that enable drip-free disconnections and prevent

air entering the lines at connection, minimizing downtime and enhancing operator safety. These easy-to-use, lightweight couplings are chemically resistant and offer broad chemical compatibility

#### APPLICATIONS FOR NS4 COUPLINGS INCLUDE:

- Large-format printing
- Blood analysis
- Inkjet printing
- Many more

#### Specifications • • •

PRESSURE: Vacuum to 120 psi, 8.3 bar

**TEMPERATURE:** 32°F to 150°F (0°C to 65°C)

#### MATERIALS:

Main body components and valves:

Glass-filled polypropylene

Thumb latch: Glass-filled polypropylene

Springs: 316 stainless steel Valve spring: 316 stainless steel External spring: 316 stainless steel O-rings: EPDM (FKM optional) Soft-touch overmold: TPV\*

**COLOR:** Gray with dark gray overmold

LUBRICANTS: Krytox® PFPE (inert)

TUBING SIZE: 1/8" to 3/8", 3.2mm to 9.5mm

#### INTERFACE:

Multi-drop 2-ft. (60 cm) cable with RS-232 or USB connector

#### **OPERATING VOLTAGE:**

RS-232: 8V - 24V

USB: 5V

**POWER CONSUMPTION: 1.2W maximum** 

#### RF COMMUNICATION RANGE:

Approximately 1"

#### RFID TAG:

ISO 15693, 13.56MHz, 112 bytes programmable

#### **Coupling Readers • POLYPROPYLENE**



# **Coupling Inserts • POLYPROPYLENE**



IN-LINE HOSE BARB (RS-232 Interface) 1/8" ID 1/4" ID 3/8" ID

9.5mm

METRIC EQ. 3 2mm 6.4mm

iNS4DT2200200 iNS4DT2200400 iNS4DT2200600

For product dimensions visit cpcworldwide com

B

All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters. Couplings are pictured with valves unless otherwise noted

<sup>\*</sup>The overmold material is known as TPV (thermoplastic vulcanizate). This TPV is an alloy of polypropylene thermoplastic and fully vulcanized EPDM rubber. The material is typically resistant to water, acids and bases

# OUR MOST COMPACT COUPLING

IdentiQuik® Subminiature Coupling (SMC)



IdentiQuik Subminiature Couplings (SMC) offer a reliable and secure alternative to luer-type connections, combining twist-to-connect simplicity with the security of RFID technology. Free rotation eliminates kinking and accidental connection during use.

#### APPLICATIONS FOR SMC COUPLINGS INCLUDE:

Air lines

- Reagents
- Surgical irrigation lines
- Many more

#### Specifications • • •

PRESSURE: Vacuum to 100 psi, 6.9 bar

**TEMPERATURE:** 32°F to 150°F (0°C to 65°C)

#### MATERIALS:

Main components and valves: Acetal, Polypropylene

**Locking sleeves:** Acetal

Valves: Acetal

Valve spring: 316 stainless steel

**0-rings:** Buna-N with acetal or EPDM with

polypropylene

#### **COLOR (MAIN COMPONENTS):**

Natural white (acetal), almond (polypropylene)

#### **TUBING SIZE:**

1/16" to 1/8" ID, 1.6mm to 3.2mm ID NOTE: For 1/16" ID tubing, maximum tube OD is .150".

INTERFACE: Multi-drop 2-ft. (60 cm) cable with RS-232 or USB connector

#### **OPERATING VOLTAGE:**

RS-232: 8V − 24V • USB: 5V

**POWER CONSUMPTION:** 1.2W maximum

#### RF COMMUNICATION RANGE:

Approximately 1"

#### RFID TAG:

ISO 15693, 13.56MHz, 112 bytes programmable

# **Coupling Bodies** • ACETAL



IN-LINE **HOSE BARB** 

**TUBING SIZE** 1/16" ID 1/8" ID

METRIC EQ. 1.6mm ID 3.2mm ID

iSMFT0103 iSMFT0203

iSMFDT0103 iSMFDT0203

.76/.91 .76 .91/.91

#### • POLYPROPYLENE



TERMINATION IN-LINE HOSE BARB

TUBING SIZE 1/16" ID 1/8" ID

METRIC EQ. 1.6mm ID 3.2mm ID

STRAIGHT THRU iSMFT0100

B .76

# **Coupling Inserts** • ACETAL



IN-LINE HOSE BARB (Acetal Locking

TUBING SIZE 1/16" ID 1/8" ID

METRIC EQ. 1.6mm ID 3.2mm ID

STRAIGHT THRU iSMMT0103 iSMMT0203

.76 .76

#### POLYPROPYLENE

All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters Couplings are pictured with valves unless otherwise noted.



TERMINATION IN-LINE HOSE BARB (Acetal Locking

**TUBING SIZE** 1/16" ID 1/8" ID

METRIC EQ. 1.6mm ID 3.2mm ID

STRAIGHT THRU iSMMT0100 iSMMT0200

.76

.76

# UNMATCHED CUSTOMIZATION **CAPABILITIES**



#### For more than 35 years, CPC engineers have worked with

customers worldwide to design application-specific solutions that resolve their fluid handling-related challenges and improve product performance. We can do the same for you.

#### COMMON IDENTIQUIK® COUPLINGS CUSTOM APPLICATIONS INCLUDE:

- Weight, pressure, temperature and flow sensing
- Bulkhead mount readers
- Multiplexing technology



# CONNECTIONS for LIFE SCIENCES MEDICAL & BIOPROCESSING



# SMC SERIES



# Specifications • • •

#### PRESSURE:

Vacuum to 100 psi, 6.9 bar

#### TEMPERATURE:

-40°F to 250°F (-40°C to 121°C)

#### **MATERIALS:**

Main components: Polycarbonate, USP

Class VI, ADCF

Locking sleeves: Polycarbonate, USP

Class VI, ADCF

O-rings: Buna-N, USP Class V, ADCF

#### COLOR:

Main components: Purple tint

#### STERILIZATION:

Gamma: Up to 50 kGy irradiation

#### **TUBING SIZES:**


1/16" to 1/8" ID, 1.6mm to 3.2mm ID

**The SMC** is a twist-to-connect coupling that provides a reliable and more secure alternative to luer-type connections. It also allows for the tubing to rotate freely when connected. This important feature prevents both kinked tubing and accidental disconnection during use.

FEATURES	BENEFITS	
Twist to connect	Prevents accidental disconnects	
Free coupling rotation	Eliminates kinked tubing	
Safer; prevents misconnections	Does not mate to luers	

# These graphs are intended to give you a general idea of the performance capabilities of each product line. The shaded area of each graph represents the operating range of the product family, i.e., upper and lower values are shown. Therefore, depending on the exact coupling configurations selected, you can reasonably expect values to fall within the shaded area.

#### SMC Air Flow, 100 psig Inlet Pressure



**SMC Water Flow** 

.95 1.89 2.84 3.78 4.73 5.67 6.62

.75 1.00 1.25 1.50 1.75 Flow (gpm)

# Liquid Flow Rate Information for Couplings

The chart below shows the flow rate for CPC couplings. Each coupling was tested with water at 70°F (21°C). To determine flow rates for specific coupling configurations use the formula at the right.


# $C_{\nu}$ values for subminiature couplings

BODIES	INSERTS	SMM01	SMM02
	SMF01	.03	.03
	SMF02	.03	.19



- **Q** = Flow rate in gallons per
- C<sub>v</sub> = Average coefficient across various flow rates (see chart)
- **ΔP**= Pressure drop across coupling (psi)
- S = Specific gravity of liquid

# **SMC SERIES DIMENSIONS**



# **Coupling Bodies • POLYCARBONATE**



# **Coupling Inserts • POLYCARBONATE**



# **Coupling Set • POLYCARBONATE**





**The SRC connector** is a small bore connector that eliminates the potential for misconnections with luer fittings. The intuitive design is simple for end users to operate and provides an audible "click" for added assurance and a secure, leak-free connection. This new connector features a unique elastomeric seal with a smooth, gap-free flow path that provides better flow characteristics and eliminates stagnant flow areas.

FEATURES	BENEFITS		
Prevents misconnections	Won't connect to IV lines and other non-compatible applications		
Reliable elastomeric seals	Leak-free connections		
Audible "click" and positive feel	Operators know when a secure, reliable connection is made		
Allows tubing rotation	Prevents accidental disconnect or kinked tubing		

# Specifications • • •

#### PRESSURE:

Vacuum to 25 psi, 1.7 bar

#### TEMPERATURE:

32°F to 110°F (0°C to 43°C) (rated for shipment uncoupled up to 160°F)

#### MATERIALS:

Main components: Polypropylene with overmolded TPE\* pads, USP Class VI, ADCF

#### **INTEGRAL SEALS:**

Overmolded TPE\*, USP Class VI, ADCF

#### COLOR:

White with blue pads


#### **STERILIZATION:**

EtO, e-beam or gamma

#### **TUBING SIZES:**


1/8" to 3/16" ID, 3.2mm to 4.8mm ID

\*The overmold material TPV (themoplastic vulcanizate) is used with the SRC polypropylene couplings. TPV is an alloy of polypropylene thermoplastic and fully vulcanized EPDM rubber. TPV is typically resistant to water, acids and bases. The overmold material TPE (thermoplastic elastomer) is used with the NS4 ABS couplings. TPE is a blend of additives and copolymers in a special formulation that forms extremely durable bonds to the ABS substrate, while offering the traditional properties of soft-touch overmold.



This graph is intended to give you a general idea of the

# **SRC SERIES DIMENSIONS**



#### **Panel Dimensions**

	PANEL	PANEL THICKNESS	PANEL	PANEL NUT
	OPENING	MAXMIN.	NUT HEX	THREAD
COUPLING INSERTS	7/16	.2103	1/2	7/16-24UNS

# A = Height/Diameter B = Total Length

# **Coupling Bodies**

POLYPROPYLENE

Add "CL" to the end of the plastic SRC part numbers below for cleanroom manufactured versions.



TERMINATION IN-LINE HOSE BARB **TUBING SIZE** 1/8" ID 3/16" ID METRIC EQ. 3.2mm ID 4.8mm ID

STRAIGHT TH SRC1702 SRC1703 .62 .62 B 1.18 1.18

#### **Coupling Inserts**

#### POLYPROPYLENE



TERMINATION IN-LINE HOSE BARB TUBING SIZE 1/8" ID 3/16" ID METRIC EQ. 3.2mm ID 4.8mm ID STRAIGHT THRU SRC2202 SRC2203 .38 .38 B 1.19 1.19

#### • CHROME-PLATED BRASS



TERMINATION PANEL MOUNT HOSE BARB **TUBING SIZE** 1/8" ID 3/16" ID

METRIC EQ. 3.2mm ID 4.8mm ID STRAIGHT THRU SRC4202CB SRC4203CB .57 .57

1.66 1.66


#### Accessories



*PART NO.* SRC3L **DESCRIPTION**Tethered cap and plug for SRC body and insert

MATERIALS
White medical-grade TPE

# MPC SERIES CONNECTOR



MPC Series couplings add ease of use and security to critical fluid handling applications. Choose from a full line of connectors and configurations, including pressure sealing caps and plugs, in sizes to fit 1/8" to 3/8" tubing. MPC couplings offer optional locking sleeves to further guard against accidental disconnects. In addition, coupling halves can be rotated when connected to reduce tube kinks.

FEATURES	BENEFITS	
Ergonomic thumb latch	Easy to operate — even with gloved hands	
USP Class VI materials	Meet biocompatibility requirements	
Sterilizable by autoclave, Et0, e-beam, or gamma	Reusable, yet economical enough to allow disposability	
Parting line-free hose barb	Eliminates potential leak path	
ADCF-free materials	Meet BSE/TSE requirements	

# **Product Dimensions**



#### Specifications • • •

#### OII3

#### PRESSURE:

Vacuum to 60 psi, 4.1 bar

#### **TEMPERATURE:**

#### ABS:

-40°F to 160°F (-40°C to 71°C)

#### Polycarbonate:

-40°F to 250°F (-40°C to 121°C)

#### Polysulfone:

-40°F to 300°F (-40°C to 149°C)

#### **MATERIALS:**

#### Main components:

ABS (white), USP Class VI Polycarbonate (purple tint), USP Class VI, ADCF Polysulfone (amber tint), USP Class VI, ADCF

#### Locking sleeves:

Polysulfone (white), USP Class VI, ADCF

#### Thumb Latches:

Polycarbonate (white), USP Class VI, ADCF Polysulfone (amber tint), USP Class VI, ADCF

**O-rings:** Silicone (clear), platinum-cured, USP Class VI, ADCF and Buna-N (black), USP Class V

#### STERILIZATION:

Gamma: Up to 50 kGy irradiation

#### Autoclave:

**Polycarbonate:** Up to 250°F (121°C), 30 minutes, up to 10 repetitions. Sterilize uncoupled only.


**Polysulfone:** Up to 270°F (132°C) for 60 minutes, up to 25 repetitions. Sterilize uncoupled only.

#### **TUBING SIZES:**

1/8" to 3/8" ID (3.2mm to 9.5mm)

**WARNING:** Pressure, temperature, chemicals, and operating environment can affect the performance of couplings. It is the customer's responsibility to test the suitability of CPC's products in their own application conditions.


This graph is intended to give you a general idea of the performance capabilities of each product line. The shaded area of the graph represents the operating range of the product family, i.e. upper and lower values are shown. Therefore, depending on the exact coupling configurations selected, you can reasonably expect values to fall within the shaded area.



# MPC SERIES CONNECTOR DIMENSIONS

# **Coupling Bodies**

ABS



1

1	
	TERMINATION
200	IN-LINE
100	HOSE BARB

TUBING	METRIC EQ.	FLOW	STRAIGHT THRU
1/4" ID	6.4mm ID	.21"	MPC17004T
3/8" ID	9.5mm ID	.29"	MPC17006T

B
1.30 (33.0
1.30 (33.0

1.30 (33.0)

1.30 (33.0)

1.30 (33.0)

1.30 (33.0)


.60 (15.2)

.60 (15.2)

.60 (15.2)

.60 (15.2)

#### POLYCARBONATE




#### • POLYSULFONE


600	TERMINATION IN-LINE HOSE BARB	<b>TUBING</b> 1/8" ID 1/4" ID 3/8" ID	<b>METRIC EQ.</b> 3.2mm ID 6.4mm ID 9.5mm ID	FLOW .09" .21" .29"	STRAIGHT THRU MPC17002T39 MPC17004T39 MPC17006T39	.96 (24.4) .96 (24.4) .96 (24.4)	B 1.10 (27.9) 1.10 (27.9) 1.10 (27.9)
OFF STREET	IN-LINE HOSE BARB WITH LOCK	1/8" ID 1/4" ID 3/8" ID	3.2mm ID 6.4mm ID 9.5mm ID	.09" .21" .29"	MPCK17002T39 MPCK17004T39 MPCK17006T39	.99 (25.2) .99 (25.2) .99 (25.2)	1.10 (27.9) 1.30 (33.0) 1.30 (33.0)

# Coupling Inserts


• ABS



#### POLYCARBONATE



#### POLYSULFONE



# MPC SERIES CONNECTOR (cont.)



#### **Mating Parts**

SEALING CAP	
MPC32003	
MPC32039	

MPC30003M

MPC30039M

O-RING Silicone Seal USP Class VI Silicone Seal USP Class VI

DESCRIPTION

Leash plug for MPC body

Leash cap for MPC insert

MPCK32003

MPCK32039

PART NO. MPC30L

# SEALING CAP W/LOCK

MATERIAL Polycarbonate Polysulfone

.96 (24.4) .99 (25.2) A .75 (19.1)

.75 (19.1)

1.30 (33.0) 1.30 (33.0)

B

B 1.24 (31.5) 1.24 (31.5)

**MATERIALS** 

MATERIAL

Polycarbonate

Polysulfone

Soft, flexible, medical-grade PVC Soft, flexible, medical-grade PVC

# Liquid Flow Rate Information for Couplings

The chart below shows the flow rate for CPC couplings. Each coupling was tested with water at 70°F (21°C). To determine flow rates for specific coupling configurations use the formula at the right.

#### C<sub>v</sub> VALUES FOR MPC SERIES COUPLINGS

BODIES INSERTS	MPC22002T03	MPC22004T03	MPC22006T03
MPC17002T03	0.11	-	0.21
MPC17004T03	-	2.8	2.8
MPC17006T03	0.14	2.8	5.5

- ${f Q}~=~{\sf Flow}$  rate in gallons per
- C<sub>v</sub> = Average coefficient across various flow rates (see chart)
- ΔP = Pressure drop across coupling (psi)
- S = Specific gravity of liquid

MPC Series nates with Back-to-Sanitary Series (pages 112-115).



# BACK-TO-BACK SERIES CONNECTOR



MPC/MPX Back-to-Back Adapters give end users the flexibility of connecting single-use systems that feature identical coupling connections at the end of their tubing. Combining both MPC and MPX couplings provides a reducing option for users who need to transition between tubing diameters ranging from 1/8" to 1/2".

FEATURES	BENEFITS
Compatible with MPC and MPX Series inserts	Easy conversion to industry standard connections or single-use systems
Tubing reduction option	Allows easy transition between multiple size tubing from 1/8" to 1/2" ID
Ergonomic thumb latches	Easy to operate - even with gloved hands
ADCF-free materials	Meet BSE/TSE requirements

# Back-To-Back Insert Adapters

#### POLYSULFONE



# Back-To-Back Body Adapters

#### • POLYCARBONATE

450	PART NO.	ТҮРЕ	A	B
35.460	MPC17C1703	MPC to MPC	.96 (24.5)	1.81 (46.0)
The state of	MPX17X1703	MPX to MPX	1.28 (32.5)	2.44 (62.0)
Charles of the Control of the Contro	MPC17X1703	MPC to MPX	1.28 (32.5)	2.13 (54.1)

#### • POLYSULFONE



All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters

#### Specifications • • •



#### PRESSURE:

Vacuum to 60 psi, 4.1 bar

#### **TEMPERATURE:**

#### Polycarbonate:

-40°F to 250°F (-40°C to 121°C)

#### Polysulfone:

-40°F to 300°F (-40°C to 149°C)

#### **MATERIALS:**

#### Main Components:

Polycarbonate (purple tint), USP Class VI, ADCF Polysulfone (amber tint), USP Class VI, ADCF

#### Thumb Latches:

Polycarbonate (white), USP Class VI, ADCF Polysulfone (amber tint), USP Class VI, ADCF

#### 0-rings:

Silicone (clear), platinum-cured, USP Class VI, ADCF

#### **STERILIZATION:**

Gamma: Up to 50 kGy irradiation


Polycarbonate: Up to 250°F (121°C) for 30 minutes, up to 10 repetitions.


Sterilize uncoupled only.

Polysulfone: Up to 270°F (132°C) for 60 minutes, up to 25 repetitions. Sterilize uncoupled only.

WARNING: Pressure, temperature, chemicals, and operating environment can affect the performance of couplings. It is the customer's responsibility to test the suitability of CPC's products in their own application conditions.







A = Height/Diameter

B = Total Length



MPX Series couplings add ease of use and security to your most critical fluid handling applications. Choose from a full line of connectors and configurations, including pressure sealing caps and plugs in sizes to fit 3/8" and 1/2" tubing. MPX couplings offer optional locking sleeves to further guard against accidental disconnects. In addition, coupling halves can be rotated when connected reducing tube kinks.

FEATURES	BENEFITS
Ergonomic thumb latch	Easy to operate — even with gloved hands
USP Class VI materials	Meet biocompatibility requirements
Sterilizable by autoclave, Et0, e-beam, or gamma	Reusable, yet economical enough to allow disposability
Parting line-free hose barb	Eliminates potential leak path
ADCF-free materials	Meet BSE/TSE requirements

# Specifications • • •

#### PRESSURE:

Vacuum to 60 psi, 4.1 bar

#### **TEMPERATURE:**

#### Polycarbonate:

-40°F to 250°F (-40°C to 121°C)

#### Polysulfone:

-40°F to 300°F (-40°C to 149°C)

#### **MATERIALS:**

#### Main components:

Polycarbonate (purple tint), USP Class VI, ADCF Polysulfone (amber tint), USP Class VI, ADCF

#### Locking sleeves:

Polysulfone (white), USP Class VI, ADCF

#### Thumb Latches:

Polycarbonate (white), USP Class VI, ADCF Polysulfone (amber tint), USP Class VI, ADCF

**0-rings:** Silicone (clear), platinum-cured, USP Class VI, ADCF

#### STERILIZATION:

Gamma: Up to 50 kGy irradiation

#### Autoclave:


Polycarbonate: Up to 250°F (121°C), 30 minutes, up to 10 repetitions. Sterilize uncoupled

Polysulfone: Up to 270°F (132°C) for 60 minutes, up to 25 repetitions. Sterilize uncoupled


#### **TUBING SIZES:** 3/8" to 1/2" ID (9.5mm to 12.7mm)

WARNING: Pressure, temperature, chemicals, and operating environment can affect the performance of couplings. It is the customer's responsibility to test the suitability of CPC's products in their own application conditions.

#### > This graph is intended to give you a general idea of the performance capabilities of each product line. The shaded area of the graph represents the operating range of the product family, i.e. upper and lower values are shown. Therefore, depending on the exact coupling configurations selected, you can reasonably expect values to fall within the shaded area.



# MPX SERIES CONNECTOR DIMENSIONS



A = Height/Diameter B = Total Length

#### **Coupling Bodies**

#### • POLYCARBONATE



TERMINATION IN-LINE HOSE BARB	<b>TUBING</b> 1/2" ID	METRIC EQ. 12.7mm ID	<b>FLOW</b> .50"	STRAIGHT THRU MPX17803	A 1.28 (32.5)	B 1.96 (49.8)
IN-LINE HOSE BARB WITH LOCK	1/2" ID	12.7mm ID	.50"	MPXK17803	1.28 (32.5)	1.96 (49.8)

#### POLYSULFONE



ERMINATION V-LINE OSE BARB	<b>TUBING</b> 1/2" ID	METRIC EQ. 12.7mm ID	<b>FLOW</b> .50"	STRAIGHT THRU MPX17839	A 1.28 (32.5)	B 1.96 (49.8)
N-LINE HOSE BARB VITH LOCK	1/2" ID	12.7mm ID	.50"	MPXK17839	1.28 (32.5)	1.96 (49.8)

# **Coupling Inserts**

#### • POLYCARBONATE



-	TERMINATION	TUBING	METRIC EQ.	FLOW	STRAIGHT THRU	O-RING	A	B
	IN-LINE	3/8" ID	9.5mm ID	.38"	MPX22603M	Silicone Seal USP Class VI	.85 (21.6)	1.90 (48.3)
	HOSE BARB	1/2" ID	12.7mm ID	.50"	MPX22803M	Silicone Seal USP Class VI	.85 (21.6)	1.90 (48.3)

Silicone Seal USP Class VI

#### POLYSULFONE



(pages 112-115).

TERMINATION	TUBING	METRIC EQ.	<b>FLOW</b>	STRAIGHT THRU	O-RING	A	B
IN-LINE	3/8" ID	9.5mm ID	.38"	MPX22639M	Silicone Seal USP Class VI	.85 (21.6)	1.90 (48.3)
HOSE BARB	1/2" ID	12.7mm ID	.50"	MPX22839M	Silicone Seal USP Class VI	.85 (21.6)	1.90 (48.3)

# **Mating Parts**

MPX30039M



SEALING CAP MPX32003	SEALING CAP W/LOCK MPXK32003	<i>MATERIAL</i> Polycarbonate	A 1.28 (32.5)	B 1.67 (42.4)
MPX32039	MPXK32039	Polysulfone	1.28 (32.5)	1.67 (42.4)
SEALING PLUG MPX30003M	<b>O-RING</b> Silicone Seal USP Class VI	<b>MATERIAL</b> Polycarbonate	A 1.10 (27.9)	B 1.66 (42.2)

All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters.

Polysulfone

1.10 (27.9)

1.66 (42.2)

# SANITARY SERIES CONNECTOR



# Specifications • • •

#### PRESSURE:

Vacuum to 60 psi, 4.1 bar

#### TEMPERATURE:

-40°F to 300°F (-40°C to 149°C)

#### **MATERIALS:**

#### Main components:

Polysulfone (amber tint), USP Class VI, ADCF

**0-rings:** Silicone (clear), platinum-cured,

USP Class VI, ADCF

#### STERILIZATION:

Gamma: Up to 50 kGy irradiation

**Autoclave:** Up to 270°F (132°C) for 60 minutes, up to 25 repetitions. Sterilize uncoupled only.

#### **TERMINATION SIZES:**

3/4", 1" and 1-1/2" sanitary

WARNING: Pressure, temperature, chemicals, and operating environment can affect the performance of couplings. It is the customer's responsibility to test the suitability of CPC's products in their own application conditions.

**Sanitary couplings** attach directly to 3/4", 1" and 1-1/2" sanitary terminations to provide greater flexibility for integrating components into single-use or hybrid (single-use to stainless) process systems. Standard bag systems with quick couplings can be easily connected to equipment with sanitary terminations, while single-use cartridge filters can be converted to incorporate quick couplings for greater system modularity.

FEATURES

BENEFITS

3/4", 1" and 1-1/2" sanitary terminations

Install to equipment with sanitary gaskets and sanitary clamps

Compatible with MPC and MPX
Series couplings

Quick and easy connections to industry standard plastic couplings on single-use bags and tube sets

Integral coupling adaptor


Provides flexibility to easily convert sanitary terminations on filter cartridge or equipment

ADCF-free materials

Meet BSE/TSE requirements

Note:
Mates with MPC and
MPX bodies, inserts
and sealing caps and
plugs and Back-to-Back
Adapters
(pages 106-111).

# SANITARY SERIES CONNECTOR DIMENSIONS



# **Coupling Bodies • POLYSULFONE**

		PART NO.	SIZE	A	B	C
	200	MPC3301239	3/4"	.98 (24.9)	1.40 (35.6)	1.0 (25.4)
	1000	MPX3301239	3/4"	1.28 (32.5)	1.70 (43.2)	1.0 (25.4)
The state of the s		MPC3301639	1"	1.50 (38.1)	1.40 (35.6)	1.50 (38.1)

# **Coupling Inserts • POLYSULFONE**

<b>PART NO.</b> MPC44012T39M MPC44024T39M	<b>SIZE</b> 3/4" 1-1/2"	<b>O-RING</b> Silicone Seal USP Class VI Silicone Seal USP Class VI	B 1.40 (35.6) 1.40 (35.6)
MPX44012T39M MPX44024T39M	3/4" 1-1/2"	Silicone Seal USP Class VI Silicone Seal USP Class VI	1.71 (43.4) 1.71 (43.4)

All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters.

Note: Sanitary couplings are compatible with both stainless steel and plastic clamps. Clamps and gaskets are referenced for illustration and are not available through CPC.





# SANIQUIK<sup>™</sup> SERIES CONNECTOR



# Specifications • • •

#### PRESSURE:

Vacuum to 60 psi, 4.1 bar

#### TEMPERATURE:

-40°F to 300°F (-40°C to 149°C)

#### **MATERIALS:**

Main component: 316L stainless steel
O-rings: Silicone (clear), platinum-cured,
USP Class VI, ADCF


#### **STERILIZATION:** Autoclave

#### **TERMINATION SIZES:**

3/4" and 1-1/2" sanitary


**WARNING:** Pressure, temperature, chemicals, and operating environment can affect the performance of couplings. It is the customer's responsibility to test the suitability of CPC's products in their own application conditions.

**CPC's SaniQuik<sup>TM</sup> connection** answers the question of how to integrate single-use components with your existing stainless processing equipment. This integral sanitary termination attaches to hard-plumbed systems with tri-clover clamps. Once attached it permits quick and easy connection to single-use bag systems, manifolds or tube sets with CPC disposable coupling bodies. SaniQuik connections reduce sanitary gasket replacement, enabling cost-effective media transfer solutions for feeding, harvesting or sampling applications.



FEATURES	BENEFITS	
3/4" and 1-1/2" sanitary standard terminations	Connect to hard plumbed systems with sanitary gaskets and sanitary clamps	
Compatible with MPC and MPX Series couplings	Quick and easy connections to industry standard plastic couplings on single-use bag and tube sets	
Integral coupling adaptor	Disconnecting coupling reduces sanitary gasket replacement	
ADCF-free materials	Meet BSE/TSE requirements	

# SANIQUIK<sup>™</sup> SERIES CONNECTOR DIMENSIONS





#### **Connections** • 316L STAINLESS



DESCRIPTION
SILICONE SEAL
USP CLASS VI

PART NO.
SQCC221212M
SQCC222424M
SQCX221212M
SQCX222416M
SQCX222424M

MATING COUPLING MPC Series MPC Series MPX Series MPX Series MPX Series SANITARY SIZE 3/4" 1-1/2" 3/4" 1-1/2"

.98" (24.9) 1.98" (50.3) .98" (24.9) 1.98" (50.3) B 1.39" (35.3) 1.50" (38.1) 1.43" (39.1) 1.50" (38.1)

1.50" (38.1)

**Accesories** • SILICONE (CLEAR)

DESCRIPTION
PLATINUM-CURED
USP CLASS VI
REPLACEMENT
SEALS

**PART NO.** 2260100 2260200

MATING SANIQUIK

SQCC221212M, SQCC222424M

SQCX221212M, SQCX222416M, SQCX222424M

All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters.


BORE

3/4"

1-1/2"

3/4"

1-1/2"



# MPU SERIES CONNECTOR

# Specifications • • •

#### PRESSURE:

Vacuum to 35 psi, 2.4 bar

#### TEMPERATURE:

-40°F to 300°F (-40°C to 149°C)

#### **MATERIALS:**

#### Main components:

Polysulfone (amber tint), USP Class VI, ADCF

**O-rings:** Silicone (clear), platinum-cured, USP Class VI, ADCF

#### **STERILIZATION:**

**Gamma:** Up to 50 kGy irradiation


Autoclave: Up to 270°F (132°C) for 60 minutes, up to 25 repetitions. Sterilize uncoupled only.


#### TUBING SIZES: 3/4" ID (19.0 mm)

WARNING: Pressure, temperature, chemicals, and operating environment can affect the performance of couplings. It is the customer's responsibility to test the suitability of CPC's products in their own application conditions.

> These graphs are intended to give you a general idea of the performance capabilities of each product line. The shaded area of each graph represents the operating

range of the product family, i.e., upper and lower values are shown. Therefore, depending on the exact coupling configurations selected, you can reasonably expect values to fall within the shaded area. MPU Water Flow






# The MPU's twist-to-connect design features an easy-to-use locking mechanism that guards against accidental disconnects and provides a reliable, secure connection. A 3/4" hose barb provides smooth, rapid media transfer.

FEATURES	BENEFITS	
3/4" hose barb	Facilitates rapid fill and empty of bioprocessing bags	
Locking feature	Guards against accidental disconnects	
Sharp barb end	Minimizes fluid turbulence and dead space	
Shrouded, leak-free seal & smooth, internal flow path	Protect valuable fluids and eliminate potential to contaminate fluid path	
ADCF-free materials	Meet BSE/TSE requirements	

# MPU SERIES CONNECTOR DIMENSIONS



A = Height/Diameter **B** = Total Length

# **Coupling Bodies • POLYSULFONE**



IN-LINE HOSE BARB

3/4" ID

METRIC EQ. 19.0mm ID

STRAIGHT THRU MPU171239

1.75 (44.5) 2.37 (60.2)

# **Coupling Inserts • POLYSULFONE**



**TERMINATION** IN-LINE HOSE BARB

3/4" ID

METRIC EQ. 19.0mm ID

STRAIGHT THRU MPU221239M

O-RING

Silicone Seal USP Class VI

1.56 (39.6) 2.88 (73.2)

#### **Accessories** • POLYSULFONE



SEALING CAP MPU32039		MATERIAL Polysulfone USP Class VI	A 1.75 (44.5)	B .79 (20.1)
SEALING PLUG	<b>O-RING</b>	MATERIAL Polysulfone USP Class VI	A	B
MPU30039M	Silicone Seal		1.56 (39.6)	1.38 (35.1)

All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters.

# **HFC39 SERIES CONNECTOR**



# Specifications • • •

#### PRESSURE:

Vacuum to 125 psi, 8.6 bar

#### TEMPERATURE:

-40°F to 280°F (-40°C to 138°C)

#### **MATERIALS:**

#### Main components:

Polysulfone (amber tint), USP Class VI, ADCF

**0-rings:** Silicone (clear), platinum-cured,

USP Class VI, ADCF

Springs: 316L stainless steel

#### **STERILIZATION:**

**Gamma:** Up to 50 kGy irradiation. Sterilize coupled or uncoupled.

**Autoclave:** Up to 270°F (132°C) for 60 minutes. Up to 25 repetitions for uncoupled units and up to one repetition for coupled units.

#### **TERMINATION SIZES:**


1/4", 3/8" and 1/2" ID hose barb (6.4mm, 9.5mm and 12.7mm)


# **HFC39 Series sterile disconnect couplings** prevent external organisms from entering into the media flow path upon disconnection. Automatic shutoff valves close off the flow path aseptically protecting valuable media while also eliminating the need for pinch clamps and tube welders. The easy-to-use thumb latch

FEATURES	BENEFITS	
Simple one-step disconnection	Maintains media sterility in each half by preventing external organisms from entering the flow path	
Automatic shutoff valves	Stop flow and eliminate need for pinch clamps	
CPC Click	Audible confirmation of secure connection	
Lightweight	Easy integration with single-use assemblies	
ADCF-free materials	Meet BSE/TSE requirements	


design provides a secure, leak-free connection and enables one-hand disconnects.

These graphs are intended to give you a general idea of the performance capabilities of each product line. The shaded area of the graph represents the operating range of the product family, i.e. upper and lower values are shown. Therefore, depending on the exact coupling configurations selected, you can reasonably expect values to fall within the shaded area.





# HFC39 SERIES DIMENSIONS



# **Coupling Bodies • POLYSULFONE**



#### **Coupling Inserts • POLYSULFONE**



# **Mating Parts**



All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters. Couplings are pictured with valves unless otherwise noted.

For validation reports, visit cpcworldwide.com/ bio and for extractables data contact info@cpcworldwide.com

# ASEPTIQUIK® S SERIES CONNECTORS



# AseptiQuik® S Small Format 1/8" and 1/4" Connectors

provide quick and easy sterile connections for small-flow applications, even in nonsterile environments. The "FLIP-CLICK-PULL" design of AseptiQuik S enables users to easily transfer small volumes of media with less risk of operator error. The connector's genderless and robust design provides reliable performance without the need for clamps, fixtures or tube welders. Biopharmaceutical manufacturers can now make 1/8" and 1/4" sterile connections with the quality and market availability they expect from the leader in single-use connection technology.

FEATURES	BENEFITS	
Genderless design	Eases single-use systems specifications with one part number for both halves	
FLIP-CLICK-PULL	Intuitive three-step connection process reduces risk of operator error	
Integrated pull tab covers	Provide pre-assembly protection and ensure simultaneous removal of both membranes	
Robust construction	Repeatable and reliable performance with no additional hardware required	
CPC Click	Audible confirmation of assembly	

# Specifications • • •

#### PRESSURE:

Up to 60 psi, 4.1 bar

#### TEMPERATURE:

39°F to 104°F (4°C to 40°C)

#### STERILIZATION:

Gamma: Up to 50kGy irradiation

AutoClave High Temp (HT) Version: Up to 266°F (130°C) for 60 minutes

#### **TERMINATION SIZES:**

1/8", 1/4" and 3/8" ID hose barb (3.2mm, 6.4mm and 9.5mm), 3/4" sanitary and MPC insert

#### MATERIALS:

#### Main Components:

Polycarbonate (white), USP Class VI, ADCF

#### Pull Tabs/Caps:

Polycarbonate (blue, standard version), USP Class VI, ADCF

Polycarbonate (white, HT version), USP Class VI. ADCF

#### Seals:


Silicone (clear), platinum-cured, USP Class VI, ADCF

#### Membrane:


Polyethylene (standard version), USP Class VI, ADCF


Hydrophobic polyethersulfone (HT version), USP Class VI, ADCF, PTFE strip sticker

> This graph is intended to give you a general idea of the performance capabilities of each product line. The shaded area of the graph represents the operating range of the product family, i.e. upper and lower values are shown. Therefore, depending on the exact coupling configurations selected, you can reasonably expect values to fall within the shaded area.



# ASEPTIQUIK® S SERIES DIMENSIONS







# **Coupling Bodies**

#### • POLYCARBONATE with Blue Pull Tabs



TERMINATION	
1/8" HOSE BARB	
1/4" HOSE BARB	
3/8" HOSE BARB	
3/4" SANITARY	
MPC INSERT	

PART NO.	
AQS17002	
AQS17004	
AQS17006	
AQS33012	
AQS17MPC	

В	G
1.30	0.50
1.45	0.65
1.45	0.65
1.60	0.80
1.49	0.69

#### POLYCARBONATE HT with White Pull Tabs



TERMINATION 1/8" HOSE BARB 1/4" HOSE BARB 3/8" HOSE BARB 3/4" SANITARY MPC INSERT

<b>PART NO.</b> AQS17002HT AQS17004HT AQS17006HT	2.25 2.25 2.25 2.25
AQS33012HT AQS17MPCHT	2.25 2.25 2.25

2.25 2.25 2.25 2.25

2.25 2.25

В	C
1.30	0.5
1.45	0.6
1.45	0.6
1.60	0.8
1.49	0.6

# FLIP-CLICK-PULL Assembly Procedure





Unsnap and flip down the protective pull tab covers on each connector half.





Align the connector halves with the pull tabs hanging down. Slide the two halves together.

until you hear an audible "CPC Click."

"CLICK"



To complete the connection, simply snap the pull tabs together by pushing on the CPC logos and Independently squeeze each side of the connector pull the membranes from the connector.

**120** cpcworldwide.com • 800-444-2474

# ASEPTIQUIK® G SERIES CONNECTORS



Genderless AseptiQuik® G Connectors enable quick and easy sterile connections, even in non-sterile environments. The easy-to-use genderless design simplifies system integration and minimizes the risk of operator error. The connectors' robust construction provides enhanced user confidence and reliable performance without the need for clamps, fixtures or tube welders. Biopharmaceutical manufacturers benefit from a full range of interchangeable 1/4" to 3/4" flow solutions with the quality and market availability they expect from the leader in single-use connection technology.

FEATURES	BENEFITS
Genderless design	Eases integration of single-use systems with one part number for both halves
Robust construction	Repeatable and reliable performance with no additional hardware required
FLIP-CLICK-PULL	Innovative three-step connection process reduces risk of operator error
Integrated pull tab covers	Pull tabs act as protective cover reducing part complexity and ensure simultaneous removal of both membranes
CPC Click	Audible confirmation of assembly

# Specifications • • •

#### PRESSURE:

Up to 60 psi, 4.1 bar Up to 75 psi, 5.1 bar for 48 hours

#### TEMPERATURE:

34°F to 104°F (1°C to 40°C)

#### STERILIZATION:

Gamma: Up to 50kGy irradiation

Autoclave High Temp (HT) Version: Up to 266°F (130°C) for 60 minutes

#### **TERMINATION SIZES:**

1/4", 3/8", 1/2", 3/4" ID hose barb (6.4mm, 9.5mm, 12.7mm, 19.0mm) and 3/4" sanitary

#### MATERIALS:

#### Main Components:

Polycarbonate (white), USP Class VI, ADCF

#### Pull Tabs/Caps:

Polycarbonate (blue, standard version), USP Class VI, ADCF Polycarbonate (white, HT version), USP Class VI, ADCF

Silicone (clear), platinum-cured, USP Class VI, ADCF

#### Membrane:

Polyethylene (standard version), USP Class VI, ADCF Hydrophobic polyethersulfone (HT version), USP Class VI, ADCF, PTFE strip sticker


tubing sizes from


This graph is intended to give you a general idea of the performance capabilities of each product line. The shaded area of the graph represents the operating range of the product family, i.e. upper and lower values are shown. Therefore, depending on the exact coupling configurations selected, you can reasonably expect values to fall within the shaded area.

# AseptiQuik G Flow Flow (lpm) 227.1 454.2 681.3 908.4 1135.5 120.0 180.0 Flow (gpm)

# ASEPTIQUIK® G SERIES DIMENSIONS







#### **Coupling Bodies**

• POLYCARBONATE with Blue Pull Tabs



TERMINATION	PART NO.	A	В	C
1/4" HOSE BARB	AQG17004	2.62	1.62	0.82
3/8" HOSE BARB	AQG17006	2.62	1.62	0.82
1/2" HOSE BARB	AQG17008	2.62	1.95	1.15
3/4" HOSE BARB	AQG17012	2.62	2.36	1.56
3/4" SANITARY	AQG33012	2.62	1.66	0.86

#### • POLYCARBONATE HT with White Pull Tabs



TERMINATION	PART NO.	A	В	C
1/4" HOSE BARB	AQG17004HT	2.62	1.62	0.82
3/8" HOSE BARB	AQG17006HT	2.62	1.62	0.82
1/2" HOSE BARB	AQG17008HT	2.62	1.95	1.15
3/4" HOSE BARB	AQG17012HT	2.62	2.36	1.56
3/4" SANITARY	AQG33012HT	2.62	1.66	0.86

Validation reports are available at pcworldwide.com/bio

# FLIP-CLICK-PULL Assembly Procedure



Unsnap and flip down the protective pull tab

covers on each connector half.





Align the connector halves with the pull tabs hanging down. Slide the two halves together. Independently squeeze each side of the connector until you hear an audible "CPC Click."

To complete the connection, simply snap the pull tabs together by pushing on the CPC logos and pull the membranes from the connector.

**122** cpcworldwide.com • 800-444-2474

# ASEPTIQUIK® C SERIES CONNECTORS



AseptiQuik® C Connectors provide quick and easy sterile connections, even in non-sterile environments. AseptiQuik's "CLICK-PULL-TWIST" design enables users to transfer media easily with less risk of operator error. The connector's robust design provides reliable performance without the need for clamps, fixtures or tube welders. Biopharmaceutical manufacturers can make sterile connections with the quality and market availability they expect from the leader in single-use connection technology.

#### **FEATURES BENEFITS** CLICK-PULL-TWIST design Intuitive three-step connection process reduces risk of operator error Membrane pull tabs Ensure simultaneous and secure removal of both membranes Repeatable and reliable performance with Robust construction no additional hardware required Integrated lock ring Secures final connection preventing disassembly Audible confirmation of completed CPC Click assembly steps

Open access through multiple supply

chain partners

# Specifications • • •

#### PRESSURE:

Up to 60 psi, 4.1 bar

#### TEMPERATURE:

39°F to 104°F (4°C to 40°C)

#### TYPICAL FLOW RATE:

Cv = 14.4 max

#### **STERILIZATION:**

Gamma: Up to 50kGy irradiation

Autoclave High Temp (HT) Version:

Up to 266°F (130°C) for 30 minutes

#### **TERMINATION SIZES:**

3/8" (9.5mm) and 1/2" (12.7mm) ID hose barb, and 3/4" sanitary

#### **MATERIALS:**

#### Main Components:

Polycarbonate (white), USP Class VI, ADCF

#### Lock Ring:

Polycarbonate (blue), USP Class VI, ADCF

#### **Pull Tabs:**

Polycarbonate (blue, standard version), USP Class VI, ADCF

Polycarbonate (white, HT version), USP Class VI, ADCF

#### Caps:

Polypropylene (clear), USP Class VI, ADCF


#### Seals:

Silicone (clear), platinum-cured, USP Class VI. ADCF

#### Membrane:


Polyethylene (standard version), USP Class VI. ADCF Hydrophobic polyethersulfone (HT version), USP Class VI, ADCF

This graph is intended to give you a general idea of the performance capabilities of each product line. The shaded area of the graph represents the operating range of the product family, i.e. upper and lower values are shown. Therefore, depending on the exact coupling configurations selected, you can reasonably expect values to fall within the shaded area.




# ASEPTIQUIK® C SERIES DIMENSIONS





#### Coupling Bodies

POLYCARBONATE with Blue Pull Tabs



#### POLYCARBONATE HT with White Pull Tabs




#### Coupling Inserts

• POLYCARBONATE with Blue Pull Tabs



#### POLYCARBONATE HT with White Pull Tabs



# CLICK-PULL-TWIST Assembly Procedure



Align male and female couplings, push together until "CPC Click" confirmation

Slight rotation of lock ring may be required for proper alignment prior to connection



Snap membrane pull tabs together and pull from connector.



Twist the blue lock ring clock-wise until "CPC Click" and rib to arrow alignment confirms secure final connection.

Market availability

# ASEPTIQUIK® DC SERIES CONNECTORS



AseptiQuik® DC Connectors are the first all-in-one single-use connection technology to offer both a sterile connect and sterile disconnect. With the AseptiQuik DC Connector, manufacturers can make a quick and easy sterile connection and disconnection, even in non-sterile environments.

AseptiQuik DC's intuitive "CLICK-PULL-TWIST" design enables users to transfer media easily with less risk of operator error. After transfer is complete, the connector features a simple one-step disconnection that maintains media sterility by preventing external organisms from entering into the media flow path. The connector's robust design and automatic shutoff valves provide reliable performance without the need for sanitary clamps, fixtures or tube welders.

FEATURES	BENEFITS
CLICK-PULL-TWIST design	Intuitive three-step connection process reduces risk of operator error
Simple one-step disconnection	Maintains media sterility in each half by preventing external organisms from entering the flow path
Membrane pull tabs	Ensure simultaneous and secure removal of both membranes
Robust construction	Repeatable and reliable performance with no additional hardware required
CPC Click	Audible confirmation of completed assembly steps
Market availability	Open access through multiple supply chain partners

# Specifications • • •

#### PRESSURE:

Up to 60 psi, 4.1 bar

#### TEMPERATURE:

39°F to 104°F (4°C to 40°C)

#### STERILIZATION:

**Gamma:** Up to 50kGy irradiation Autoclave High Temp (HT) Version: Up to 266°F (130°C) for 30 minutes

#### **TERMINATION SIZES:**

1/4". 3/8" and 1/2" ID hose barb (6.4mm, 9.5mm and 12.7mm)

#### **MATERIALS:**

Main Components: Polycarbonate (white). USP Class VI. ADCF

Lock Ring: Polycarbonate (blue), USP Class VI, ADCF

#### Pull Tabs:

Polycarbonate (blue, standard version), USP Class VI, ADCF Polycarbonate (white, HT version), USP Class VI. ADCF

Caps: Polypropylene (clear), USP Class VI, ADCF

Seals: Silicone (clear), platinum-cured, USP Class VI, ADCF


#### Membrane:

Polyethylene (standard version), USP Class VI, ADCF Hydrophobic polyethersulfone (HT version), USP Class VI, ADCF

**Springs:** 316L stainless steel


AseptiQuik halves: **IQCDC** inserts mate with standard AQC bodies **AQCDC** bodies mate wth standard AQC

This graph is intended to give you a general idea of the performance capabilities of each product line. The shaded area of the graph represents the operating range of the product family, i.e. upper and lower values are shown. Therefore, depending on the exact coupling configurations selected, you can reasonably expect values to fall within the shaded area.



# ASEPTIQUIK® DC SERIES DIMENSIONS





#### **Coupling Bodies** • POLYCARBONATE



**TERMINATION** 1/4" HOSE BARB 3/8" HOSE BARB 1/2" HOSE BARB

AQCDC17004 AQCDC17006 AQCDC17008

1/4" 3/8" 2.36 (59.9) 2.36 (59.9) 2.36 (59.9)

5.33 (135.4) 5.33 (135.4) 5.62 (142.7)

0.60 (15.2) 0.60 (15.2) 0.89 (22.6)

#### **Coupling Inserts • POLYCARBONATE**



1/4" HOSE BARB 3/8" HOSE BARB 1/2" HOSE BARB PART NO. AQCDC22004 AQCDC22006 AQCDC22008 **FLOW** 1/4" 3/8" 3/8"

A 2.50 (63.5) 2.50 (63.5) 2.50 (63.5)

5.41 (137.4) 5.41 (137.4) 5.70 (144.8)

0.60 (15.2) 0.60 (15.2) 0.89 (22.6)

# CLICK-PULL-TWIST Assembly Procedure



Align male and female couplings, push together until audible "CPC Click" confirmation.

Slight rotation of lock ring may be required for proper alignment prior to connection.



Snap membrane pull tabs together and pull from connector.



Twist the blue lock ring clock-wise until "CPC Click" and rib to arrow alignment confirms secure final connection.

# To Disconect



When fluid transfer is finished, press the thumb latch down to complete the sterile disconnection. Both halves will remain sterile.

**126** cpcworldwide.com • 800-444-2474

# ASEPTIQUIK® X SERIES CONNECTORS



AseptiQuik® X Large Format 1" Connectors provide quick and easy sterile connections for high flow applications, even in non-sterile environments. AseptiQuik X's "TWIST-PULL-TWIST" design enables users to quickly transfer large volumes of media easily with less risk of operator error. The connector's robust design provides reliable performance without the need for clamps, fixtures or tube welders. Biopharmaceutical manufacturers can now make 1" sterile connections with the quality and market availability they expect from the leader in single-use connection technology.

FEATURES	BENEFITS
TWIST-PULL-TWIST Design	Intuitive three-step connection process reduces risk of operator error
Membrane pull tabs	Ensure simultaneous and secure removal of both membranes
Robust construction	Repeatable and reliable performance with no additional hardware required
Integrated lock ring	Secures final connection preventing disassembly
CPC Click	Audible confirmation of completed assembly steps

#### Specifications • • •

#### PRESSURE:

Up to 60 psi, 4.1 bar

#### TEMPERATURE:

39°F to 104°F (4°C to 40°C)

#### STERILIZATION:

**Gamma:** Up to 50kGy irradiation Autoclave High Temp (HT) Version: Up to 266°F (130°C) for 30 minutes

#### **TERMINATION SIZES:**

3/4" and 1" ID hose barb (19.0mm and 25.4mm). 1-1/2" sanitary

#### **MATERIALS:**

#### Main Components:

PVDF (white), USP Class VI, ADCF

#### Lock Ring:

PVDF (blue), USP Class VI, ADCF

#### **Pull Tabs:**

Polycarbonate (blue, standard version), USP Class VI, ADCF

Polycarbonate (white, HT version), USP Class VI, **ADCF** 

#### Caps:

Polypropylene (clear), USP Class VI, ADCF


#### Seals:

Silicone (clear), platinum-cured, USP Class VI, ADCF

#### Membrane:


Polyethylene (standard version), USP Class VI, ADCF Hydrophobic polyethersulfone (HT version), USP Class VI, ADCF, PTFE strip sticker


This graph is intended to give you a general idea of the performance capabilities of each product line. The shaded area of the graph represents the operating range of the product family, i.e. upper and lower values are shown. Therefore, depending on the exact coupling configurations selected, you can reasonably expect values to fall within the shaded area.



# ASEPTIQUIK® X SERIES DIMENSIONS









#### Coupling Bodies

POLYCARBONATE with Blue Pull Tabs




POLYCARBONATE HT with White Pull Tabs




#### Coupling Inserts


• POLYCARBONATE with Blue Pull Tabs



#### POLYCARBONATE HT with White Pull Tabs



# TWIST-PULL-TWIST Assembly Procedure




indicator of the white body. Insert

the two halves together.

Twist the blue lock ring clockwise (approximately 75°) until audible "CPC Click" is heard. Alignment of the small blue lock ring rib indicator with the white body's rib indicator confirms final rotation is complete.







Twist the blue lock ring clockwise until the final audible "CPC Click" is heard (90° from the initial starting point). Alignment of the long blue lock ring rib indicator with the white body's rib indicator confirms initial rotation is complete.


# STEAM-THRU® SERIES CONNECTOR



Steam-Thru Connections allow a quick and easy sterile connection between stainless steel biopharmaceutical processing equipment and disposable bag and tube assemblies. The single-use design saves time and money by eliminating unnecessary cleaning procedures and reducing validation burden associated with reusable components.

FEATURES	BENEFITS	
Innovative three-port design	Allows a true steam-through SIP process which eliminates "dead legs" and the need for laminar flow hoods	
Patented valve design	Allows sterile connection and disconnection and permits high media flow rate	
Thumb latch/Tear-away sleeve	Secures valve position, provides visual indicator of process stage	
3/4" and 1-1/2" Terminations	Easily connects to process equipment	
ADCF-free materials	Meet BSE/TSE requirements	

#### These graphs are intended to give you a general idea of the performance capabilities of each product line. The shaded area of the graph represents the operating range of the product family, i.e. upper and lower values are shown. Therefore, depending on the exact coupling configurations selected, you can reasonably expect values to fall within the shaded area



# Specifications • • •

#### PRESSURE:

#### Steam position:

Up to 30 psi, 2.1 bar (Steam-Thru) Up to 35 psi, 2.4 bar (Steam-Thru II)

Flow position: Vacuum to 20 psi, 1.4 bar

#### **TEMPERATURE:**

#### Steam position:

Up to 266°F (130°C) for 60 minutes (Steam-Thru) Up to 275°F (135°C) for 60 minutes (Steam-Thru II)

Flow position: 39°F to 104°F (4°C to 40°C)

#### **MATERIALS:**

Connection: Polysulfone (amber tint), USP Class VI. ADCF

**0-rings:** Silicone (clear), platinum-cured, USP Class VI, ADCF

**Tear-away sleeve:** Polyethylene or polycarbonate (Steam-Thru only), USP Class VI

#### TYPICAL FLOW RATE:

 $C_{v} = 4.2 - 4.6$  (Steam-Thru)  $C_{v}^{v} = 5.2 - 8.0$  (Steam-Thru II)

#### **STERILIZATION:**

Gamma: Up to 50 kGy irradiation

Autoclave: Up to 265°F (129°C) for 60 minutes, up to two cycles (applies only to part numbers STC1700500-STC1700800 and STC2020000-STC2021000)


#### SIP process:

Up to 266°F (130°C) for 60 minutes (Steam-Thru) Up to 275°F (135°C) for 60 minutes (Steam-Thru

#### **TERMINATION SIZES:**

3/8" (9.5mm) to 1/2" (12.7mm) ID hose barb (Steam-Thru)


3/8" (9.5mm) to 1/2" (12.7mm) ID hose barb and 3/4" sanitary (Steam-Thru II)



# STEAM-THRU SERIES DIMENSIONS

# Steam-Thru® Configurations

Steam-Thru Connection's patented three-port design allows steam to pass directly through the lower ports to "steam on" to stainless equipment. After the SIP cycle is completed, the connector's valve is actuated, creating a sterile flow path to single-use systems.




# **Coupling Bodies • POLYSULFONE**



# Steam-Thru II Configurations

Steam-Thru II Connections offer the flexibility of "steam on" and "steam off" functionality. The innovative design allows the valve to be returned to the steam position enabling a second SIP cycle following media transfer. The "steam off" disconnection of single-use systems minimizes cross-contamination risks associated with reusable components.



# **Coupling Bodies • POLYSULFONE**



DESCRIPTION WITH **POLYCARBONATE** SLEEVE

PART NO.	TERMINATIONS	A	B	C	D	<b>E</b>
STC2020000	3/4" x 3/4" sanitary x 1/2" HB	1.42 (36.1)	6.84 (173.7)	5.93 (150.6)	2.40 (61.0)	.88 (22.4)
STC2020100	3/4" x 3/4" sanitary x 3/8" HB	1.42 (36.1)	6.76 (171.7)	5.93 (150.6)	2.40 (61.0)	.80 (20.3)
STC2020200	3/4" x 1-1/2" sanitary x 1/2" HB	1.42 (36.1)	6.84 (173.7)	5.93 (150.6)	2.40 (61.0)	.88 (22.4)
STC2020300	3/4" x 1-1/2" sanitary x 3/8" HB	1.42 (36.1)	6.76 (171.7)	5.93 (150.6)	2.40 (61.0)	.80 (20.3)
STC2020900	3/4" x 3/4" sanitary x 3/4" sanitary	1.42 (36.1)	6.60 (167.6)	5.93 (150.6)	2.40 (61.0)	.62 (15.7)
STC2021000	3/4" x 1-1/2" sanitary x 3/4" sanitary	1.42 (36.1)	6.60 (167.6)	5.93 (150.6)	2.40 (61.0)	.62 (15.7)

All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters.

# ASEPTIQUIK® STC SERIES CONNECTOR



AseptiQuik® STC Connectors integrate the AseptiQuik® sterile

connector and the Steam-Thru® II SIP connector, giving manufacturers even greater flexibility for hybrid stainless steel and single-use processing. The AseptiQuik STC connector combines either an AseptiQuik G or AseptiQuik C sterile connector with a Steam-Thru II connection that can be mounted directly to a stainless steel vessel via a sanitary termination. AseptiQuik STC connectors utilizing an AseptiQuik G sterile connector offer a genderless design, simplifying system integration at the aseptic connector end. AseptiQuik STC connectors utilizing an AseptiQuik C sterile connector allow system designers to incorporate a keyed female-to-male connection at the aseptic connector end.

The union of the two connectors into a single unit through a sanitary clamp allows an AseptiQuik sterile connection to be steamed on to stainless equipment via SIP. After the SIP cycle, a wide range of single-use systems can be connected. The SIP process can be done in advance allowing a quick and easy sterile connection to the AseptiQuik half without having to wait 30-60 minutes for SIP prior to harvest.

FEATURES	BENEFITS
AseptiQuik design	Innovative three-step connection process reduces risk of operator error
Genderless design	Eases integration of single-use systems with universal mating between AseptiQuik G Connectors
Innovative three-port steam design	Allows a true steam-through SIP process which eliminates "dead legs"
Robust construction	Repeatable and reliable performance with no additional hardware required
Patented steam valve design	Allows sterile connection and disconnection to stainless equipment and permits a high media flow rate.
Sanitary interface between the two connectors	More secure connection than tubing with cable ties
CPC Click	Audible confirmation of assembly steps
Market availability	Open access through multiple supply chain partners

#### Specifications • • •

#### PRESSURE:

**Steam position:** Up to 35 psi, 2.4 bar **Flow position:** Up to 20 psi, 1.4 bar

#### TEMPERATURE:

#### Steam position:

Up to 275°F (135°C) for 60 minutes

Flow position: 39°F to 104°F (4°C to 40°C)

#### **STERILIZATION:**

**Gamma:** Up to 50kGy irradiation

AutoClave: High Temp (HT) version: Up to 266°F (130°C) for 30 minutes (AQSTC) Up to 266°F (130°C) for 60 minutes (AQGSTC)

**Note**: A slight clockwise rotation of the clamp nut may be needed after autoclave.

#### **MATERIALS:**

#### **Main Components:**

**AseptiQuik** - Polycarbonate (white), USP Class VI, ADCF

Steam-Thru II - Polysulfone (amber tint), USP Class VI, ADCF

**Seals:** Silicone (clear), platinum-cured, USP Class VI, ADCF

#### Pull Tabs:

Polycarbonate (blue, standard version), USP Class VI, ADCF Polycarbonate (white, HT version), USP Class VI. ADCF

Caps: Polypropylene (clear), USP Class VI, ADCF

#### Membrane:


Polyethylene (standard version), USP Class VI, ADCF

Hydrophobic polyethersulfone (HT versions), USP Class VI, ADCF, PTFE strip sticker

Clamp: Nylon 66 (white), USP Class VI



# ASEPTIQUIK® STC SERIES DIMENSIONS



A = Height/Diameter

Mates with gendered male and female AseptiQuik Series.

#### **Coupling Bodies**

POLYCARBONATE with Blue Pull Tabs



TERMINATION 3/4" X 3/4" SANITARY STEAM-THRU II WITH ASEPTIQUIK BODY

3/4" X 1-1/2" SANITARY STEAM-THRU II WITH ASEPTIQUIK BODY (PICTURED)

PART NO.	A	B	©
AQSTC2330900	2.78 (70.6 mm)	9.25 (235 mm)	8.42 (213.9 mm)
AQSTC2331000	2.78 (70.6 mm)	9.25 (235 mm)	8.42 (213.9 mm)

POLYCARBONATE HT with White Pull Tabs




#### 3/4" X 3/4" SANITARY STEAM-THRU II WITH ASEPTIQUIK BODY (PICTURED)

3/4" X 1-1/2" SANITARY STEAM-THRU II WITH ASEPTIQUIK BODY

PART NO. AQSTC2330900HT	2.78 (70.6 mm)	B 9.25 (235 mm)	8.42 (213.9 mm)
AQSTC2331000HT	2.78 (70.6 mm)	9.25 (235 mm)	8.42 (213.9 mm)

# ASEPTIQUIK® G STC SERIES DIMENSIONS



A = Height/Diameter

B = Total LengthC = Actuated Length

Mates with genderless AseptiQuik G

#### Coupling Bodies

• POLYCARBONATE with Blue Pull Tabs



TERMINATION 3/4" X 3/4" SANITARY STEAM-THRU II WITH GENDERLESS ASEPTIQUIK G BODY

3/4" X 1-1/2" SANITARY STEAM-THRU II WITH GENDERLESS ASEPTIQUIK G BODY

PART NO. AQGSTC2330900	<b>A</b> 3.05	B 8.26	7.36
AQGSTC2331000	3.05	8.26	7.36

POLYCARBONATE HT with White Pull Tabs




TERMINATION 3/4" X 3/4" SANITARY STEAM-THRU II WITH GENDERLESS ASEPTIQUIK G BODY


3/4" X 1-1/2" SANITARY STEAM-THRU II WITH GENDERLESS ASEPTIQUIK G BODY


440.020000000	A 3.05	B 8.26	<b>C</b> 7.36
	3.05	8.26	7.36

# ASEPTIQUIK® STC SERIES ASSEMBLY PROCEDURE


# CLICK-PULL-TWIST Assembly Procedure







Il in a steam trap creating a "steam on" sterile connection.





Twist the blue lock ring clockwise until audible "CPC Click." Alignment of the lock ring rib with the body's arrow indicator confirms final connection.



Once the AseptiQuik assembly is complete and the steam trap has been closed, simply press the thumb latch to allow the valve to be pushed into the flow position. The "CPC Click" confirms the transition is complete.



Once the valve is locked into the flow position, you are ready for sterile fluid transfer to or from the process equipment.





After sterile media transfer is complete, simply press the thumb latch and pull the valve back into the steam position. The "CPC Click" confirms the transition is complete.




# ASEPTIQUIK® G STC SERIES ASSEMBLY PROCEDURE

# FLIP-CLICK-PULL Assembly Procedure





















# CONNECTIONS for CHEMICAL HANDLING APPLICATIONS





# CQH/CQV SERIES CONNECTOR





The CQH and CQV Series are designed for high-purity use and feature all plastic construction. Molded, virgin materials and lubricant-free design enable use in the most demanding applications. Broad chemical compatibility and cleanroom manufacturing make them ideal for use in critical wet processes.

FEATURES	BENEFITS		
100% metal free	No risk of metal contaminants or corrosion		
High flow valve design	High flow in a compact package		
Disconnect under pressure	Speeds servicing and reduces risk of injury		
Polypropylene and PVDF	Broad chemical compatibility		









#### Specifications • • •

#### **CQH PRESSURE:**

0 - 80 psig (0 - 5.5 bar)

#### **CQH TEMPERATURE:**

32°F to 225°F (0°C to 107°C)

#### **CQH MATERIALS:**

Main components: Natural, virgin polypropylene Valve o-rings: FKM (black)

External insert o-ring: Simriz® FFKM perfluoro

Valve (wetted) and thumb latch spring: PEEK® Flare nuts: PVDF

Lubricants: None used

**CQV PRESSURE:** 0 - 80 psig (0 - 5.5 bar)

#### **CQV TEMPERATURE:**

0°F to 225°F (-18°C to 107°C)

#### **CQV MATERIALS:**

Main components: Natural, virgin PVDF Valve o-rings: Chemraz FFKM perfluoro (white) External insert o-ring: Chemraz FFKM perfluoro

Valve (wetted) and thumb latch spring: PEEK

Flare nuts: PVDF Lubricants: None used

#### SPILLAGE/AIR INCLUSION:

~1.5 cc (ml)/disconnect (reconnect)

#### **TUBING SIZES:**

1/4" to 1/2" flare; 3/8" to 3/4" NPT


NOTE: CQH & CQV are also available with optional PPS springs and Chemraz FFKM or EPDM seals.



These graphs are intended to give you a general idea of the performance capabilities of each product line. The shaded area of each graph represents the operating range of the product family, i.e., upper and lower values are shown. Therefore, depending on the exact coupling configurations selected, you can reasonably expect values to fall within the shaded area.



# CQH/CQV SERIES CONNECTOR DIMENSIONS



# **CQH Coupling Bodies • POLYPROPLYENE**

(Page	TERMINATION IN-LINE	TUBING/THREAD SIZE 1/4" OD 3/8" OD 1/2" OD 3/8" Taper	TERMINATION TYPE Flare compression Flare compression Flare compression Male NPT	SHUTOFF CQHD06100104 CQHD06100106 CQHD06100108 CQHD06100206	1.44 1.44 1.44 1.44	B 3.32 3.38 3.38 2.73	1.00
93	TERMINATION PANEL MOUNT	1/4" OD 3/8" OD 1/2" OD 3/8" Taper	Flare compression Flare compression Flare compression Male NPT	CQHD06110104 CQHD06110106 CQHD06110108 CQHD06110206	1.82 1.82 1.82 1.82	3.32 3.38 3.37 2.73	1.00


# **CQH Coupling Inserts • POLYPROPLYENE**

TERMINATION IN-LINE	TUBING/THREAD SIZE 1/4" OD 3/8" OD 1/2" OD 3/8" Taper	TERMINATION TYPE Flare compression Flare compression Flare compression Male NPT	SHUTOFF CQHD06200104 CQHD06200106 CQHD06200108 CQHD06200206	1.00 1.00 1.00 1.16	B 2.49 2.55 2.56 1.90	1.00	
	3/8" Taper	Male NP1	CQHD06200206	1.16	1.90	1.00	
	3/4" Taper	Male NPT	CQHD06200212	1.16	2.26	1.00	

# **CQV Coupling Bodies • PVDF**

We do	TERMINATION IN-LINE	TUBING/THREAD SIZE 1/4" OD 3/8" OD 1/2" OD 3/8" Taper	TERMINATION TYPE Flare compression Flare compression Flare compression Male NPT	SHUTOFF CQVD06100104 CQVD06100106 CQVD06100108 CQVD06100206	1.44 1.44 1.44 1.44	B 3.31 3.38 3.37 2.70	1.00
W.	TERMINATION PANEL MOUNT	1/4" OD 3/8" OD 1/2" OD 3/8" Taper	Flare compression Flare compression Flare compression Male NPT	CQVD06110104 CQVD06110106 CQVD06110108 CQVD06110206	1.82 1.82 1.82 1.82	3.24 3.38 3.37 2.70	1.00

# **CQV** Coupling Inserts • PVDF



All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters.

# **CQG SERIES CONNECTOR**




The CQG06 Series high purity couplings feature our patented pressurebalanced, non-spill design. Molded virgin polypropylene and a 100% springless flow path provide broad chemical resistance and exceptionally high flow capacity, allowing instant disconnects (and reconnects), even under pressure.

FEATURES	BENEFITS
Non-spill design	Ultimate protection from chemicals and fumes
Pressure-balanced	Failsafe disconnect, even under pressure; easy to reconnect at high pressure
Springless flow path design	Eliminates source of metallic contaminant

These graphs are intended to give you a general idea of the performance capabilities of each product line. The shaded area of each graph represents the operating range of the product family, i.e., upper and lower values are shown. Therefore, depending on the exact coupling configurations selected, you can reasonably expect values to fall within the

**CQGD06** Pressure Range



# Specifications • • •

#### PRESSURE:

0 - 80 psig (0 - 5.5 bar)

#### TEMPERATURE:

32°F to 150°F (0°C to 66°C)

#### MATERIALS:

Main components:

Natural, virgin polypropylene

Seals: FKM (Simriz® FFKM optional)

Springs (non-wetted): Hastelloy® C

Flare nuts: PVDF Lubricants: None used **Dual Containment Nut** and Panel Mount Fitting: Virgin, natural polypropylene

Flare nuts: PVDF Panel nuts: HDPE Panel mount o-rings: FKM

#### SPILLAGE/AIR INCLUSION:

<0.1 cc (ml)/disconnect (reconnect)


PANEL MOUNT: Optional adapter kit


#### **TUBING SIZES:**

3/8" to 3/4" flare; 3/8" to 3/4" NPT

WARNING: Pressure, temperature, chemicals, and operating environment can affect the performance of couplings. It is the customer's responsibility to test the suitability of CPC's products in their own application conditions.

# **CQG SERIES CONNECTOR DIMENSIONS**





#### **Coupling Bodies • POLYPROPLYENE**

TERMINATION IN-LINE	TUBING/THREAD SIZE 3/8" OD 1/2" OD 3/4" OD 1/2" Taper	TERMINATION TYPE Flare compression Flare compression Flare compression Male NPT	SHUTOFF CQGD06100106 CQGD06100108 CQGD06100112 CQGD06100208	1.96 1.96 1.96 1.96	B 4.21 4.42 4.49 3.81	1.00	
------------------------	---	---	---	------------------------------	-----------------------------------	------	--

#### **Coupling Inserts • POLYPROPLYENE**

TERMINATION IN-LINE	TUBING/THREAD SIZE 3/8" OD 1/2" OD 3/4" OD 3/8" Taper	TERMINATION TYPE Flare compression Flare compression Flare compression Male NPT	SHUTOFF CQGD06200106 CQGD06200108 CQGD06200112 CQGD06200206‡	1.62 1.62 1.62 1.62	B 4.35 4.45 4.63 3.81	1.00	
	3/8" Taper 1/2" Taper	Male NPT Male NPT	CQGD06200206‡ CQGD06200208	1.62 1.62	3.81 3.95	1.00 1.00	
	3/4" Taper	Male NPT	CQGD06200212‡	1.62	4.15	1.00	

#### **Panel Mount Adaptor** • POLYETHYLENE

DESCRIPTION	PART NO.
Panel mount adapter kit (fits both bodies & inserts)	CQG06PMKIT01
Tanci mount adapter hit (its both bodies & inserts)	OQUOUI MINITOI

All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters

#### ChemQuik® Dual Containment Flare Nuts • POLYPROPYLENE



The ChemQuik® Dual Containment System is an easy way to "double contain" critical chemical lines, protecting plant and personnel in case a primary process line ruptures or "sweats". The system provides a protective secondary line to catch any fluid and routes it to a safe location.

These fittings work with any ChemQuik coupling with fine thread flare terminations or a common flare style fitting. In addition, the panel mount version can be mounted into a pump cabinet or other panel mount connection point. The primary line can then be routed from the coupling directly to a pump, connected to a ChemQuik coupling at the panel mount fitting.

CQPMDCNUT0612

CQPMDCNUT0812

CQPMDCNUT1216

The "weep port" serves to vent the area between the primary and secondary lines so that pressure cannot build up in case of a primary line rupture. The leaking fluid can then be routed to a containment vessel or to a leak detector.

FLARE NUT	DESCRIPTION	
CQDCNUT0408	Dual containment nut, 1/4" OD process line x	
	1/2" OD secondary containment line.	
CQDCNUT0612	Dual containment nut, 3/8" OD process line x	
	3/4" OD secondary containment line.	
CQDCNUT0812	Dual containment nut, 1/2" OD process line x	
	3/4" OD secondary containment line.	
CQDCNUT1216	Dual containment nut, 3/4" OD process line x	
	1" OD secondary containment line.	

#### FLARE FITTING CQPMDCNUT0408

Panel mount dual containment nut, 1/4" OD Process line x 1/2" OD secondary containment line. Panel mount dual containment nut. 3/8" OD Process line x 3/4" OD secondary containment line. Panel mount dual containment nut, 1/2" OD Process line x 3/4" OD secondary containment line. Panel mount dual containment nut, 3/4" OD Process line x 1" OD secondary containment line.

**140** cpcworldwide.com • 800-444-2474

shaded area.


# DRUMQUIK PRO / DRUMQUIK PUR



DrumQuik® PRO and DrumQuik PUR dispensing systems are used with common pumps to create a closed transfer system for both general and high purity chemicals. These easy-to-use products meet the strict safety requirements of chemical handling applications as well as the demands of harsh operating environments. DrumQuik dispensing systems minimize chemical exposure and facilitate the delivery of chemicals from bulk transfer containers to end use applications.

For general chemical dispensing, DrumQuik PRO offers a rugged and costeffective solution. DrumQuik PUR is a 100% metal-free, cleanroom packaged, PVDF product designed specifically for applications demanding the highest level of purity. All materials are compliant with FDA standards for food and consumables applications.

> This graph characterizes the flow of the DrumQuik PRO with the coupler vent port open to the atmosphere and a pump drawing liquid out.



# DrumQuik Dispensing Systems for Safer Chemical Management

Whether you're dispensing edible oils, condiments, detergents, specialty automotive fluids, or industrial chemicals, acids, caustics or solvents, DrumQuik PRO and the high purity DrumQuik PUR can make your process cleaner, faster, safer and smarter.

- Provide the critical connection between your chemical container and the end use process
- Create a closed transfer system to protect workers and keep a cleaner work environment

DrumQuik solutions are available for a variety of container types providing interchangeable, reusable connections for:

- Plastic or steel drums
- Chemical packs
- IBCs or totes

Note: For high purity applications DrumQuik PUR provides the assurance of compatibility with chemicals used in semiconductor and flat panel production, pharmaceutical and photovoltaic or solar cell manufacturing.

The DrumQuik PRO dispensing system has and IBCs for use with 49 CFR 173 and German BAM 11027.

> Made from virgin materials compliant with FDA and ROHS standards.



# DrumQuik® PRO and DrumQuik PUR are changing the world of chemical management with their unique closed-system design. Innovative DrumQuik systems help companies reduce costs, increase worker safety and protect the environment. From basic chemical handling to specialty high purity and pharmaceutical applications, DrumQuik dispensing systems make chemical management cleaner, faster, safer and smarter.

### Specifications • • •



#### PRESSURE:

**Uncoupled:** Vacuum to 45 psig. 3.1 bar Coupled: Vacuum to drum rating

#### TEMPERATURE:

**Drum insert:** -20°F to 120°F (-20°C to 49°C) Coupler: HDPE, 32°F to 120°F (0°C to 49°C) PVDF, 0°F to 120°F (-18°C to 49°C)

#### MATERIALS:

Drum insert & dip-tube:

Virgin polyethylene (HDPE)

### Shipping plug:

Virgin polyethylene (HDPE)

#### Coupler:

Virgin polypropylene (DrumQuik PRO) Virgin PVDF (DrumQuik PUR)

#### Spring (coupler only):

316 stainless steel, Hastelloy® C, virgin PEEK UHP, PPS or PTFE encapsulated 316 stainless steel

#### Seals:

Coupler seals: EPDM. FKM. Simriz® UHP FFKM Perfluoroelastomer

Bung seal: Polyolefin, EPDM, FKM

Shipping plug seal: EPDM, FKM

#### **Torque specifications:**

2" Buttress - 20ft-lbs (27 N-m) BCS 56x4 - 20ft-lbs (27 N-m) BCS 70mm - 22ft-lbs (30 N-m)

2" NPS (plastic drum) - 20ft-lbs (27 N-m) 2" NPS (steel drum) - 15ft-lbs (20 N-m) Shipping plug - 15inch-lbs (1.7 N-m)

NOTE: The torque specifications listed above are general recommendations which may need to be adjusted + or depending on specific application

# DrumQuik PRO

#### **FEATURES BENEFITS** Quick connect/disconnect Allows faster drum change-outs Stops flow when disconnected Automatic shutoff valve Reduces messy handling and Disposable drum insert chemical exposure Ported vent system Allows make-up air into drum or connection of blanket gas (N2,CO2) for vapor management Multiple standard thread sizes Can be used with many common steel and plastic containers

# DrumQuik PUR

FEATURES	BENEFITS
Quick connect/disconnect	Allows faster drum change-outs
All plastic construction	No metallic corrosion/contamination
Automatic shutoff valve	Stops flow when disconnected
Cleanroom manufactured	Assures high level of purity
Disposable drum insert	Minimizes chemical exposure & contamination

# **Couplers** • POLYPROPYLENE

2	PART NUMBER	DESCRIPTION
PR	DQPR0120312	3/4" Hose barb (EPDM seals, 316 SST spring)
3	DQPR0120208	1/2" Male NPT (EPDM seals, 316 SST spring)
DRUMQUIK	DQPR0120208202	1/2" Male NPT (FKM seals, 316 SST spring)
2	DQPR0120208206	1/2" Male NPT (EPDM seals, Hastelloy C spring)
	DQPR0120208207	1/2" Male NPT (FKM seals, Hastelloy C spring)
	DQPR0120712	3/4" Male BSPP (EPDM seals, 316 SST spring)
	DQPR0120712202	3/4" Male BSPP (FKM seals, 316 SST spring)
	DQPR0120712206	3/4" Male BSPP (EPDM seals, Hastelloy C spring)
	DOPRO120712207	3/// Mala RSPP (FKM spals Hastallov C spring)

## **Couplers** • PVDF (3/4" FLARE TERMINATION)

8	PART NUMBER	DESCRIPTION
UIK PL	DQPUR120112	FKM liquid contact seals, FKM vent seal, PEEK plastic spring
DRUMQUIK	DQPUR120112203	Simriz® FFKM liquid contact seals, FKM vent seal, PEEK plastic spring
_	DQPUR120112208	FKM liquid contact seals, FKM vent seal, PPS plastic spring
	DQPUR120112209	Simriz FFKM liquid contact seals, FKM vent seal, PPS plastic spring
	DQPUR120112210	FKM liquid contact seals, FKM vent seal, PTFE encapsulated 316 stainless steel (TESS) spring
	DQPUR120112211	Simriz FFKM liquid contact seals, FKM vent seal, PTFE encapsulated 316 stainless steel (TESS) spring

## **Drum Inserts** • POLYETHYLENE (HDPE)

~	PART NUMBER	DESCRIPTION
<u>٦</u>	DQPR020BCS70	70mm thread, EPDM seals (70x6 and 70x5 plastic drums)
Š	DQPR020BCS70202	70mm thread, FKM seals (70x6 and 70x5 plastic drums)
ğ	DQPR0202BUT	2" American buttress thread, EPDM seals (No. Am. Plastic drums
2	DQPR0202BUT202	2" American buttress thread, FKM seals (No. Am. Plastic drums
<u>م</u>	DQPR02056X4	BCS 56x4 thread, EPDM seals (European plastic drums)
2	DQPR02056X4202	BCS 56x4 thread, FKM seals (European plastic drums)
A	DQPR020G2NPSST	2" NPS (ISO G2) thread, EPDM seals (steel drums and IBC caps)
ਡ	DQPR020G2NPSST202	2" NPS (ISO G2) thread, FKM seals (steel drums and IBC caps)
DRUMQUIK PRO & DRUMQUIK PUR	DQPR020G2NPSPL	2" NPS (ISO G2) thread, EPDM flat bung seal (plastic drums and IBC caps)
Ī	DQPR020G2NPSPL202	2" NPS (ISO G2) thread, FKM flat bung seal (plastic drums and IBC caps)
	DQPUR202BUT	2" American buttress thread, Polyolefin (POE) bung seal, EPDM shipping plug seal
	DQPUR202BUT202	2" American buttress thread, Polyolefin (POE) bung seal, FKM shipping plug seal


# **Dip-Tubes** • POLYETHYLENE (HDPE)

PUR	PART NUMBER	DESCRIPTION
8 ح	DQPRODT0410	13.0" (330mm) long for 5 gal/25L jerry cans
PRO 8	DQPRODT0985	35.5" (905mm) long for 55 gal/200L drums
	DQPRODT1070	39.0" (990mm) long for 275 gal/1040L IBCs
Š	DQPRODT1274	47.0" (1194mm) long for 330 gal/1250L IBCs
DRUMQUIK	DQPURDT0410	13.0" (330mm) long for 5 gal/25L jerry cans
28	DQPURDT0985	35.5" (905mm) long for 55 gal/200L drums
<b>"</b>	DQPURDT1070	39.0" (990mm) long for 275 gal/1040L IBCs
	DQPURDT1274	47.0" (1194mm) long for 330 gal/1250L IBCs

# Assembled Drum Inserts & Dip-Tubes

~	PART NUMBER	DESCRIPTION
MQUIK PL	DQPUR2BUT0410	Drum insert assembly, NAMPAC square jerry can, 20L, 330mm long (from sealing surface)
DRUMQ	DQPUR2BUT0985	Drum insert assembly with universal dip-tube (drum), 905mm long (from sealing surface)
۵	DQPUR2BUT1070	Drum insert assembly with universal dip-tube (drum), 990mm long (from sealing surface)
	DQPUR2BUT1274	Drum insert assembly with universal dip-tube (IBC), 1194mm long (from sealing surface)

# DrumQuik® PRO & PUR Dispensing Connector System



#### **VAPOR MANAGEMENT OPTIONS:**

- 1. OPEN: Vent port left open to atmosphere
- 2. **REMOTE**: Vent port routed to fume hood, scrubber or vapor recovery
- 3. CHECKED: Connect vent check valve to allow make-up air in without allowing fumes to escape
- 4. BLANKET GAS: Connect a CPC quick disconnect coupling for blanket gas (N2, CO2)
- 5. FILTERED: Connect vent filter

# DRUMQUIK PRO / DRUMQUIK PUR ACCESSORIES

# **Color Coding and Keying Options**

• COLORED KEY KITS FOR PHYSICAL CODING: POLYETHYLENE

DQPRO PART NUMBER	DQPUR PART NUMBER	DESCRIPTION
DQPROKEYCPGRY00	DQPURKEYCPGRY00	Key kit for coupler, gray
DQPROKEYCPRED01	DQPURKEYCPRED01	Key kit for coupler, red
DQPROKEYCPYEL02	DQPURKEYCPYEL02	Key kit for coupler, yellow
DQPROKEYCPGRN03	DQPURKEYCPGRN03	Key kit for coupler, green
DQPROKEYCPBRN04	DQPURKEYCPBRN04	Key kit for coupler, brown
DQPROKEYCPBLU05	DQPURKEYCPBLU05	Key kit for coupler, blue
DQPROKEYCPBLK99	DQPURKEYCPBLK99	Key kit for coupler, black
DQPROKEYCPCOMBO	DQPURKEYCPCOMBO	One of every coupler key ring and cap
Key kits for coupler includ	de coupler key ring and cap.	
DQPROKEYDIRED01	DQPURKEYDIRED01	Key kit for drum insert, red
DQPROKEYDIYEL02	DQPURKEYDIYEL02	Key kit for drum insert, yellow
DQPROKEYDIGRN03	DQPURKEYDIGRN03	key kit for drum insert, green
DQPROKEYDIBRN04	DQPURKEYDIBRN04	Key kit for drum insert, brown
DQPROKEYDIBLU05	DQPURKEYDIBLU05	Key kit for drum insert, blue
DQPROKEYDIBLK99	DQPURKEYDIBLK99	Key kit for drum insert, black
DQPROKEYDICOMBO	DQPURKEYDICOMBO	One of every drum insert key ring and button plug

# **Color Coded Parts**

• COLORED CAPS & BUTTON PLUGS FOR VISUAL CODING: POLYETHYLENE

PART NUMBER	DESCRIPTION
DQPROCAPGRY00	Colored cap for coupler, gray
DQPROCAPRED01	Colored cap for coupler, red
DQPROCAPYEL02	Colored cap for coupler, yellow
DQPROCAPGRN03	Colored cap for coupler, green
DQPROCAPBRN04	Colored cap for coupler, brown
DQPROCAPBLU05	Colored cap for coupler, blue
DQPROCAPBLK99	Colored cap for coupler, black
DQPROCAPCOMBO	One each of every color cap
DQPROPLUGGRY00	Colored button plug for drum insert, gray
DQPROPLUGRED01	Colored button plug for drum insert, red
DQPROPLUGYEL02	Colored button plug for drum insert, yellow
DQPROPLUGGRN03	Colored button plug for drum insert, green
DQPROPLUGBRN04	Colored button plug for drum insert, brown
DQPROPLUGBLU05	Colored button plug for drum insert, blue
DQPROPLUGBLK99	Colored button plug for drum insert, black
DQPROPLUGCOMBO	One each of every color plug

Note: Colored caps and button plugs can be used with both DrumQuik PRO and



# DrumQuik® PRO & PUR Key Kits



# Color Coding



# DRUMQUIK PRO / DRUMQUIK PUR ACCESSORIES



Fully press dip-tube into drum insert socket.



Place the drum insert assembly in drum.



Turn shipping plug counterclockwise to remove. (Use a large, flat blade or #4 Phillips screwdriver.)



Insert coupler into drum insert assembly.



Turn the coupler locking ringclockwise until stop is reached (approximately 2 full turns); valve is opened during this process.

# Tools

PART NUMBER 2408900	<b>DESCRIPTION</b> Torque socket tool (blue, glass-filled nylon)
2479100	Bung wrench (aluminum)
2375600	Key installation tool for key kits (Shipping plug may also be used to install keys)
DQPURFT08	Flaring Tool, 1/2" tubing
DQPURFT12	Flaring Tool, 3/4" tubing

# Sample Kits (Includes one coupler, one insert and one dip-tube)

PART NUMBER	DESCRIPTION
DQPRONASAMP	2" Buttress drum insert, 905mm dip-tube, 1/2 NPT coupler (EPDM 0-rings, 316 SST spring in coupler)
DQPRONASAMP207	2" Buttress drum insert, 905mm dip-tube, 1/2 NPT coupler (FKM seals & Hastelloy® C spring in coupler)
DQPROEUSAMP	BSC 56x4 drum insert, 905mm dip-tube, 3/4 BSPP coupler (EPDM seals, 316 SST spring in coupler)
DQPR0EUSAMP207	BSC 56x4 drum insert, 905mm dip-tube, 3/4 BSPP coupler (FKM seals & Hastelloy® C spring in coupler)
DQPRONPSSTSAMP	2" NPS (ISO G2) drum insert, 905mm dip-tube, 1/2 NPT coupler (EPDM Bung 0-rings, 316 SST spring in coupler)
DQPRONPSSTSAMP207	2" NPS (ISO G2) drum insert, 905mm dip-tube, 1/2 NPT coupler (FKM seals & Hastelloy® C spring in coupler)
DQPRONPSPLSAMP	2" NPS (ISO G2) drum insert, 905mm dip-tube, 1/2 NPT coupler (EPDM Flat Gasket, 316 SST Spring)
DQPRONPSPLSAMP207	2" NPS (ISO G2) drum insert, 905mm dip-tube, 1/2 NPT coupler (FKM seals & Hastelloy® C spring in coupler)
DQPURSAMP202	Kit includes 1 Coupler (DQPUR 12 0112 (PEEK® & FKM) & 1 assembled Dip-tube (DQPUR 2BUT 0985 202)
DQPURSAMP209	Kit includes 1 Coupler (DQPUR 12 0112 209 (PPS & Simriz®) & 1 Dip-tube assy. (DQPUR 2BUT 0985 202)

# Shipping Plugs

DART WWDER	DESCRIPTION
PART NUMBER	DESCRIPTION
DQPROPUMPAD	Internal 3/4 NPT Thread pump adaptor shipping plug (HDPE, EPDM)
DQPROPUMPAD202	Internal 3/4 NPT Thread pump adaptor shipping plug (HDPE, FKM)
DQPRPUMPAD202	Extended shipping plug for venting internal pressure from
DQPROSHPLUGEX	drum (HDPE, EPDM)
DQPROSHIPLUGEX202	Extended shipping plug for venting internal pressure from drum (HDPE, FKM)
2007900	Shipping plug for drum insert (Does not include shipping plug o-ring)
2290100	Shipping plug o-ring, EPDM
2307200	Shipping plug o-ring, FKM
DQPURSP	Shipping Plug with EPDM seal
DOPURSP202	Shipping Plug with FKM seal

# Sample Kits (Includes one coupler, one insert and one dip-tube)

PART NUMBER	DESCRIPTION
DQPURDISEALPOE	Drum insert POE bung seal
DQPURDISEALFEP	Drum insert FEP encapsulated silicone bung seal
DQPURCPSEALFKM	Coupler seal kit (2 seals), FKM
DQPURCPSEALSMZ	Coupler seal kit (2 seals), Simriz®
2423400	Coupler o-ring kit, EPDM (Includes 2 o-rings)
2423500	Coupler o-ring kit, FKM (Includes 2 o-rings)
2166000	Drum insert flat gasket, 2" buttress, FDA EPDM
2166100	Drum insert flat gasket, 2" buttress, FKM
2165800	Drum insert o-ring, 2" NPS (ISO G2) EPDM for steel drums
2165900	Drum insert o-ring, 2" NPS (ISO G2) FKM for steel drums
2165800	Drum insert o-ring, BCS 56x4, EPDM (same as ISO G2 for steel)
2165900	Drum insert o-ring, BCS 56x4, FKM (same as ISO G2 for steel)
2525200	Drum insert flat gasket, 2" NPS (ISO G2) EPDM for plastic drums
2525300	Drum insert flat gasket, 2" NPS (ISO G2) FKM for plastic drums
2594900	Replacement Drum Insert Seal, 2" Buttress, FDA
DQPR070X6AD	70 X 6 Adaptor with EPDM O-Ring
DQPR070X6AD202	70 X 6 Adaptor with FKM (Viton®) O-Rings

# Accessories

PAKI NUMBEK	DESCRIPTION
TETHER/DUST CAP	
DOPROTETHER	Tether/Dust Can for use with DrumQuik PRO shinning plugs and coupler

#### **COUPLER VENT PORT OPTIONS**

2879800	Vent check valve, 1/2" ID x 1/4" male NPT (316 SST, Santoprene)
2879801	Vent check valve, 1/2" ID x 1/4" male NPT (Hastelloy® C, FKM)
DQPR0FIL01	Air filter, coupler vent port, 1/4" male NPT (polypropylene, 70 micron)
APCD24004	Vent coupling insert, 1/4" male NPT (Acetal, 316 SST, Buna-N)
APCD17004SH	Vent coupling body with shroud, 1/4" hose barb (Acetal, 316 SST, Buna-N)
APCD17006SH	Vent coupling body with shroud, 3/8" hose barb (Acetal, 316 SST, Buna-N)
EFCD24412	Vent coupling insert, 1/4" male NPT (Polypropylene, 316 SST, EPDM)
EFCD17412	Vent coupling body, 1/4" hose barb (Polypropylene, 316 SST, EPDM)
EFCD17612	Vent coupling body, 3/8" hose barb (Polypropylene, 316 SST, EPDM)
NOTF: Parts listed above a	are suggested for use in the vent port, but any 1/4" NPT CPC coupling can be used for

blanket gas applications (see page 144)

BACKFLUW PREV	ENTION OPTIONS • POLYPROPYLENE
2478900	Check valve for couplers, 1/2" female NPT x 1/2" male NPT (EPDM, 302 SST)
2478901	Check valve for couplers, 1/2" female NPT x 1/2" male NPT (FKM, 302 SST)
2478902	Check valve for couplers, 1/2" female NPT x 1/2" male NPT (EPDM, Hastelloy C)
2478903	Check valve for couplers, 1/2" female NPT x 1/2" male NPT (FKM, Hastelloy C)
2479000	Foot valve for dip-tube, 1/2" hose barb (EPDM, 316 SST)
2479001	Foot valve for dip-tube, 1/2" hose barb (FKM, 316 SST)
2479002	Foot valve for dip-tube, 1/2" hose barb (EPDM, Hastelloy C)
2479003	Foot valve for dip-tube, 1/2" hose barb (FKM, Hastelloy C)
NOTE, CPC distrib	autos the above check valves and feet valves. These valves are not produced by CPC and CPC

makes no claims or warranty regarding suitability for use in specific applications.

COUPLEK LIQUID PUKI F	ITTINGS • PULTPRUPTLENE
2402600	1/2" female NPT x 3/8" OD compression
2402700	1/2" female NPT x 1/2" OD compression
2402800	1/2" female NPT x 5/8" OD compression
2402900	1/2" female NPT x 3/8" ID hose barb
2403000	1/2" female NPT x 1/2" ID hose barb
2403100	1/2" female NPT x 3/4" ID hose barb

#### DQPUR ACCESSORY FITTINGS

DQPUR MC0404	Male Connector, PVDF, 1/4" Flare x 1/4" NPT
DQPUR MC0604	Male Connector, PVDF, 3/8" Flare x 1/4" NPT
DQPUR MC0804	Male Connector, PVDF, 1/2" Flare x1/4" NPT
DQPUR FR1208	Flare Reducer, PFA, 3/4" Female Flare x 1/2" Male Flare
DQPUR FI1216	Flare Increaser, PFA, 3/4" Female Flare x 1" Male Flare
DQPUR TFMC1208	Female flare x male connector, PFA, 3/4" Female Flare x 1/2" Male NPT
DQPUR SE0808	Sweep Elbow, PFA, 1/2" Female Flare x 1/2" Male Flare
DQPURSE1212	Sweep Elbow, PFA, 3/4" Female Flare x 3/4" Male Flare
DQPURDCPM0812	Double Containment Panel Mount Fitting, PFA, 1/2" process line x 3/4" secondary containment line
DQPURDCPM1216	Double Containment Panel Mount Fitting, PFA, 3/4" process line x 1" secondary containment line
DQPURDCN0812	Double Containment Adapter, PFA, 1/2" process line x 3/4" secondary containment line
DQPURDCN1216	Double Containment Adapter, PFA, 3/4" process line x 1" secondary containment line

plug or key with shoulder up.

#### Torque Tool



Refer to specifications for torque requirements (page 143)



#### **Backflow Check and Foot Valves**

Prevents liquid from draining back into drum when pump is off





#### **Coupler Liquid Port Fittings**





installation

# DRUMQUIK PRO / DRUMQUIK PUR ACCESSORIES



Universal 3-Port Drum Adaptor Kit (shown with HFC12 and NS4)



#### DrumQuik® 3-Port UDA Kit • POLYPROPYLENE

The 3-Port UDA features a 3/4" male NPT thread that mates with common containers, bung plugs and caps. The adaptor ports enable simultaneous chemical extraction, venting or application of blanket gas. In addition, the third port can be used to recirculate priming chemical to the container in metering pump

Optional system components including foot valves, air filters and vent check valves are also available. For your convenience, CPC offers a variety of drum adaptor plugs to facilitate easy drum connections. (See plugs and caps on page 149.)

PART NUMBER	DESCRIPTION
DQUDA3PKIT	3-port adaptor, EPDM seal, locking nut, 35.5" (905mm) dip-tube, 1/2" x 3/8" &
•	1/2" x 1/4" reducer bushings and 1/4" & 3/8" NPT plugs
DQUDA3PKIT202	3-port adaptor, FKM seal, locking nut, 35.5" (905mm) dip-tube, 1/2" x 3/8" &
	1/2" x 1/4" reducer bushings and 1/4" & 3/8" NPT plugs
DQUDA3P	3-port adaptor bulk packaged (1/2" female NPT liquid outlet port, 3/8" & 1/4"
	NPT vent/recirculation ports)
DQUDA3PNUT	Locking nut for the 3-port adaptor
2717000	EPDM seal for the 3-port adaptor
2717100	FKM seal for the 3-port adaptor
2621100	Reducer bushing, black polypropylene, 1/2" male NPT x 3/8" female NPT
2621200	Reducer bushing, black polypropylene, 1/2" male NPT x 1/4" female NPT
2621300	Pipe plug, black polypropylene, 1/4" male NPT
2621400	Pipe plug, black polypropylene, 3/8" male NPT
	Use with PLC, Non-Spill, EFC and ChemQuik Series couplings.

Note: DrumQuik® LIDA 3-Port uses DrumQuik PRO din-tubes and accessories

### DrumQuik UDA Kit • POLYPROPYLENE

The DrumQuik® UDA Kit turns ChemQuik® or general purpose (HFC or NSH) couplings into a dip-tube that can easily be threaded into the 3/4" female NPS port of common drum closures. It provides a simple and inexpensive way to dispense chemicals in pumped systems. For your convenience, CPC offers a variety of drum adaptor plugs to facilitate easy drum connections. (See plugs and caps on page 149.)

Press together to assemble			
Dip-tube	Bushing		Y
		Coupling Insert	Coupling Body

PART NUMBER	DESCRIPTION
DQUDAKIT	Adaptor & DQPRODT1274 47" long dip-tube



Universal Drum Adaptor Kit (shown with NSH Series)

used with many standard CPC quick disconnect couplings

## **BottleQuik® 2-Port Kit UDA • POLYPROPYLENE**

The BottleQuik® is intended for use on reagent bottles with SP400-38mm threads. It can be used in either the upright position or inverted and has two 3/8" female NPT ports for liquid and vent ports.

DQUDACAP2PKITPP

2-Port adaptor CAP for 38mm SP400 bottles & dip-tubes, 17.75" (451mm) long, FKM o-ring. 3/8" female NPT liquid port, 3/8" female NPT vent port, 3/8" female NPT vent

Use above kits with CQG, CQH and CQV Series couplings.

## DrumQuik Asian Drum Adaptors • POLYPROPYLENE

The DrumQuik® Asian drum adaptors allow ChemQuik® or general purpose (HFC or NSH) couplings to be easily connected to the unique threads common in many drum closures (with integral dip-tubes) manufactured in Asia, e.g., Kodama, Accelo, Stella, Dung Woo, etc. A 3/4" female NPT inner thread can accept any coupling or fitting with a 3/4" male NPT termination.

#### PART NUMBER DESCRIPTION

DQADAPPO212 • TYPE A Adaptor, FKM seals, 3/4" female NPT inner port x 3/4" NPS male thread to fit drum

DQADAPP0622 • TYPE B Adaptor, FKM seals, 3/4" female NPT inner port x 1-3/8" TYPE B male buttress thread to fit drum plug.

# **DrumQuik Plugs and Caps**

For your convenience, CPC offers several sizes of common bung plugs and caps that have the internal 3/4" female NPS thread, into which the DrumQuik Universal Drum Adaptors/ChemQuik couplings are threaded.

#### PART NUMBER

DQMBUNGPP2BUT12 DQMBUNGPP2NPS12 DQMBUNGPP56X412 DOMBLINGPP70X612 DOMCAPPP70mm DQMCAPPP63MM

Bung closure, polypropylene, EPDM o-ring, 2" buttress x 3/4" female NPS Bung closure, polypropylene, EPDM o-ring, 2" NPS x 3/4" female NPS Bung closure, polypropylene, EPDM o-ring, BCS56x4 x 3/4" female NPS Bung closure, polypropylene, EPDM o-ring, BCS70x6, 3/4" female NPS Cap closure, polypropylene, EPDM o-ring, 70mm x 3/4" female NPS Cap closure, polypropylene, EPDM o-ring, 63mm, 3/4" female NPS

All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters. **NOTE:** The DQ UDA kits are designed to work with CPC couplings that have 3/4" NPT male terminations. The DQ UDA 3-Port kits are designed to work with a 1/2" NPT male on the liquid port and a 1/4" NPT male on the vent port. A 1/2" male NPT to 1/4" female NPT reducer bushing may be required to accommodate smaller NS4, PLC and PMC Series couplings. BottleQuik features two 3/8" female NPT ports. NOTE: DrumQuik Adaptors are intended for end-use applications not requiring compliance with UN/DOT standards for use in transportation of hazardous materials



BottleQuik 2-Port UDA Bottle Adaptor Kit (shown with CQG and CQH)



For use with closure plugs with a 3/4" NPS

center port.



For use with closure plugs with a 1-3/8" buttress center port.





3/4" NPT inner port incorporates a face seal to assure a leak-free connection



DQ UDA kits can be used with a wide range of closures.



# Specifications • • •

#### PRESSURE:

20" Hg vacuum to 15 psi, 1.0 bar (body only)

#### TEMPERATURE:

32°F to 160°F (0°C to 71°C)

#### MATERIALS:

#### Main body components and valves:

Polypropylene

#### Springs:

Body: 316 stainless steel Cap: 302 stainless steel

**UDC caps:** Acetal, polypropylene

#### O-rings (food grade):

UDC: EPDM, silicone, Buna-N

iUDC: EPDM. silicone

CIP adapter: Acetal

COLOR: Light gray with dark gray latch

CLOSURE SIZE: SP-400 thread/38mm

snap-in

#### **TUBING SIZES:**

3/8" to 3/4" ID, 9.5mm to 19.0mm ID

# under NSF/ANSI standard 169

# The patented 3/8" flow Universal Dispensing Coupler

(UDC) provides a universal connection to a 38mm threaded or snap-in fitment neck. Make instant connections to bulk packaging systems including bag-in-box (BIB), flexible and rigid packaging styles. Automatic flush face valves minimize costly or dangerous spillage. An ergonomic design and a large, shrouded thumb latch pad produce a coupling that is easy to grip and simple to operate.

The UDC Series is also available with optional RFID (Radio Frequency Identification) capability (see page 95).

FEATURES	BENEFITS	
Universal 38mm cap	Fits industry standard fitment necks	
Flush face/non-spill valves	Minimize spillage	
Integral terminations	Fewer leak points, shorter assemblies, faster installations	
Shrouded thumb latch	Protects against accidental disconnects	
RFID technology (iUDC only)	Helps maintain warranty, integrity, brand protection	



# **UDC SERIES DIMENSIONS**

# **Coupling Bodies** • POLYPROPYLENE



# Caps • ACETAL



TERMINATION	DESCRIF
THREAD-ON	Without
38MM VALVED	With thre

DESCRIPTION	SILICONE O-RING
Without thread seal	95800
With thread seal	95801

EPDM O-RING 97400 97401

A	В
.70	1.47
.70	1.47

### • POLYPROPYLENE



TERMINATION SNAP-IN 38MM VALVED
THREAD-ON 38MM VALVED

<b>DESCRIPTION Note:</b> Snap-in cl	SILICONE O-RING 9461300 osure designed to fit Scholle 900 bag necks.	<b>EPDM O-RING</b> 9500000	A 1.15	B 1.77	
With thread seal		9462000	1.70	1.48	

BUNA-N O-RING

96400

96401

# Clean-in-Place (CIP) Adapter • ACETAL



IN-LINE

1/2" HB (no o-ring) 96000

A 1.15

SHUTOFF

iUDCDR17

iUDCDR17

iUDCDR17

iUDCDR1600401

iUDCDR1600601

iUDCDR1600801

B 1.56

# **Coupling Readers • POLYPROPYLENE**



IN-LINE HOSE BARB
PANEL MOUNT

HOSE BARB	
PANEL MOUNT HOSE BARB	

TUBING SIZE	METRIC EQ.	SHUTOFF (RS-232
1/4" ID	6.4mm ID	iUDCDR1700400
3/8" ID	9.5mm ID	iUDCDR1700600
1/2" ID	12.7mm ID	iUDCDR1700800
5/8" ID	15.9mm ID	iUDCDR1701000
3/4" ID	19mm ID	iUDCDR1701200
***************************************		
1/4" ID	6.4mm ID	iUDCDR1600400
3/8" ID	9.5mm ID	iUDCDR1600600

(USB)	A	B
700401	2.00	3.49
700601	2.00	3.49
700801	2.00	3.49
	2.00	3.49
	2.00	2.40

3.49 3.49 Refer to cpcworldwide.com

for dimensions

# Caps • ACETAL



**TERMINATION** THREAD-ON 38MM VALVED

**DESCRIPTION** Without thread seal With thread seal

1/2" ID

METRIC EQ. 38mm 38mm

12.7mm ID

SHIITOFF iUDCDTCN3803 iUDCDTCS3803

SHUTOFF (RS-232)

iUDCDR1600800

A 1.70

B 1.47 1.47

#### POLYPROPYLENE



TERMINATION SNAP-IN 38MM VALVED	DESCRIPTION	METRIC EQ. 38mm	SHUTOFF iUDCDTSN0000	A 1.62	B 1.76
THREAD-ON	Without thread seal	38mm	iUDCDTCN3800	1.70	1.48
38MM VALVED	With thread seal	38mm	iUDCDTCS3800	1.70	1.48

All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters. Couplings are pictured with valves unless otherwise noted

# PUNCTURE SEAL SERIES CONNECTOR



The 1/4" flow Puncture Seal dispensing system provides clean and safe liquid dispensing from bag-in-box (BIB) or other types of flexible packaging. The system consists of a Puncture Seal closure cap and a specially designed quick disconnect coupling body. The closure makes a positive thread or snap-in seal with a 38mm fitment neck. The Puncture Seal coupling body automatically breaks the membrane seal on the closure cap when connected. To switch to a new container, simply depress the CPC thumb latch to disconnect the system and then connect the Puncture Seal coupling body to a fresh container.

#### **FEATURES BENEFITS** Easily breaks the Puncture Seal Puncture Seal coupling Polypropylene and Greater chemical resistance polyethylene materials Puncture Seal membrane Eliminates shipping caps and o-rings CPC thumb latch One-hand connection and disconnection

# Specifications • • •

#### PRESSURE:

120 psi, 8.3 bar (coupling body)

#### TEMPERATURE:

32°F to 160°F (0°C to 71°C)

#### **MATERIALS:**

Main components and valves: Polypropylene 38mm closure: Polyethylene (LDPE)

Thumb latch: Stainless steel Valve spring: 316 stainless steel External spring and pin: Stainless steel


**0-rings:** EPDM

CLOSURE SIZE: SP-400 38mm

#### **TUBING SIZES:**

1/4" and 3/8" ID, 6.4mm and 9.5mm ID


# Closure Detail



# DID YOU KNOW?

### CPC'S PUNCTURE SEAL CLOSURES CAN BE USED TO DISPENSE VIRTUALLY ANY LIQUID FROM FLEXIBLE PACKAGES?

The Puncture Seal makes transporting and dispensing cleaning concentrates, dry cleaning chemicals and soaps cleaner and safer. Printing inks, photo developers and dying agents can also be dispensed quickly and easily, without messy drips or spills. The Puncture Seal can also be used with beverage and condiment dispensing systems, just look for the NSF-listed part numbers.



# **Coupling Bodies • POLYPROPYLENE**



## **Closure** • POLYETHYLENE

THREAD-ON


VALVED

Note: Snap-in closure designed to fit Scholle 900 bag necks.



**TERMINATION** PART NO. LENGTH HEIGHT PSC38mm 1.12 38MM VALVED SNAP-IN 38MM 496200 1.88 1.62







**FitQuik® Connectors** from CPC are high quality fittings for leak-free tubing connections. These precision molded fittings are designed to eliminate tubing leak points in applications such as medical devices, analytical instrumentation or air-driven equipment.

Cost-competitive FitQuik Connectors are offered in many standard configurations and in custom designs to meet unique customer specifications.

**Designed to fit.** Quality FitQuik® Connectors are instrumentation - grade fittings that ensure a secure fit in even the most demanding applications. These fittings feature precise hose barbs and thread forms to ensure the right fit every time.

**Materials available:** Black and white nylon, polypropylene, PVDF and polycarbonate (available in medical-grade and USP Class VI, ADCF)

# TUBE TO TUBE FITTINGS

# Tube to Tube Fittings - Straight

	PART NO. HS2 HS230 HS231 HS270 HS291	DESCRIPTION  Polypropylene Straight Fitting, 1/16" Hose Barbs White Nylon Straight Fitting, 1/16" Hose Barbs Black Nylon Straight Fitting, 1/16" Hose Barbs PVDF Straight Fitting, 1/16" Hose Barb Purple Tint Polycarbonate Straight Fitting, 1/16" Hose Barbs
	HS3 HS330 HS331 HS370 HS391	Polypropylene Straight Fitting, 3/32" Hose Barbs White Nylon Straight Fitting, 3/32" Hose Barbs Black Nylon Straight Fitting, 3/32" Hose Barbs PVDF Straight Fitting, 3/32" Hose Barb Purple Tint Polycarbonate Straight Fitting, 3/32" Hose Barbs
	HS4 HS430 HS431 HS470 HS491	Polypropylene Straight Fitting, 1/8" Hose Barbs White Nylon Straight Fitting, 1/8" Hose Barbs Black Nylon Straight Fitting, 1/8" Hose Barbs PVDF Straight Fitting, 1/8" Hose Barb Purple Tint Polycarbonate Straight Fitting, 1/8" Hose Barbs
Part P	HS5 HS530 HS531 HS570 HS591	Polypropylene Straight Fitting, 5/32" Hose Barbs White Nylon Straight Fitting, 5/32" Hose Barbs Black Nylon Straight Fitting, 5/32" Hose Barbs PVDF Straight Fitting, 5/32" Hose Barb Purple Tint Polycarbonate Straight Fitting, 5/32" Hose Barbs
	HS6 HS630 HS631 HS670 HS691	Polypropylene Straight Fitting, 3/16" Hose Barbs White Nylon Straight Fitting, 3/16" Hose Barbs Black Nylon Straight Fitting, 3/16" Hose Barbs PVDF Straight Fitting, 3/16" Hose Barb Purple Tint Polycarbonate Straight Fitting, 3/16" Hose Barbs
	HS8 HS830 HS831 HS870 HS891	Polypropylene Straight Fitting, 1/4" Hose Barbs White Nylon Straight Fitting, 1/4" Hose Barbs Black Nylon Straight Fitting, 1/4" Hose Barbs PVDF Straight Fitting, 1/4" Hose Barb Purple Tint Polycarbonate Straight Fitting, 1/4" Hose Barbs
	HS12 HS1230 HS1231 HS1270	Polypropylene Straight Fitting, 3/8" Hose Barbs White Nylon Straight Fitting, 3/8" Hose Barbs Black Nylon Straight Fitting, 3/8" Hose Barbs PVDF Straight Fitting, 3/8" Hose Barbs

Tube to Tube Fittings - Elbow

HE230

HE231

HE330 HE331

HE370

# Tube to Tube Fittings - Elbow (cont.)

~	PART NO. HE4 HE430 HE431 HE470 HE491	DESCRIPTION Polypropylene Elbow Fitting, 1/8" Hose Barbs White Nylon Elbow Fitting, 1/8" Hose Barbs Black Nylon Elbow Fitting, 1/8" Hose Barbs PVDF Elbow Fitting, 1/8" Hose Barbs Purple Tint Polycarbonate Elbow Fitting, 1/8" Hose Barbs
1	HE5 HE530 HE531 HE570 HE591	Polypropylene Elbow Fitting, 5/32" Hose Barbs White Nylon Elbow Fitting, 5/32" Hose Barbs Black Nylon Elbow Fitting, 5/32" Hose Barbs PVDF Elbow Fitting, 5/32" Hose Barbs Purple Tint Polycarbonate Elbow Fitting, 5/32" Hose Barbs
	HE6 HE630 HE631 HE670 HE691	Polypropylene Elbow Fitting, 3/16" Hose Barbs White Nylon Elbow Fitting, 3/16" Hose Barbs Black Nylon Elbow Fitting, 3/16" Hose Barbs PVDF Elbow Fitting, 3/16" Hose Barbs Purple Tint Polycarbonate Elbow Fitting, 3/16" Hose Barbs
	HE8 HE830 HE831 HE870 HE891	Polypropylene Elbow Fitting, 1/4" Hose Barbs White Nylon Elbow Fitting, 1/4" Hose Barbs Black Nylon Elbow Fitting, 1/4" Hose Barbs PVDF Elbow Fitting, 1/4" Hose Barbs Purple Tint Polycarbonate Elbow Fitting, 1/4" Hose Barbs
	HE12 HE1230 HE1231 HE1270	Polypropylene Elbow Fitting, 3/8" Hose Barbs White Nylon Elbow Fitting, 3/8" Hose Barbs Black Nylon Elbow Fitting, 3/8" Hose Barbs PVDF Elbow Fitting, 3/8" Hose Barbs

FitQuik polycarbonate fittings are compatible with autoclave sterilization. See sterilization guide on page 163.

# Tube to Tube Fittings - Tee

HT2

7	HT230 HT231 HT270 HT291	White Nylon Tee Fitting, 1/16" Hose Barbs Black Nylon Tee Fitting, 1/16" Hose Barbs PVDF Tee Fitting, 1/16" Hose Barbs Purple Tint Polycarbonate Tee Fitting, 1/16" Hose Barbs
7	HT3 HT330 HT331 HT370 HT391	Polypropylene Tee Fitting, 3/32" Hose Barbs White Nylon Tee Fitting, 3/32" Hose Barbs Black Nylon Tee Fitting, 3/32" Hose Barbs PVDF Tee Fitting, 3/32" Hose Barbs Purple Tint Polycarbonate Tee Fitting, 3/32" Hose Barbs
7	HT4 HT430 HT431 HT470 HT491	Polypropylene Tee Fitting, 1/8" Hose Barbs White Nylon Tee Fitting, 1/8" Hose Barbs Black Nylon Tee Fitting, 1/8" Hose Barbs PVDF Tee Fitting, 1/8" Hose Barbs Purple Tint Polycarbonate Tee Fitting, 1/8" Hose Barbs
7	HT5 HT530 HT531 HT570 HT591	Polypropylene Tee Fitting, 5/32" Hose Barbs White Nylon Tee Fitting, 5/32" Hose Barbs Black Nylon Tee Fitting, 5/32" Hose Barbs PVDF Tee Fitting, 5/32" Hose Barbs Purple Tint Polycarbonate Tee Fitting, 5/32" Hose Barbs

Polypropylene Tee Fitting, 1/16" Hose Barbs

#### Note: Parts not shown to scale. For customer drawings and technical specifications please visit cpcworldwide.com

Polypropylene Elbow Fitting, 1/16" Hose Barbs

Purple Tint Polycarbonate Elbow Fitting, 1/16" Hose Barbs

Purple Tint Polycarbonate Elbow Fitting, 3/32" Hose Barbs

White Nylon Elbow Fitting, 1/16" Hose Barbs Black Nylon Elbow Fitting, 1/16" Hose Barbs

Polypropylene Elbow Fitting, 3/32" Hose Barbs White Nylon Elbow Fitting, 3/32" Hose Barbs

Black Nylon Elbow Fitting, 3/32" Hose Barbs

PVDF Elbow Fitting, 3/32" Hose Barbs

PVDF Elbow Fitting, 1/16" Hose Barbs

# Tube to Tube Fittings - Tee (cont.)

	PART NO. HT6 HT630 HT631 HT670 HT691	DESCRIPTION Polypropylene Tee Fitting, 3/16" Hose Barbs White Nylon Tee Fitting, 3/16" Hose Barbs Black Nylon Tee Fitting, 3/16" Hose Barbs PVDF Tee Fitting, 3/16" Hose Barbs Purple Tint Polycarbonate Tee Fitting, 3/16" Hose Barbs
7	HT8 HT830 HT831 HT870 HT891	Polypropylene Tee Fitting, 1/4" Hose Barbs White Nylon Tee Fitting, 1/4" Hose Barbs Black Nylon Tee Fitting, 1/4" Hose Barbs PVDF Tee Fitting, 1/4" Hose Barbs Purple Tint Polycarbonate Tee Fitting, 1/4" Hose Barbs
7	HT12 HT1230 HT1231 HT1270	Polypropylene Tee Fitting, 3/8" Hose Barbs White Nylon Tee Fitting, 3/8" Hose Barbs Black Nylon Tee Fitting, 3/8" Hose Barbs PVDF Tee Fitting, 3/8" Hose Barbs

PART NO. HY2 HY230 HY231 HY270 HY291	DESCRIPTION  Polypropylene Y Fitting, 1/16" Hose Barbs White Nylon Y Fitting, 1/16" Hose Barbs Black Nylon Y Fitting, 1/16" Hose Barbs PVDF Y Fitting, 1/16" Hose Barbs Purple Tint Polycarbonate Y Fitting, 1/16" Hose Barbs
HY3	Polypropylene Y Fitting, 3/32" Hose Barbs
HY330	White Nylon Y Fitting, 3/32" Hose Barbs
HY331	Black Nylon Y Fitting, 3/32" Hose Barbs
HY370	PVDF Y Fitting, 3/32" Hose Barbs
HY391	Purple Tint Polycarbonate Y Fitting, 3/32" Hose Barbs
HY4	Polypropylene Y Fitting, 1/8" Hose Barbs
HY430	White Nylon Y Fitting, 1/8" Hose Barbs
HY431	Black Nylon Y Fitting, 1/8" Hose Barbs
HY470	PVDF Y Fitting, 1/8" Hose Barbs
HY491	Purple Tint Polycarbonate Y Fitting, 1/8" Hose Barbs
HY5	Polypropylene Y Fitting, 5/32" Hose Barbs
HY530	White Nylon Y Fitting, 5/32" Hose Barbs
HY531	Black Nylon Y Fitting, 5/32" Hose Barbs
HY570	PVDF Y Fitting, 5/32" Hose Barbs
HY591	Purple Tint Polycarbonate Y Fitting, 5/32" Hose Barbs
HY6	Polypropylene Y Fitting, 3/16" Hose Barbs
HY630	White Nylon Y Fitting, 3/16" Hose Barbs
HY631	Black Nylon Y Fitting, 3/16" Hose Barbs
HY670	PVDF Y Fitting, 3/16" Hose Barbs
HY691	Purple Tint Polycarbonate Y Fitting, 3/16" Hose Barbs
HY8	Polypropylene Y Fitting, 1/4" Hose Barbs
HY830	White Nylon Y Fitting, 1/4" Hose Barbs
HY831	Black Nylon Y Fitting, 1/4" Hose Barbs
HY870	PVDF Y Fitting, 1/4" Hose Barbs
HY891	Purple Tint Polycarbonate Y Fitting, 1/4" Hose Barbs
HY12	Polypropylene Y Fitting, 3/8" Hose Barbs
HY1230	White Nylon Y Fitting, 3/8" Hose Barbs
HY1231	Black Nylon Y Fitting, 3/8" Hose Barbs
HY1270	PVDF Y Fitting, 3/8" Hose Barbs

# Tube to Tube Fittings - Straight Reducer

านมธ เบ	IUDG II	ittiigs - Straight Neudoei
	PART NO.	DESCRIPTION
	HSR32	Polypropylene Reducer Fitting, 3/32" to 1/16" Hose Barb
	MSR3230	White Nylon Reducer Fitting, 3/32" to 1/16" Hose Barb
	HSR3231	Black Nylon Reducer Fitting, 3/32" to 1/16" Hose Barb
	HSR3270	PVDF Straight Reducer Fitting, 3/32" to 1/16" Hose Barb
	HSR3291	Purple Tint Polycarbonate Reducer Fitting, 3/32" to 1/16"
		Hose Barb
***************************************		D
0.0	HSR42	Polypropylene Reducer Fitting, 1/8" to 1/16" Hose Barb
	HSR4230	White Nylon Reducer Fitting, 1/8" to 1/16" Hose Barb
The same of the sa	HSR4231	Black Nylon Reducer Fitting, 1/8" to 1/16" Hose Barb
	HSR4270	PVDF Straight Reducer Fitting, 1/8" to 1/16" Hose Barb
	HSR4291	Purple Tint Polycarbonate Reducer Fitting, 1/8" to 1/16" Hose Barb
	HSR43	Polypropylene Reducer Fitting, 1/8" to 3/32" Hose Barb
	HSR4330	White Nylon Reducer Fitting, 1/8" to 3/32" Hose Barb
The state of the s	HSR4331	Black Nylon Reducer Fitting, 1/8" to 3/32" Hose Barb
	HSR4370	PVDF Straight Reducer Fitting, 1/8" to 3/32" Hose Barb
	HSR4391	Purple Tint Polycarbonate Reducer Fitting, 1/8" to 3/32"
		Hose Barb
••••	повео	Delugrapulana Paduaar Eitting 5/2011 to 1/101 Head Doub
60	HSR52	Polypropylene Reducer Fitting, 5/32" to 1/16" Hose Barb White Nylon Reducer Fitting, 5/32" to 1/16" Hose Barb
1	HSR5230	White Nylon Reducer Fitting, 5/32" to 1/16" Hose Barb Black Nylon Reducer Fitting, 5/32" to 1/16" Hose Barb
	HSR5231 HSR5270	PVDF Straight Reducer Fitting, 5/32" to 1/16" Hose Barb
Salar Sa	HSR5270	Purple Tint Polycarbonate Reducer Fitting, 5/32" to 1/16"
	HONJEJI	Hose Barb
***************************************		
	HSR53	Polypropylene Reducer Fitting, 5/32" to 3/32" Hose Barb
	MSR5330	White Nylon Reducer Fitting, 5/32" to 3/32" Hose Barb
and the same	HSR5331	Black Nylon Reducer Fitting, 5/32" to 3/32" Hose Barb
	HSR5370	PVDF Straight Reducer Fitting, 5/32" to 3/32" Hose Barb
The same of the sa	HSR5391	Purple Tint Polycarbonate Reducer Fitting, 5/32" to 3/32"
		Hose Barb
***************************************	LICDE A	Delugranulana Dadusay Fitting F/201 to 1/01 Hose Dayb
	HSR54 HSR5430	Polypropylene Reducer Fitting, 5/32" to 1/8" Hose Barb
	HSR5431	White Nylon Reducer Fitting, 5/32" to 1/8" Hose Barb Black Nylon Reducer Fitting, 5/32" to 1/8" Hose Barb
	HSR5470	PVDF Straight Reducer Fitting, 5/32" to 1/8" Hose Barb
	HSR5491	Purple Tint Polycarbonate Reducer Fitting, 5/32" to 1/8"
	110110431	Hose Barb
***************************************		
40	HSR62	Polypropylene Reducer Fitting, 3/16" to 1/16" Hose Barb
1	HSR6230	White Nylon Reducer Fitting, 3/16" to 1/16" Hose Barb
	HSR6231	Black Nylon Reducer Fitting, 3/16" to 1/16" Hose Barb
-	HSR6270	PVDF Straight Reducer Fitting, 3/16" to 1/16" Hose Barb
	HSR6291	Purple Tint Polycarbonate Reducer Fitting, 3/16" to 1/16"
		Hose Barb
_	HSR63	Polypropylene Reducer Fitting, 3/16" to 3/32" Hose Barb
	HSR6330	White Nylon Reducer Fitting, 3/16" to 3/32" Hose Barb
	HSR6331	Black Nylon Reducer Fitting, 3/16" to 3/32" Hose Barb
and the same of th	HSR6370	PVDF Straight Reducer Fitting, 3/16" to 3/32" Hose Barb
	HSR6391	Purple Tint Polycarbonate Reducer Fitting, 3/16" to 3/32"
		Hose Barb
A CONTRACTOR	HSR64	Polypropylene Reducer Fitting, 3/16" to 1/8" Hose Barb
	HSR6430	White Nylon Reducer Fitting, 3/16" to 1/8" Hose Barb
The state of the s	HSR6431	Black Nylon Reducer Fitting, 3/16" to 1/8" Hose Barb
	HSR6470	PVDF Straight Reducer Fitting, 3/16" to 1/8" Hose Barb
	HSR6491	Purple Tint Polycarbonate Reducer Fittings, 3/16" to 1/8" Hose Barb
\	.,X.,	nuse Dally
X/X/X	HSR65	Polypropylene Reducer Fitting, 3/16" to 5/32" Hose Barb
	HSR6530	White Nylon Reducer Fitting, 3/16" to 5/32" Hose Barb
	HSR6531	Black Nylon Reducer Fitting, 3/16" to 5/32" Hose Barb
	HSR6570	PVDF Straight Reducer Fitting 3/16" to 5/32" Hose Barh

HSR6570 PVDF Straight Reducer Fitting, 3/16" to 5/32" Hose Barb

Purple Tint Polycarbonate Reducer Fittings, 3/16" to 5/32"

HSR6591

Hose Barb

Smart fluid handling to take you forward, faster. 157

# TUBE TO TUBE FITTINGS

# Tube to Tube Fittings - Straight Reducer



Polypropylene Reducer Fitting, 1/4" to 1/8" Hose Barb White Nylon Reducer Fitting, 1/4" to 1/8" Hose Barb Black Nylon Reducer Fitting, 1/4" to 1/8" Hose Barb PVDF Straight Reducer Fitting, 1/4" to 1/8" Hose Barb Purple Tint Polycarbonate Reducer Fittings, 1/4" to 1/8"



#### HSR85 HSR8530 HSR8570 HSR8591

HSR84

HSR8430

HSR8431

HSR8470

Polypropylene Reducer Fitting, 1/4" to 5/32" Hose Barb White Nylon Reducer Fitting, 1/4" to 5/32" Hose Barb Black Nylon Reducer Fitting, 1/4" to 5/32" Hose Barb PVDF Straight Reducer Fitting, 1/4" to 5/32" Hose Barb Purple Tint Polycarbonate Reducer Fittings, 1/4" to 5/32" Hose Barb



# HSR8630 HSR8691

Polypropylene Reducer Fitting, 1/4" to 3/16" Hose Barb White Nylon Reducer Fitting, 1/4" to 3/16" Hose Barb Black Nylon Reducer Fitting, 1/4" to 3/16" Hose Barb PVDF Straight Reducer Fitting, 1/4" to 3/16" Hose Barb Purple Tint Polycarbonate Reducer Fittings, 1/4" to 3/16" Hose Barb



Polypropylene Reducer Fitting, 3/8" to 1/4" Hose Barb HSR12830 White Nylon Reducer Fitting, 3/8" to 1/4" Hose Barb HSR12831 Black Nylon Reducer Fitting, 3/8" to 1/4" Hose Barb HSR12870 PVDF Straight Reducer Fitting, 3/8" to 1/4" Hose Barb

# **PLUGS AND** THREAD TO THREAD

# Threaded Plugs



MP30 MP31

Polypropylene Threaded Plug, Tapered 10-32 White Nylon Threaded Plug, Tapered 10-32 Black Nylon Threaded Plug, Tapered 10-32 PVDF Threaded Plug Fitting, Tapered 10-32

Polypropylene Threaded Plug, 1/16" NPT

White Nylon Threaded Plug, 1/16" NPT





Black Nylon Threaded Plug, 1/16" NPT PVDF Threaded Plug, 1/16" NPT Polypropylene Threaded Plug, 1/8" NPT White Nylon Threaded Plug, 1/8" NPT





N8P31

Polypropylene Threaded Plug, 1/4" NPT White Nylon Threaded Plug, 1/4" NPT Black Nylon Threaded Plug, 1/4" NPT PVDF Theaded Plug Fitting, 1/4" NPT

Black Nylon Threaded Plug, 1/8" NPT

PVDF Threaded Plug Fitting, 1/8" NPT

# Thread to Thread -Straight

#### DESCRIPTION

N32 N3230 N3231 N3270

Polypropylene Straight Fitting, Tapered 10-32 Threads White Nylon Straight Fitting, Tapered 10-32 Threads Black Nylon Straight Fitting, Tapered 10-32 Threads PVDF Straight Fitting, Tapered 10-32 Threads

# THREAD TO TUBE FITTINGS

# Thread to Tube Fittings - Straight



#### DESCRIPTION

Polypropylene Straight Fitting, 1/4-28 UNF to 1/16" Hose Barb White Nylon Straight Fitting, 1/4-28 UNF to 1/16" Hose Barb Black Nylon Straight Fitting, 1/4-28 UNF to 1/16" Hose Barb PVDF Straight Fitting, 1/4-28 UNF to 1/16" Hose Barb



KS330 KS331 KS370

KS2

KS230

KS231

KS270

Polypropylene Straight Fitting, 1/4-28 UNF to 3/32" Hose Barb White Nylon Straight Fitting, 1/4-28 UNF to 3/32" Hose Barb Black Nylon Straight Fitting, 1/4-28 UNF to 3/32" Hose Barb PVDF Straight Fitting, 1/4-28 UNF to 3/32" Hose Barb



KS4 Polypropylene Straight Fitting, 1/4-28 UNF to 1/8" Hose Barb KS430 White Nylon Straight Fitting, 1/4-28 UNF to 1/8" Hose Barb KS431 Black Nylon Straight Fitting, 1/4-28 UNF to 1/8" Hose Barb KS470 PVDF Straight Fitting, 1/4-28 UNF to 1/8" Hose Barb



KS5 KS530 KS531 KS570

Polypropylene Straight Fitting, 1/4-28 UNF to 5/32" Hose Barb White Nylon Straight Fitting, 1/4-28 UNF to 5/32" Hose Barb Black Nylon Straight Fitting, 1/4-28 UNF to 5/32" Hose Barb PVDF Straight Fitting, 1/4-28 UNF to 5/32" Hose Barb



MS2 MS230 MS231 MS270

Polypropylene Straight Fitting, Tapered 10-32 to 1/16" Hose Barb White Nylon Straight Fitting, Tapered 10-32 to 1/16" Hose Barb Black Nylon Straight Fitting, Tapered 10-32 to 1/16" Hose Barb PVDF Straight Fitting, Tapered 10-32 to 1/16" Hose Barb



MS330 MS331 MS370

Polypropylene Straight Fitting, Tapered 10-32 to 3/32" Hose Barb White Nylon Straight Fitting, Tapered 10-32 to 3/32" Hose Barb Black Nylon Straight Fitting, Tapered 10-32 to 3/32" Hose Barb PVDF Straight Fitting, Tapered 10-32 to 3/32" Hose Barb





Polypropylene Straight Fitting, Tapered 10-32 to 1/8" Hose Barb White Nylon Straight Fitting, Tapered 10-32 to 1/8" Hose Barb Black Nylon Straight Fitting, Tapered 10-32 to 1/8" Hose Barb PVDF Straight Fitting, Tapered 10-32 to 1/8" Hose Barb



MS5 MS530 MS531 MS570

Polypropylene Straight Fitting, Tapered 10-32 to 5/32" Hose Barb White Nylon Straight Fitting, Tapered 10-32 to 5/32" Hose Barb Black Nylon Straight Fitting, Tapered 10-32 to 5/32" Hose Barb PVDF Straight Fitting, Tapered 10-32 to 5/32" Hose Barb



MS630 MS631

Polypropylene Straight Fitting, Tapered 10-32 to 3/16" Hose Barb White Nylon Straight Fitting, Tapered 10-32 to 3/16" Hose Barb Black Nylon Straight Fitting, Tapered 10-32 to 3/16" Hose Barb PVDF Straight Fitting, Tapered 10-32 to 3/16" Hose Barb



Polypropylene Straight Fitting, 1/16" NPT to 1/16" Hose Barb N2S230 White Nylon Straight Fitting, 1/16" NPT to 1/16" Hose Barb Black Nylon Straight Fitting, 1/16" NPT to 1/16" Hose Barb N2S231 PVDF Straight Fitting, 1/16" NPT to 1/16" Hose Barb



N2S330 N2S331

Polypropylene Straight Fitting, 1/16" NPT to 3/32" Hose Barb White Nylon Straight Fitting, 1/16" NPT to 3/32" Hose Barb Black Nylon Straight Fitting, 1/16" NPT to 3/32" Hose Barb PVDF Straight Fitting, 1/16" NPT to 3/32" Hose Barb



N2S430 N2S431

Polypropylene Straight Fitting, 1/16" NPT to 1/8" Hose Barb White Nylon Straight Fitting, 1/16" NPT to 1/8" Hose Barb Black Nylon Straight Fitting, 1/16" NPT to 1/8" Hose Barb PVDF Straight Fitting, 1/16" NPT to 1/8" Hose Barb



N2S530 N2S531 N2S570

Polypropylene Straight Fitting, 1/16" NPT to 5/32" Hose Barb White Nylon Straight Fitting, 1/16" NPT to 5/32" Hose Barb Black Nylon Straight Fitting, 1/16" NPT to 5/32" Hose Barb PVDF Straight Fitting, 1/16" NPT to 5/32" Hose Barb





Polypropylene Straight Fitting, 1/8" NPT to 1/16" Hose Barb White Nylon Straight Fitting, 1/8" NPT to 1/16" Hose Barb Black Nylon Straight Fitting, 1/8" NPT to 1/16" Hose Barb PVDF Straight Fitting, 1/8" NPT to 1/16" Hose Barb

# Thread to Tube Fittings - Straight (cont.)

DESCRIPTION

Polypropylene Straight Fitting, 1/8" NPT to 3/32" Hose Barb N4S330 White Nylon Straight Fitting, 1/8" NPT to 3/32" Hose Barb Black Nylon Straight Fitting, 1/8" NPT to 3/32" Hose Barb PVDF Straight Fitting, 1/8" NPT to 3/32" Hose Barb

Polypropylene Straight Fitting, 1/8" NPT to 1/8" Hose Barb

White Nylon Straight Fitting, 1/8" NPT to 1/8" Hose Barb

Black Nylon Straight Fitting, 1/8" NPT to 1/8" Hose Barb



N4S530 N4S531

Polypropylene Straight Fitting, 1/8" NPT to 5/32" Hose Barb White Nylon Straight Fitting, 1/8" NPT to 5/32" Hose Barb Black Nylon Straight Fitting, 1/8" NPT to 5/32" Hose Barb PVDF Straight Fitting, 1/8" NPT to 5/32" Hose Barb

Polypropylene Straight Fitting, 1/8" NPT to 3/16" Hose Barb

PVDF Straight Fitting, 1/8" NPT to 1/8" Hose Barb



N4S630 White Nylon Straight Fitting, 1/8" NPT to 3/16" Hose Barb N4S631 Black Nylon Straight Fitting, 1/8" NPT to 3/16" Hose Barb PVDF Straight Fitting, 1/8" NPT to 3/16" Hose Barb



N4S8 Polypropylene Straight Fitting, 1/8" NPT to 1/4" Hose Barb White Nylon Straight Fitting, 1/8" NPT to 1/4" Hose Barb N4S830 N4S831 Black Nylon Straight Fitting, 1/8" NPT to 1/4" Hose Barb PVDF Straight Fitting, 1/8" NPT to 1/4" Hose Barb



N8S470

Polypropylene Straight Fitting, 1/4" NPT to 1/8" Hose Barb White Nylon Straight Fitting, 1/4" NPT to 1/8" Hose Barb Black Nylon Straight Fitting, 1/4" NPT to 1/8" Hose Barb PVDF Straight Fitting, 1/4" NPT to 1/8" Hose Barb



Polypropylene Straight Fitting, 1/4" NPT to 5/32" Hose Barb White Nylon Straight Fitting, 1/4" NPT to 5/32" Hose Barb N8S530 N8S531 Black Nylon Straight Fitting, 1/4" NPT to 5/32" Hose Barb PVDF Straight Fitting, 1/4" NPT to 5/32" Hose Barb



N8S631 N8S670

Polypropylene Straight Fitting, 1/4" NPT to 3/16" Hose Barb White Nylon Straight Fitting, 1/4" NPT to 3/16" Hose Barb Black Nylon Straight Fitting, 1/4" NPT to 3/16" Hose Barb PVDF Straight Fitting, 1/4" NPT to 3/16" Hose Barb



N8S830 N8S831

Polypropylene Straight Fitting, 1/4" NPT to 1/4" Hose Barb White Nylon Straight Fitting, 1/4" NPT to 1/4" Hose Barb Black Nylon Straight Fitting, 1/4" NPT to 1/4" Hose Barb PVDF Straight Fitting, 1/4" NPT to 1/4" Hose Barb



N8S1030 N8S1031

Polypropylene Straight Fitting, 1/4" NPT to 5/16" Hose Barb White Nylon Straight Fitting, 1/4" NPT to 5/16" Hose Barb Black Nylon Straight Fitting, 1/4" NPT to 5/16" Hose Barb PVDF Straight Fitting, 1/4" NPT to 5/16" Hose Barb



N8S1230 N8S1231

Polypropylene Straight Fitting, 1/4" NPT to 3/8" Hose Barb White Nylon Straight Fitting, 1/4" NPT to 3/8" Hose Barb Black Nylon Straight Fitting, 1/4" NPT to 3/8" Hose Barb PVDF Straight Fitting, 1/4" NPT to 3/8" Hose Barb

# Thread to Tube Fittings - Elbow





DESCRIPTION ME230

Polypropylene Elbow Fitting, Tapered 10-32 to 1/16" Hose Barb White Nylon Elbow Fitting, Tapered 10-32 to 1/16" Hose Barb Black Nylon Elbow Fitting, Tapered 10-32 to 1/16" Hose Barb PVDF Elbow Fitting, Tapered 10-32 to 1/16" Hose Barbs



ME330 MF331

Polypropylene Elbow Fitting, Tapered 10-32 to 3/32" Hose Barb White Nylon Elbow Fitting, Tapered 10-32 to 3/32" Hose Barb Black Nylon Elbow Fitting, Tapered 10-32 to 3/32" Hose Barb PVDF Elbow Fitting, Tapered 10-32 to 3/32" Hose Barbs

# Thread to Tube Fittings - Elbow (cont)



ME430

Polypropylene Elbow Fitting, Tapered 10-32 to 1/8" Hose Barb White Nylon Elbow Fitting, Tapered 10-32 to 1/8" Hose Barb Black Nylon Elbow Fitting, Tapered 10-32 to 1/8" Hose Barb PVDF Elbow Fitting, Tapered 10-32 to 1/8" Hose Barbs

Polypropylene Elbow Fitting, Tapered 10-32 to 5/32" Hose Barb



ME530 White Nylon Elbow Fitting, Tapered 10-32 to 5/32" Hose Barb ME531 Black Nylon Elbow Fitting, Tapered 10-32 to 5/32" Hose Barb MF570 PVDF Elbow Fitting, Tapered 10-32 to 5/32" Hose Barbs Polypropylene Elbow Fitting, 1/8" NPT to 1/8" Hose Barb NAFA



N4F470 PVDF Elbow Fitting, 1/8" NPT to 1/8" Hose Barb Polypropylene Elbow Fitting, 1/8" NPT to 3/16" Hose Barb White Nylon Elbow Fitting, 1/8" NPT to 3/16" Hose Barb N4E631 Black Nylon Elbow Fitting, 1/8" NPT to 3/16" Hose Barb

PVDF Elbow Fitting, 1/8" NPT to 3/16" Hose Barb

White Nylon Elbow Fitting, 1/8" NPT to 1/8" Hose Barb

Black Nylon Elbow Fitting, 1/8" NPT to 1/8" Hose Barb



Polypropylene Elbow Fitting, 1/8" NPT to 1/4" Hose Barb N4E830 White Nylon Elbow Fitting, 1/8" NPT to 1/4" Hose Barb N4F831 Black Nylon Elbow Fitting, 1/8" NPT to 1/4" Hose Barb N4F870 PVDF Elbow Fitting, 1/8" NPT to 1/4" Hose Barb



Polypropylene Elbow Fitting, 1/4" NPT to 1/4" Hose Barb White Nylon Elbow Fitting, 1/4" NPT to 1/4" Hose Barb N8F830 Black Nylon Elbow Fitting, 1/4" NPT to 1/4" Hose Barb N8F831 N8F870 PVDF Elbow Fitting, 1/4" NPT to 1/4" Hose Barb



Polypropylene Elbow Fitting, 1/4" NPT to 3/8" Hose Barb White Nylon Elbow Fitting, 1/4" NPT to 3/8" Hose Barb Black Nylon Elbow Fitting, 1/4" NPT to 3/8" Hose Barb PVDF Elbow Fitting, 1/4" NPT to 3/8" Hose Barb

# Thread to Tube Fittings - Tee



DESCRIPTION MT230

MT231

MT531

Polypropylene Tee Fitting, Tapered 10-32 to 1/16" Hose Barbs White Nylon Tee Fitting, Tapered 10-32 to 1/16" Hose Barbs Black Nylon Tee Fitting, Tapered 10-32 to 1/16" Hose Barbs PVDF Tee Fitting, Tapered 10-32 to 1/16" Hose Barbs



Polypropylene Tee Fitting, Tapered 10-32 to 3/32" Hose Barbs MT330 White Nylon Tee Fitting, Tapered 10-32 to 3/32" Hose Barbs MT331 Black Nylon Tee Fitting, Tapered 10-32 to 3/32" Hose Barbs PVDF Tee Fitting, Tapered 10-32 to 3/32" Hose Barbs



Black Nylon Tee Fitting, Tapered 10-32 to 1/8" Hose Barbs PVDF Tee Fitting, Tapered 10-32 to 1/8" Hose Barbs MT530

Polypropylene Tee Fitting, Tapered 10-32 to 5/32" Hose Barbs White Nylon Tee Fitting, Tapered 10-32 to 5/32" Hose Barbs Black Nylon Tee Fitting, Tapered 10-32 to 5/32" Hose Barbs PVDF Tee Fitting, Tapered 10-32 to 5/32" Hose Barbs

Polypropylene Tee Fitting, Tapered 10-32 to 1/8" Hose Barbs

White Nylon Tee Fitting, Tapered 10-32 to 1/8" Hose Barbs

The unique shrouded hose barb on our 1/16" fittings provides additional stability to prevent breakage in applications where small tubing connections might be pulled or twisted during use.

Note: Parts not shown to scale. For customer drawings and technical specifications please visit cpcworldwide.com

# **BLOOD PRESSURE**

# Blood Pressure Fittings - Male Bayonet

PART NO. DESCRIPTION

X	BC330 BC331	White Nylon Male Bayonet Fitting to 3/32" Hose Barb Black Nylon Male Bayonet Fitting to 3/32" Hose Barb
D	BC430 BC431	White Nylon Male Bayonet Fitting to 1/8" Hose Barb Black Nylon Male Bayonet Fitting to 1/8" Hose Barb
×	BC530 BC531	White Nylon Male Bayonet Fitting to 5/32" Hose Barb Black Nylon Male Bayonet Fitting to 5/32" Hose Barb
D	BC630 BC631	White Nylon Male Bayonet Fitting to 3/16" Hose Barb Black Nylon Male Bayonet Fitting to 3/16" Hose Barb

# Blood Pressure Fittings - Screw-type

	P	
 • • • • • •	• • • • •	

DESCRIPTION GS320

White Acetal Screw-type Fitting, 5/16" UNF to 3/32" Hose Barb GS420 White Acetal Screw-type Fitting, 5/16" UNF to 1/8" Hose Barb White Acetal Screw-type Fitting, 5/16" UNF to 5/32" Hose Barb GS520



GM430M2

FOR SMC AND SMF1 BLOOD

PRESSURE CONNECTORS,

SEE PAGE 18-21.

White Acetal Screw-type Fitting, 5/16" UNF to 3/16" Hose Barb White Nylon with Polypropylene and TPE rotating nut. Buna-N seal. 5/16" UNF rotating internal thread X 1/8" HB White Nylon with Polypropylene and TPE rotating nut. Buna-N seal. 5/16" UNF rotating internal thread X 5/32" HB White Nylon with Polypropylene and TPE rotating nut, Buna-N seal. 5/16" UNF rotating internal thread X 3/16" HB

# LUERS

# Male Luers



**DESCRIPTION** 



Polypropylene Male Luer to 1/16" Hose Barb White Nylon Male Luer to 1/16" Hose Barb LM2130 LM2131 Black Nylon Male Luer to 1/16" Hose Barb LM2170 PVDF Male Luer to 1/16" Hose Barb Purple Tint Polycarbonate Male Luer to 1/16" Hose Barb LM2191

Male luers feature an integral locking ring!



LM31 Polypropylene Male Luer to 3/32" Hose Barb LM3130 White Nylon Male Luer to 3/32" Hose Barb LM3131 Black Nylon Male Luer to 3/32" Hose Barb PVDF Male Luer to 3/32" Hose Barb LM3170 Purple Tint Polycarbonate Male Luer to 3/32" Hose Barb LM3191



LM4130 LM4131 IM4170 I M4191

Polypropylene Male Luer to 1/8" Hose Barb White Nylon Male Luer to 1/8" Hose Barb Black Nylon Male Luer to 1/8" Hose Barb PVDF Male Luer to 1/8" Hose Barb Purple Tint Polycarbonate Male Luer to 1/8" Hose Barb



I M51 LM5130 LM5131 LM5170 LM5191

Polypropylene Male Luer to 5/32" Hose Barb White Nylon Male Luer to 5/32" Hose Barb Black Nylon Male Luer to 5/32" Hose Barb PVDF Male Luer to 5/32" Hose Barb Purple Tint Polycarbonate Male Luer to 5/32" Hose Barb



LM6130 LM6131 LM6170 LM6191 LM8170

Polypropylene Male Luer to 3/16" Hose Barb White Nylon Male Luer to 3/16" Hose Barb Black Nylon Male Luer to 3/16" Hose Barb PVDF Male Luer to 3/16" Hose Barb Purple Tint Polycarbonate Male Luer to 3/16" Hose Barb PVDF Male Luer to 1/4" Hose Barb

#### Female Luers

**DESCRIPTION** 



Polypropylene Female Luer to 1/16" Hose Barb White Nylon Female Luer to 1/16" Hose Barb Black Nylon Female Luer to 1/16" Hose Barb PVDF Female Luer to 1/16" Hose Barb Purple Tint Polycarbonate Female Luer to 1/16" Hose Barb



LF31 LF3130 LF3131 LF3170 LF3191

Polypropylene Female Luer to 3/32" Hose Barb White Nylon Female Luer to 3/32" Hose Barb Black Nylon Female Luer to 3/32" Hose Barb PVDF Female Luer to 3/32" Hose Barb Purple Tint Polycarbonate Female Luer to 3/32" Hose Barb



LF4130 LF4131 LF4191

Polypropylene Female Luer to 1/8" Hose Barb White Nylon Female Luer to 1/8" Hose Barb Black Nylon Female Luer to 1/8" Hose Barb PVDF Female Luer to 1/8" Hose Barb Purple Tint Polycarbonate Female Luer to 1/8" Hose Barb



LF51 LF5130 LF5131 LF5170 LF5191

Polypropylene Female Luer to 5/32" Hose Barb White Nylon Female Luer to 5/32" Hose Barb Black Nylon Female Luer to 5/32" Hose Barb PVDF Female Luer to 5/32" Hose Barb Purple Tint Polycarbonate Female Luer to 5/32" Hose Barb



LF61 LF6130 LF6131 LF6170 LF6191

Polypropylene Female Luer to 3/16" Hose Barb White Nylon Female Luer to 3/16" Hose Barb Black Nylon Female Luer to 3/16" Hose Barb PVDF Female Luer to 3/16" Hose Barb Purple Tint Polycarbonate Female Luer to 3/16" Hose Barb



LF8130 LF8131 LF8170

Polypropylene Female Luer to 1/4" Hose Barb White Nylon Female Luer to 1/4" Hose Barb Black Nylon Female Luer to 1/4" Hose Barb PVDF Female Luer to 1/4" Hose Barb Purple Tint Polycarbonate Female Luer to 1/4" Hose Barb

#### Panel Mount Male Luers



PMLM21 Polypropylene Panel Mount Male Luer to 1/16" Hose Barb PMLM2130 White Nylon Panel Mount Male Luer to 1/16" Hose Barb PMLM2131 Black Nylon Panel Mount Male Luer to 1/16" Hose Barb PMLM2170 PVDF Panel Mount Male Luer to 1/16" Hose Barb



Polypropylene Panel Mount Male Luer to 3/32" Hose Barb PMLM3130 White Nylon Panel Mount Male Luer to 3/32" Hose Barb PMLM3131 Black Nylon Panel Mount Male Luer to 3/32" Hose Barb PMLM3170 PVDF Panel Mount Male Luer to 3/32" Hose Barb



PMLM41 Polypropylene Panel Mount Male Luer to 1/8" Hose Barb PMLM4130 White Nylon Panel Mount Male Luer to 1/8" Hose Barb PMLM4131 Black Nylon Panel Mount Male Luer to 1/8" Hose Barb PMLM4170 PVDF Panel Mount Male Luer to 1/8" Hose Barb



PMNUTSST Stainless Steel 1/4-28 UNF Nut PMNUT30 White Nylon Panel Mount 1/4-28 UNF Nut

**Note:** Panel mount luers require a ½" diameter panel opening



### Panel Mount Female Luers



PART NO. DESCRIPTION Polypropylene Panel Mount Female Luer to 1/16" Hose Barb PMLF2130 White Nylon Panel Mount Female Luer to 1/16" Hose Barb PMLF2131 Black Nylon Panel Mount Female Luer to 1/16" Hose Barb PMLF2170 PVDF Panel Mount Female Luer to 1/16" Hose Barb



Polypropylene Panel Mount Female Luer to 3/32" Hose Barb PMLF3130 White Nylon Panel Mount Female Luer to 3/32" Hose Barb PMLF3131 Black Nylon Panel Mount Female Luer to 3/32" Hose Barb PMLF3170 PVDF Panel Mount Female Luer to 3/32" Hose Barb



Polypropylene Panel Mount Female Luer to 1/8" Hose Barb PMLF4130 White Nylon Panel Mount Female Luer to 1/8" Hose Barb PMLF4131 Black Nylon Panel Mount Female Luer to 1/8" Hose Barb PMLF4170 PVDF Panel Mount Female Luer to 1/8" Hose Barb



PMNUTSST Stainless Steel 1/4-28 UNF Nut PMNUT30 White Nylon Panel Mount 1/4-28 UNF Nut

# Luer Caps and Plugs



LCAP30 LCAP31

**DESCRIPTION** Polypropylene Luer Cap White Nylon Luer Can Black Nylon Luer Cap PVDF Luer Cap



I PI IIG LPLUG30 LPLUG31 LPLUG70

Polypropylene Luer Plug White Nylon Luer Plug Black Nylon Luer Plug PVDF Luer Plug

rings can be used for visual line identification or to prevent connector rotation in panel mount applications

Color-coded

# Male Slip Luers

SLM2130

Polypropylene Male Slip Luer to 1/16" Hose Barb White Nylon Male Slip Luer to 1/16" Hose Barb Black Nylon Male Slip Luer to 1/16" Hose Barb PVDF Male Slip Luer to 1/16" Hose Barb Purple Tint Polycarbonate Male Slip Luer to 1/16" Hose Barb



Polypropylene Male Slip Luer to 3/32" Hose Barb White Nylon Male Slip Luer to 3/32" Hose Barb Black Nylon Male Slip Luer to 3/32" Hose Barb SI M3170 PVDF Male Slip Luer to 3/32" Hose Barb Purple Tint Polycarbonate Male Slip Luer to 3/32" Hose Barb



Polypropylene Male Slip Luer to 1/8" Hose Barb White Nylon Male Slip Luer to 1/8" Hose Barb SI M4130 Black Nylon Male Slip Luer to 1/8" Hose Barb SI M4170 PVDF Male Slip Luer to 1/8" Hose Barb SLM4191 Purple Tint Polycarbonate Male Slip Luer to 1/8" Hose Barb



Polypropylene Male Slip Luer to 5/32" Hose Barb SLM51 SLM5130 White Nylon Male Slip Luer to 5/32" Hose Barb Black Nylon Male Slip Luer to 5/32" Hose Barb PVDF Male Slip Luer to 5/32" Hose Barb SLM5191 Purple Tint Polycarbonate Male Slip Luer to 5/32" Hose Barb




Polypropylene Male Slip Luer to 3/16" Hose Barb SLM6130 White Nylon Male Slip Luer to 3/16" Hose Barb Black Nylon Male Slip Luer to 3/16" Hose Barb SLM6131 PVDF Male Slip Luer to 3/16" Hose Barb Purple Tint Polycarbonate Male Slip Luer to 3/16" Hose Barb

### Luer Accessories



Polypropylene Stationary Luer Lock Ring LMSL30 White Nylon Stationary Luer Lock Ring LMSL31 Black Nylon Stationary Luer Lock Ring Yellow Nylon Stationary Luer Lock Ring LMSL32 LMSL33 Blue Nylon Stationary Luer Lock Ring LMSI34 Red Nylon Stationary Luer Lock Ring







LMSR Polypropylene Rotating Luer Lock Ring LMSR30 White Nylon Rotating Luer Lock Ring Black Nylon Rotating Luer Lock Ring Yellow Nylon Rotating Luer Lock Ring LMSR32 Blue Nylon Rotating Luer Lock Ring Red Nylon Rotating Luer Lock Ring LMSR34



Green Nylon Rotating Luer Lock Ring LMSR36 Orange Nylon Rotating Luer Lock Ring LMSR70 PVDF Rotating Luer Lock Ring Purple Tint Polycarbonate Rotating Luer Lock Ring Polypropylene Luer Color Coded Ring



PMRL30 White Nylon Luer Color Coded Ring PMRI 31 Black Nylon Luer Color Coded Ring Yellow Nylon Luer Color Coded Ring PMRI 32 PMRI 33 Blue Nylon Luer Color Coded Ring Red Nylon Luer Color Coded Ring Green Nylon Luer Color Coded Ring PMRI 35 PMRI 36 Orange Nylon Luer Color Coded Ring PVDF Panel Mount Luer Color Coded Ring




PMNUTSST Stainless Steel 1/4-28 UNF Nut


Note: Parts not shown to scale. For customer drawings and technical specifications please visit cpcworldwide.com

# **ACCESSORIES**

### **Accessories**



#### Manifolds • SPECIFICATIONS • Material: Aluminum • Finish: Black anodized



HFC312 (no leash)

### JG® Push-To-Connect Accessories

DESCRIPTION IN-LINE LOCKING CLIPS, RED	<b>TUBING</b> 1/4" Tube OD 3/8" Tube OD	<b>PART NO.</b> 1011900 1012000	<b>DESCRIPTION</b> In-Line Collet Covers, Color Coded	<b>TUBING</b> 1/4" Tube OD 3/8" Tube OD	<b>PART NO.</b> 1012101 to 06 1012201 to 06	
	_	0 0	suffix codes: 01=Black, 02=Red, See pages 42 and 43 for JG coupl	, , , , , , , , , , , , , , , , , , , ,	05=Blue, 06=Green	

**†NOTE**: CPC's Ferruleless Polytube Fitting terminations do not require ferrules to achieve a secure connection, which makes them easier to use and reuse. PTF fittings are designed for semi-rigid tubing, i.e., polyethylene, nylon, polyurethane, etc. \*NOTE: Gaskets are not included with 10-32 PTF style fittings. To prevent leakage with these fittings, we recommend using LOCTITE® Removable Threadlocker 242 or Permanent Threadlocker 262.

# STERILIZATION AND DISINFECTANT METHODS

		METHOD							
KEY	Y = Excellent, recommended material for this sterilization method  N = No, not recommended  N/A = Not Applicable	FORMALIN	NFECTAI	ETHYL ALCOHOL	ETHYLENE OXIDE	AUTOCLAVE	E-BEAM IRRADIATION SO WITDGE WAYS	GAMMA IRRADIATION SAWAGE	DRY HEAT
MATERIAL	METALS 302 STAINLESS STEEL 316 STAINLESS STEEL CHROME-PLATED BRASS-CDA 360	Y Y N	Y Y	Y Y Y	Y Y Y	Y Y Y	Y Y Y	ү ү ү	y Y Y
	POLYMERS  ABS  ACETAL  LDPE  NYLON  POLYCARBONATE  POLYPROPYLENE (excludes)  POLYSULFONE	N Y Y Y Y Y	N/A Y Y N Y Y Y Y Y	Y Y Y N Y Y	Y Y Y Y Y Y	N Y N N Y N	Y N Y N Y Y	Y N Y N Y Y	N N N N Y N
	ELASTOMERS NITRILE/BUNA-N SILICONE KALREZ® EPR/EPDM FKM/VITON®	Y Y Y Y	Y Y Y Y	Y Y Y Y	N Y Y Y N	N Y Y Y N	Y Y N Y	Y Y N Y	Y Y Y Y

### Sterilization Methods

**DISINFECTANTS:** 70°F (20°C), Formalin, ethyl alcohol, etc. Sterilize coupled or uncoupled.

**ETHYLENE OXIDE, ETO:** Four hours, 100% Et0 @ 110°F (43°C), up to five repetitions, coupled or uncoupled.

#### **AUTOCLAVE:**

Acetal:  $250^{\circ}$ F ( $121^{\circ}$ C), 30 minutes max, up to 10 repetitions. Sterilize uncoupled only.

**Polycarbonate**: 250°F (121°C), 30 minutes, up to 10 repetitions. Sterilize uncoupled only.

**Polysulfone**: 270°F (132°C) for 60 minutes, up to 25 repetitions. Sterilize uncoupled only.

**Steam-Thru Connection**: 265°F (129°C) for 30 minutes, up to two cycles. (part number specific)

 $\mbox{HFC39}\mbox{:}~270\mbox{\,°F}~(132\mbox{\,°C})$  for 60 minutes, up to 25 repetitions for uncoupled units and up to one repetition for coupled units.

**Electron Beam (E-Beam)**: Maximum cumulative exposure of 50 kilograys. Sterilize coupled or uncoupled.

**Gamma:** Maximum cumulative exposure of 50 kilograys. Sterilize coupled or uncoupled.

**Dry Heat**:  $250^{\circ}$ F ( $121^{\circ}$ C), 12 hours, no pressure. Sterilize uncoupled only.

#### SIP process:

Up to 266°F (130°C) for 60 minutes (Steam-Thru Connection)

Up to 275°F (135°C) for 60 minutes (Steam-Thru II Connection)

Unsure of
the chemical
compatibility for your
application? Consult
the ChemCheck tool on
cpcworldwide.com or contac
CPC customer service to
receive product samples
for in-house testing
and verification.

Smart fluid handling to take you forward, faster. 163

# REGULATORY AND COMPLIANCE

#### ISO 13485:2003 Certification

ISO 13485:2003 is recognized by regulators around the world as a good basis for addressing medical device design and manufacturing regulatory requirements. It allows us to enhance product safety by proactively identifying and managing product and project risks. Becoming ISO 13485:2003 certified has allowed us to better control the consistency of manufactured products.

### ISO 9001:2008 Certification

ISO 9001:2008 is a standard which assures consistency of a product ordered by customers. Organizations having ISO 9001:2008 certification have demonstrated compliance to the ISO 9001:2008 requirements by an independent registration authority. CPC's Quality Management System has been approved and certified under the ISO 9001:2008 standard.

# Cleanroom Manufacturing

CPC manufactures certain Life Sciences and Chemical Management product lines in a cleanroom certified by an external testing service to meet or exceed ISO Class 7 (10,000) at 0.5 mm per ISO 14644 and the former Federal Standard 209E. Certification data is available upon request.

#### FDA and USDA

The U.S. Food and Drug Administration publishes, through the Code of Federal Regulations, standardized criteria which govern the acceptability of materials used in food contact. The U.S. Department of Agriculture publishes similar standards that mirror FDA criteria. Neither agency approves or disapproves products for particular applications. Most of CPC's products are made using resins that comply with applicable FDA or USDA standards. When necessary, the standard o-ring seals are replaced with specific, recognized materials.

#### REACH

REACH is the Regulation for Registration, Evaluation, Authorisation and Restriction of Chemicals. It entered into force on 1st June 2007 to streamline and improve the former legislative framework on chemicals of the European Union (EU). REACH places greater responsibility on industry to manage the risks that chemicals may pose to the health and the environment. CPC publishes a list of CPC products that are compliant with the EU regulation 1907/2006.

# Regulation of Hazardous Substances

The RoHS Directive stands for "the restriction of the use of certain hazardous substances in electrical and electronic equipment." This Directive bans the placing on the EU market of new electrical and electronic equipment containing more than agreed levels of lead, cadmium, mercury, hexavalent chromium, polybrominated biphenyl (PBB) and polybrominated diphenyl ether (PBDE) flame retardants.

#### NSF

NSF International, based in Ann Arbor, Michigan, develops and publishes consensual criteria that govern the acceptability of materials and equipment used in food and beverage processing. They also do testing to verify the performance of materials or devices to their published criteria. CPC lists many of its product lines under the criteria of NSF/ANSI Standards, which governs components used in food and beverage contact applications.





**COLDER PRODUCTS COMPANY 1001 WESTGATE DRIVE** ST. PAUL, MINNESOTA 55114 U.S.A.

PHONE: +1 651-645-0091 FAX. +1 651-645-5404 **TOLL FREE:** 800-444-2474 info@cpcworldwide.com cpcworldwide.com

COLDER PRODUCTS COMPANY GMBH **OSTRING 19** D-63762 GROSSOSTHEIM **GERMANY** 

PHONE: +49-6026-9973-0 +49-6026-9973-173 FAX:

cpcgmbh@cpcworldwide.com cpcworldwide.com

**COLDER PRODUCTS COMPANY LIMITED** FLAT B, 29/F, WEST GATE TOWER, 7 WING HONG STREET, CHEUNG SHA WAN KOWLOON, HONG KONG

PHONE: +852-2987-5272 +852-2987-2509 FAX:

asiapacific@cpcworldwide.com cpcworldwide.com







workmanship and materials for a period of 12 months from the date of sale by CPC to its initial customer (regardless of any subsequent sale of the products). This warranty is void if the product is misused, altered, tampered with or is installed or used in a manner that is inconsistent with CPC's written recommendations, specifications and/or instructions, or fails to perform due to normal wear and tear. CPC does not warrant the suitability of the product for any particular application. Determining product application suitability is solely the customer's responsibility. CPC is not liable for special, indirect, incidental, consequential or other damages including, but not limited to, loss, damage, personal injury, or any other expense directly or indirectly arising from the use of or inability to use its products either separately or in combination with other products. ALL OTHER WARRANTIES EXPRESS OR IMPLIED, WHETHER ORAL, WRITTEN OR IN ANY OTHER FORM, INCLUDING BUT NOT LIMITED TO WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ARE EXPRESSLY EXCLUDED. The sole and exclusive remedy under this warranty is limited, at the option of CPC, to replacement of the defective product or an account credit in the amount of the original selling price. All

CPC WARRANTY STATEMENT: CPC (Colder Products Company) warrants its products against defects in

describing the product's application and performance, unless otherwise authorized in writing by CPC. CPC PATENT STATEMENT: CPC takes pride in its innovative quick disconnect coupling and fittings solutions, many of which have been awarded United States and international patents. CPC has a strong tradition of

allegedly defective CPC products must be returned prepaid transportation to CPC, together with information

leadership in the quick disconnect market, and aggressively pursues and protects its proprietary information and intellectual property. In cases where it is practical and has a benefit to its customers, CPC has licensed its proprietary technology. Please contact CPC to discuss your unique needs.

CPC TRADEMARK STATEMENT: AseptiQuik®, BottleQuik®, BreakAway®, ChemQuik®, DrumQuik®, FitQuik®, IdentiQuik®, Nu-Seal®, SnapQuik®, Steam-Thru®, Softube® are registered trademarks with the U.S. Patent & Trademark Office. All other trademarks or service marks are property of their respective owners. WARNING: Due to the wide variety of possible fluid media and operating conditions, unintended consequences may result from the use of this product, all of which are beyond the control of CPC. It is the user's responsibility to carefully determine and test for compatibility for use with their application. All such risks shall be assumed by the buyer

COPYRIGHT © 2014 COLDER PRODUCTS COMPANY. All rights reserved. Colder Products Company, Colder Products and CPC are registered trademarks with the US Patent & Trademark Office.

WARNING: Due to the wide variety of possible fluid media and operating conditions, unintended consequences may result from the use of this product, all of which are beyond the control of CPC It is the user's responsibility to carefully determine and test for compatibility for use with their application. All such risks shall be assumed

CAT2020 12/15 WC/15K PRINTED IN U.S.A.